


Analysis of Sulfur-Containing Metabolites of Asparagus in Urine by SBSE-GCxGC-TOFMS

Life Science and Chemical Analysis Centre
LECO Corporation, Saint Joseph, Michigan USA



Columns:

10 m x 0.18 mm x 0.2 μm Rtx-5
1.0 m x 0.10 mm x 0.1 μm Rtx-17

TOFMS:

35 to 450 m/z @ 200 spectra/s

For questions on this analysis
e-mail us at: life_science@leco.com

For a complete listing of snapshots
and application notes visit us on the
web at www.leco.com

