



PAL3 Autosampler Series II Accessories & Smart Consumables



Content

| | |
|--|----|
| Tools & Syringes | 4 |
| Modules | 11 |
| Tray Holders, Stacks and Racks | 15 |
| Smart Syringe & Needle Portfolio | 19 |
| Smart SPME Arrow / Smart SPME Fiber | 26 |
| μSPE Options | 32 |
| Accessories | 34 |
| Consumables - Vials, Vial Caps and Seals | 37 |
| Spare Parts A-Z | 40 |



Tools & Syringes

Liquid Tools & Syringes

| Part No. | Description | Compatible with: | |
|---------------|---|------------------|-----|
| | | RTC | RSI |
| PAL3-TH-D757 | Universal Liquid Syringe Tool for syringes of 0.5µL, 1.2µL, 5µL, 10µL, 25µL, 50µL or 100µL with a needle length of 57mm Syringes not included, requires firmware version 1.3.0 or higher | X | X |
| PAL3-TH-D785 | Universal Liquid Syringe Tool for syringes of 0.5µL, 1.2µL, 5µL, 10µL, 25µL, 50µL or 100µL with a needle length of 85mm Syringes not included, requires firmware version 1.3.0 or higher | X | X |
| PAL3-TH-D857 | Universal Liquid Syringe Tool for syringes of 250µL, 500µL or 1000µL with a needle length of 57mm Syringes not included, requires firmware version 1.3.0 or higher | X | X |
| PAL3-TH-D885 | Universal Liquid Syringe Tool for syringes of 250µL, 500µL or 1000µL with a needle length of 85mm Syringes not included, requires firmware version 1.3.0 or higher | X | X |
| PAL3-TH-D1857 | Liquid Syringe Tool for a syringe of 2500µL, 5000µL, 10000µL with a needle length of 57mm Syringe not included, requires firmware version 1.3.0 or higher | X | X |



Recommended GC Syringes:

| Part No. | Description | Tool compatibility |
|-------------------|--|--------------------|
| SF0-57-B-26P-CO | 0.5µL Smart Syringe with fixed needle for Tool D7/57: Needle length 57mm, gauge 26P, scale length 27mm, Point Style conical | D7/57 |
| SE1-B-57-23P-CO | 1µL Smart Syringe body for exchangeable needles for Tool D8/57, Plunger in Needle, Needle length 57mm, gauge 23P, scale length 27mm, Point Style conical, Needle Flange Type NIP_1 | D8/57 |
| SF1-57-B-26P-CO | 1.2µL Smart Syringe with fixed needle for Tool D7/57: Needle length 57mm, gauge 26P, scale length 36mm, Point Style conical | D7/57 |
| SF1-85-B-26P-CO | 1.2µL Smart Syringe with fixed needle for Tool D7/85: Needle length 85mm, gauge 26P, scale length 36mm, Point Style conical | D7/85 |
| SF10-57-M-26S-CO | 10uL Smart Syringe with fixed needle for Tool D7/57: Needle length 57mm, metal plunger, gauge 26S, scale length 54mm, Point Style conical | D7/57 |
| SF10-57-M-23S-CO | 10uL Smart Syringe with fixed needle for Tool D7/57: Needle length 57mm, metal plunger, gauge 23S, scale length 54mm, Point Style conical | D7/57 |
| SF10-57-T-26S-CO | 10uL Smart Syringe with fixed needle for Tool D7/57: Needle length 57mm, PTFE plunger, gauge 26S, scale length 60mm, Point Style conical | D7/57 |
| SF10-57-T-23S-CO | 10uL Smart Syringe with fixed needle for Tool D7/57: Needle length 57mm, PTFE plunger, gauge 23S, scale length 60mm, Point Style conical | D7/57 |
| SF10-85-M-26S-CO | 10uL Smart Syringe with fixed needle for Tool D7/85: Needle length 85mm, metal plunger, gauge 26S, scale length 54mm, Point Style conical | D7/85 |
| SF10-85-T-26S-CO | 10uL Smart Syringe with fixed needle for Tool D7/85: Needle length 85mm, PTFE plunger, gauge 26S, scale length 60mm, Point Style conical | D7/85 |
| SF25-57-T-26S-CO | 25uL Smart Syringe with fixed needle for Tool D7/57: Needle length 57mm, PTFE plunger, gauge 26S, scale length 60mm, Point Style conical | D7/57 |
| SF25-85-T-26S-CO, | 25uL Smart Syringe with fixed needle for Tool D7/85: Needle length 85mm, PTFE plunger, gauge 26S, scale length 60mm, Point Style conical | D7/85 |
| SF50-57-T-26S-CO | 50uL Smart Syringe with fixed needle for Tool D7/57: Needle length 57mm, PTFE plunger, gauge 26S, scale length 60mm, Point Style conical | D7/57 |
| SF50-85-T-26S-CO | 50uL Smart Syringe with fixed needle for Tool D7/85: Needle length 85mm, PTFE plunger, gauge 26S, scale length 60mm, Point Style conical | D7/85 |

CAUTION: SYRINGES FOR LABORATORY USE ONLY - NOT FOR HUMAN IN VIVO USE

| Part No. | Description | Tool compatibility |
|-------------------|---|--------------------|
| SF100-57-T-23S-CO | 100uL Smart Syringe with fixed needle for Tool D7/57: Needle length 57mm, PTFE plunger, gauge 23S, scale length 60mm, Point Style conical | D7/57 |
| SF100-57-T-26S-CO | 100uL Smart Syringe with fixed needle for Tool D7/57: Needle length 57mm, PTFE plunger, gauge 26S, scale length 60mm, Point Style conical | D7/57 |
| SF100-85-T-26S-CO | 100uL Smart Syringe with fixed needle for Tool D7/85: Needle length 85mm, PTFE plunger, gauge 26S, scale length 60mm, Point Style conical | D7/85 |
| SF250-57-T-26-CO | 250uL Smart Syringe with fixed needle for Tool D8/57: Needle length 57mm, PTFE plunger, gauge 26, scale length 60mm, Point Style conical | D8/57 |
| SF250-85-T-26-CO | 250uL Smart Syringe with fixed needle for Tool D8/85: Needle length 85mm, PTFE plunger, gauge 26, scale length 60mm, Point Style conical | D8/85 |
| SF500-57-T-26-CO | 500uL Smart Syringe with fixed needle for Tool D8/57: Needle length 57mm, PTFE plunger, gauge 26, scale length 60mm, Point Style conical | D8/57 |
| SF1000-57-T-23-CO | 1000uL Smart Syringe with fixed needle for Tool D8/57: Needle length 57mm, PTFE plunger, gauge 23, scale length 60mm, Point Style conical | D8/57 |

Syringe Bodies and exchangeable Needles can be found on page 22-25

Recommended LC Syringes

| Part No. | Description | Tool compatibility |
|--------------------|--|--------------------|
| SF10-57-T-22S-FL | 10uL Smart Syringe with fixed needle for Tool D7/57: Needle length 57mm, PTFE plunger, gauge 22S, scale length 60mm, Point Style flat | D7/57 |
| SFX25-57-T-22S-FL | 25uL Smart Syringe with fixed needle for Tool D7/57: Needle length 57mm, PTFE plunger, gauge 22S, scale length 60mm, Point Style flat, x-Type | D7/57 |
| SFX100-57-T-22-FL | 100uL Smart Syringe with fixed needle for Tool D7/57: Needle length 57mm, PTFE plunger, gauge 22, scale length 60mm, Point Style flat, x-Type | D7/57 |
| SFX100-57-T-22S-FL | 100uL Smart Syringe with fixed needle for Tool D7/57: Needle length 57mm, PTFE plunger, gauge 22S, scale length 60mm, Point Style flat, x-Type | D7/57 |
| SF250-57-T-22-FL | 250uL Smart Syringe with fixed needle for Tool D8/57: Needle length 57mm, PTFE plunger, gauge 22, scale length 60mm, Point Style flat | D8/57 |
| SF500-57-T-22-FL | 500uL Smart Syringe with fixed needle for Tool D8/57: Needle length 57mm, PTFE plunger, gauge 22, scale length 60mm, Point Style flat | D8/57 |
| SF1000-57-T-22-FL | 1000uL Smart Syringe with fixed needle for Tool D8/57: Needle length 57mm, PTFE plunger, gauge 22, scale length 60mm, Point Style flat | D8/57 |

Syringe Bodies and exchangeable Needles can be found on page 22-25

Recommended Prep Syringes

| Part No. | Description | Tool compatibility |
|--------------------|---|--------------------|
| SC2500-57-T-19-FL | 2500uL Smart Syringe with fixed needle for Tool D18/57: Needle length 57mm, PTFE plunger, gauge 19, scale length 60mm, Point Style flat | D18/57 |
| SC2500-57-T-22-FL | 2500uL Smart Syringe with fixed needle for Tool D18/57: Needle length 57mm, PTFE plunger, gauge 22, scale length 60mm, Point Style flat | D18/57 |
| SC2500-57-PE-22-FL | 2500uL Smart Syringe with fixed needle for Tool D18/57: Needle length 57mm, Polyethylene plunger, gauge 22, scale length 60mm, Point Style flat | D18/57 |
| SC5000-57-T-19-FL | 5000uL Smart Syringe with fixed needle for Tool D18/57: Needle length 57mm, PTFE plunger, gauge 19, scale length 60mm, Point Style flat | D18/57 |
| SC5000-57-T-22-FL | 5000uL Smart Syringe with fixed needle for Tool D18/57: Needle length 57mm, PTFE plunger, gauge 22, scale length 60mm, Point Style flat | D18/57 |

Syringe Bodies and exchangeable Needles can be found on page 22-25

CAUTION: SYRINGES FOR LABORATORY USE ONLY - NOT FOR HUMAN IN VIVO USE

Headspace Tools & Syringes

| Part No. | Description | Compatible with: | |
|----------------|--|------------------|-----|
| | | RTC | RSI |
| PAL3-TH-HS1000 | Headspace Syringe Tool for a 1000µL syringe with a needle length of 65mm Syringe not included | X | X |
| PAL3-TH-HS2500 | Headspace Syringe Tool for a 2500µL syringe with a needle length of 65mm Syringe not included | X | X |
| PAL3-TH-HS5000 | Headspace Syringe Tool for a 5000µL syringe with a needle length of 65mm Syringe not included | X | X |




Recommended GC Syringes:

| Part No. | Description | Tool compatibility |
|------------------------|--|--------------------|
| SH1000-65-T-23-SP | 1000µL Smart Headspace Syringe with fixed needle for Tool HS1000, PTFE Plunger, Needle length 65mm, Gauge 23, Point Style Sideport | HS1000 |
| SH2500-65-T-23-SP | 2500µL Smart Headspace Syringe with fixed needle for Tool HS2500, PTFE Plunger, Needle length 65mm, Gauge 23, Point Style Sideport | HS2500 |
| SH5000-65-T-23-SP | 5000µL Smart Headspace Syringe with fixed needle for Tool HS5000, PTFE Plunger, Needle length 65mm, Gauge 23, Point Style Sideport | HS5000 |
| Gluefree option | | |
| SHGF2250-T-23-SP | 2250µl Smart Glue Free Headspace Syringe with fixed needle for Tool HS2500, PTFE plunger, needle gauge 23, Point Style Sideport | HS2500 |
| SHGF4500-T-23-SP | 4500µl Smart Glue Free Headspace Syringe with fixed needle for Tool HS5000, PTFE plunger, needle gauge 23, Point Style Sideport | HS5000 |

Syringe Bodies and exchangeable Needles can be found on page 22-25

LCMS Tool

| Part No. | Description | Compatible with: | |
|-----------------|--|------------------|-----|
| | | RTC | RSI |
| PAL3-TH-SLCMS-P | <p>LCMS-P Tool consisting of:</p> <ul style="list-style-type: none"> 1 pc LCMS-P Tool incl. Smart Syringe, Needle and FEP Loop 1 pc Pump Module and Connection Cable 1 pc PFA Transfer tube 1/16" 2 pc Solvent bottle 1000ml 2 pc Aspirating tube 1/8" <p>Attention: LCMS Wash Module not included. Requires Check Valve Kit. Compatible with all PAL RTC/PAL RSI models from FW version 2.1 or higher. Requires firmware version 3.0 or higher and PAL RTC/PAL RSI Series II for Smart functionality</p> | X | X |



Recommended Syringes and Needles

| Part No. | Description | Tool compatibility |
|-----------------|--|--------------------|
| SLCMS-100-PE | 100µl Smart LCMS Syringe for LCMS Tool, PE Plunger, Gastight, Scale Length 60 mm, Thread 1/4-28 UNF, Removable Needle not included | LCMS |
| SLCMS-100-T | 100µl Smart LCMS Syringe for LCMS Tool, PTFE Plunger, Gastight, Scale Length 60 mm, Thread 1/4-28 UNF, Removable Needle not included | LCMS |
| PAL3-LCMS-Ndl | Needle for PAL LCMS Tool, gauge 22 PST 3, length 51 mm (3 pcs per pack, incl. Needle Retaining Nut) | LCMS |
| PAL3-LCMS-Ndl-C | Coated Needle for PAL LCMS Tool, gauge 22 PST 3, length 51 mm (3 pcs per pack, incl. Needle Retaining Nut) | LCMS |
| PAL3-LCMS-Ndl-T | Titanium Needle for PAL LCMS Tool, gauge 22 PST 3, length 51 mm (incl. Needle Retaining Nut) | LCMS |

CAUTION: SYRINGES FOR LABORATORY USE ONLY - NOT FOR HUMAN IN VIVO USE

Pipette Tools

| Part No. | Description | Compatible with: | |
|-----------------|---|------------------|-----|
| | | RTC | RSI |
| PAL3-TH-Pip1000 | Pipette tool for 1000uL Tips To use with PAL Disposable Tips: PAL3-PipTip-1000 and PAL3-Pip-Tip-1000F Tips and Racks for tips are not included requires firmware version 2.4.x or higher and PAL RTC | X | |
| PAL3-TH-Pip200 | Pipette tool for 200uL Tips To use with PAL Disposable Tips: PAL3-PipTip-200 and PAL3-Pip-Tip-200F Tips and Racks for tips are not included requires firmware version 2.4.x or higher and PAL RTC | X | |



Recommended Pipette Tips:

| Part No. | Description | Tool compatibility |
|-------------------|--|--------------------|
| PAL3-PipTip-1000 | Pipette Tips 1000uL without Filter 10 pcs with 96 tips each | TH-Pip1000 |
| PAL3-PipTip-1000F | Pipette Tips 1000uL with Filter 10 pcs with 96 tips each | TH-Pip1000 |
| PAL3-PipTip-200 | Pipette Tips 200uL without Filter 10 pcs with 96 tips each | TH-Pip200 |
| PAL3-PipTip-200F | Pipette Tips 200uL with Filter 10 pcs with 96 tips each | TH-Pip200 |

Recommended Rack for Pipette Tips:

| Part No. | Description |
|-------------------|---|
| PAL3-Rack-Pip-Tip | Tip rack for 96 PAL pipette tips Tips not included |

Dilutor Tool

| Part No. | Description | Compatible with: | |
|-----------------|--|------------------|-----|
| | | RTC | RSI |
| PAL3-TH-Dilutor | Tool Holder for Dilutor Option Attention: Does not include needles and tubing requires firmware version 2.1.0 or higher and Dilutor Option | X | X |




Recommended syringes:

| Part No. | Description | Tool compatibility |
|------------------|---|--------------------|
| PAL3-SYDIL-21834 | 100µL Dilutor Syringe, Scale 60mm, PTFE Plunger, ZDV (Zero Dead Volume) | Dilutor |
| PAL3-SYDIL-21835 | 1mL Dilutor Syringe, Scale 60mm, PTFE Plunger, ZDV (Zero Dead Volume) | Dilutor |
| PAL3-SYDIL-21836 | 5mL Dilutor Syringe, Scale 60mm, PTFE Plunger, ZDV (Zero Dead Volume) | Dilutor |
| PAL3-SYDIL-21837 | 10mL Dilutor Syringe, Scale 60mm, PTFE Plunger, ZDV (Zero Dead Volume) | Dilutor |



Modules

| Part No. | Description | Compatible with: | | | |
|-----------------------|---|------------------|-----|-----|---|
| | | RTC | RSI | LSI | |
| PAL3-Agitator | Agitator for the incubation and agitation of up to 6x 20mL vials | X | X | X |  |
| PAL3-BarcodeReader | Barcode Reader with 2 active laser scanners for all standard vials containing 1 dimension (1-D) barcode in a horizontal orientation | X | X | X |  |
| PAL3-Centrifuge-2mL | Centrifuge for the centrifugation of up to 8x 2mL vials Maximum speed of 7200 rpm. this corresponds to a RCF of 5000 x g | X | X | X |  |
| PAL3-Centrifuge-Combi | Centrifuge for the centrifugation of 2x 10mL vials, 2x 20mL vials or up to 4x 2mL vials Maximum speed of 4800 rpm. this corresponds to a RCF of 2000 x g | X | X | X | |
| PAL3-Decapper | Decapper Module for opening and closing of initially closed 2mL, 10mL or 20mL vials PALsystem Vials with re-cappable vial caps are highly recommended for multiple opening and closing procedures within a workflow Requires FW Version 2.4 or higher | X | | |  |
| PAL3-DilutorMulti | Dilutor Module, consisting of: 1 pc Dilutor module 1 pc Connection cable 1 pc Dilutor Tool 1 pc Dilutor syringe 1mL 1 pc PFA Transfer tube 1/16" 5 pc Solvent bottle 1000ml 5 pc Aspirating tube 1/8" requires firmware version 2.1.0 or higher | X | X | |  |
| PAL3-DilutorSingle | Dilutor Module, consisting of: 1 pc Dilutor module 1 pc Connection cable 1 pc Dilutor Tool 1 pc Dilutor syringe 1mL 1 pc PFA Transfer tube 1/16" 1 pc Solvent bottle 1000ml 1 pc Aspirating tube 1/8" requires firmware version 2.1.0 or higher | X | X | |  |
| PAL3-FlowCell | Flow-Cell Module to accommodate up to 6 flow cells, consisting of: 1 pc Holder 88mm 1 pc Flow-Cell Stainless Steel 316 with connectors Swagelok 1/4" 5 pcs Inlet Septa | X | X | |  |

| Part No. | Description | Compatible with: | | | |
|----------------------|---|------------------|-----|-----|---|
| | | RTC | RSI | LSI | |
| PAL3-HeatexStirrer | Heatex Stirrer for intensive heating and stirring of 1x 20mL vial without stir bar Specific Heatex Stirrer insert for 10mL vial available (not included) requires firmware version 2.3 or higher (use in combination with existing PAL3-Agitator) | X | X | |  |
| PAL3-MHE-Module | Multiple Headspace Extraction Module, consisting of MHE module, MHE tool and MHE needle Gauge 22/5 | X | X | |  |
| PAL3-SPME-ArrowCond | SPME Conditioning Module for the conditioning of SPME Arrows and SPME fibers Includes a second position for conditioning a replacement SPME arrow or SPME fiber requires firmware version 2.3 or higher | X | X | |  |
| PAL3-ValveDrive-R | Valve Drive Module for Rheodyne valve (without valve) consisting of: 1 pc Valve Drive Module 1 pc PAL Bus connection cable requires firmware version 1.4.0 or higher | X | X | X |  |
| PAL3-ValveDrive-V | Valve Drive Module for Valco valve (without valve) consisting of: 1 pc Valve Drive module 1 pc PAL Bus connection cable requires firmware version 1.4.0 or higher | X | X | X |  |
| PAL3-VortexMixer | Vortex Mixer for the intensive mixing of 1x 2mL, 10 or 20mL vials with an agitation speed between 0-2000 rpm | X | X | X |  |
| PAL3-ParkStation-RTC | Park Station for the storage and robotic tool change of up to three PALtools like Liquid, HS, SPME, consisting of: 1 pc Park Station 1 pc Electrical connection cable | X | | |  |
| PAL3-SolventModule | Solvent Module for 3x 100mL solvent bottles | X | X | X |  |
| PAL3-WashFast | Fast Wash Module for two solvents (Needle Gauge 19 to 26), consisting of: 1 pc Wash Module with two wash positions and one Waste position 1 pc Pump Module with two solvent pumps 2 pc 1 liter solvent bottles with tubings and frits 1 pc PALbus connection cable requires firmware version 2.0.1 or higher | X | X | X |  |

| Part No. | Description | Compatible with: | | | |
|------------------|--|------------------|-----|-----|---|
| | | RTC | RSI | LSI | |
| PAL3-WashFast-HF | HF Fast Wash Module for two solvents (Needle Gauge 19 to 26), consisting of: 1 pc Wash Module with two wash positions and one Waste position 1 pc Pump Module with two solvent pumps 2 pc 1 liter solvent bottles with tubings and frits 1 pc PALbus connection cable Requires firmware version 2.4.x or higher | X | X | X |  |
| PAL3-WashLarge | Large Wash Module for 2x 100mL solvent vials and one waste position | X | X | X |  |
| PAL3-WashLCMS | Wash Module for LCMS Tool. with two wash positions and one waste position Requires FW 2.0.x or higher | X | X | |  |
| PAL3-WashStd | Standard Wash Module for up to 4x 10mL solvent vials and 1x 10mL waste vial | X | X | X |  |
| PAL3-Terminal | Terminal for PAL3 Autosampler Systems | X | X | X |  |



Tray Holders, Stacks and Racks

| Part No. | Description | Compatible with: | | | |
|-----------------------|--|------------------|-----|-----|---|
| | | RTC | RSI | LSI | |
| PAL3-LiqCooler | Tray Holder and sample tray for housing up to 32 vials of 10/20ml in the temperature range between 4 °C- 70°C External circulator bath for heating/cooling not included | X | X | X |  |
| PAL3-Trayholder | Tray Holder for up to 3 different MTP, DW plates or sample racks VT12, VT15, VT54, VT70 or 1 sample rack R60 | X | X | X |  |
| PAL3-Trayplate | Tray Plate for 3 microtiter / deepwell plates or racks (VT15, VT54 or VT70) | X | X | X |  |
| PAL3-PeltierStack12MT | Peltier Stack 12MT 4 °C- 40 °C for 12pcs Standard Microtiter Plates Including Power Supply and cables requires firmware version 1.5.0 or higher | X | X | X |  |
| PAL3-PeltierStack2DW | Peltier Stack 2 DW 4°C- 40°C for 2pcs Standard Microtiter-, Deepwell Plates, VT15, VT54 or VT70 racks Including Power Supply and cables requires firmware version 1.4.0 or higher | X | X | X |  |
| PAL3-PeltierStack6DW | Peltier Stack 6DW 4 °C- 40 °C for 6pcs Standard Microtiter-, Deepwell Plates, VT15, VT54 or VT70 racks Including Power Supply and cables requires firmware version 1.5.0 or higher | X | X | X |  |
| PAL3-Stack12MT | Stack 12MT for 6pcs Standard Microtiter plates requires firmware version 2.2.0 or higher | X | X | X |  |
| PAL3-Stack6DW | Stack 6DW for 6pcs Standard Microtiter-, Deepwell Plates, VT15, VT54 or VT70 racks requires firmware version 2.2.0 or higher | X | X | X |  |
| PAL3-Kit-VT15 | Sample Rack for 15 vials 10/20mL. Vials are not included, Pk of 3 pcs | X | X | X |  |
| PAL3-Kit-VT54 | Sample Rack for 54 vials 2mL. Vials are not included, Pk of 3 pcs | X | X | X | |
| PAL3-Rack-R60 | Sample Rack for 60 vials of 10/20mL, Aluminum Vials not included | X | X | X | |
| PAL3-Rack-VT12 | Sample Rack for 12 vials of 40mL. Vials are not included | X | X | X | |
| PAL3-Rack-VT15 | Sample Rack for 15 vials of 10/20mL. Vials are not included | X | X | X | |
| PAL3-Rack-VT54 | Sample Rack for 54 vials of 2mL. Vials are not included | X | X | X | |
| PAL3-Rack-VT70 | Sample Rack for 70 vials of 1mL. Vials are not included | X | X | X | |



PAL3 ITEX

ITEX Kit

| Part No. | Description | Compatible with: | |
|----------------|--|------------------|-----|
| | | RTC | RSI |
| PAL3-SITEX-Kit | PAL3 ITEX Tool, including: 1 pc Smart ITEX Syringe 1.3mL (SITEX-1300-T) 2 pc ITEX Trap TENAX TA 80/100 mesh (PAL3-ITEX-Trap-TXTA100) Compatible with all PAL RTC/PAL RSI models from FW version 2.3. or higher Requires firmware version 3.0 or higher and PAL RTC/PAL RSI Series II for Smart functionality | X | X |

ITEX Traps

| Part No. | Description | Tool compatibility |
|-------------------------|--|--------------------|
| PAL3-ITEX-Trap-Carb-S3 | PAL3 ITEX Trap Carbosieve SIII -60/ 80 mesh | ITEX |
| PAL3-ITEX-Trap-Carb1000 | PAL3 ITEX Trap Carboxen 1000 -60/ 80 mesh | ITEX |
| PAL3-ITEX-Trap-CarbC | PAL3 ITEX Trap Carbopack - C 80/100 mesh | ITEX |
| PAL3-ITEX-Trap-Mix-1 | PAL3 ITEX Trap Mix1 Tenax TA : Carboxen1000 (Bottom: Top - 1: 1) | ITEX |
| PAL3-ITEX-Trap-Mix-2 | PAL3 ITEX Trap Mix2 Carbopack C : Carbosieve SIII (Bottom: Top - 2: 1) | ITEX |
| PAL3-ITEX-Trap-Mix-3 | PAL3 ITEX Trap Mix3 Tenax GR : Carbosieve SIII (Bottom: Top - 1: 1) | ITEX |
| PAL3-ITEX-Trap-TxGR | PAL3 ITEX Trap Tenax GR - 80/100 mesh | ITEX |
| PAL3-ITEX-Trap-TXTA100 | PAL3 ITEX Trap Tenax TA 80/100 mesh | ITEX |

ITEX Syringe

| Part No. | Description | Tool compatibility |
|--------------|--|--------------------|
| SITEX-1300-T | Smart ITEX Syringe 1300µl for Tool ITEX, ITEX Trap not included. | Itex |

Smart Syringe & Needle Portfolio



Liquid Smart Syringes - Fixed Needle

| Part No. | Description | Tool compatibility |
|--------------------|--|--------------------|
| SF0-57-B-26P-CO | 0.5µL Smart Syringe with fixed needle for Tool D7/57: Needle length 57mm, gauge 26P, scale length 27mm, Point Style conical | D7/57 |
| SF1-57-B-26P-CO | 1.2µL Smart Syringe with fixed needle for Tool D7/57: Needle length 57mm, gauge 26P, scale length 36mm, Point Style conical | D7/57 |
| SF1-85-B-26P-CO | 1.2µL Smart Syringe with fixed needle for Tool D7/85: Needle length 85mm, gauge 26P, scale length 36mm, Point Style conical | D7/85 |
| SF10-57-M-23S-CO | 10uL Smart Syringe with fixed needle for Tool D7/57: Needle length 57mm, metal plunger, gauge 23S, scale length 54mm, Point Style conical | D7/57 |
| SF10-57-M-26S-CO | 10uL Smart Syringe with fixed needle for Tool D7/57: Needle length 57mm, metal plunger, gauge 26S, scale length 54mm, Point Style conical | D7/57 |
| SF10-57-T-22S-FL | 10uL Smart Syringe with fixed needle for Tool D7/57: Needle length 57mm, PTFE plunger, gauge 22S, scale length 60mm, Point Style flat | D7/57 |
| SF10-57-T-23S-CO | 10uL Smart Syringe with fixed needle for Tool D7/57: Needle length 57mm, PTFE plunger, gauge 23S, scale length 60mm, Point Style conical | D7/57 |
| SF10-57-T-26S-CO | 10uL Smart Syringe with fixed needle for Tool D7/57: Needle length 57mm, PTFE plunger, gauge 26S, scale length 60mm, Point Style conical | D7/57 |
| SF10-85-M-26S-CO | 10uL Smart Syringe with fixed needle for Tool D7/85: Needle length 85mm, metal plunger, gauge 26S, scale length 54mm, Point Style conical | D7/85 |
| SF10-85-T-26S-CO | 10uL Smart Syringe with fixed needle for Tool D7/85: Needle length 85mm, PTFE plunger, gauge 26S, scale length 60mm, Point Style conical | D7/85 |
| SF10-85-M-23S-CO | 10uL Smart Syringe with fixed needle for Tool D7/85: Needle length 85mm, metal plunger, gauge 23S, scale length 54mm, needle tip conical | D7/85 |
| SF25-57-T-26S-CO | 25uL Smart Syringe with fixed needle for Tool D7/57: Needle length 57mm, PTFE plunger, gauge 26S, scale length 60mm, Point Style conical | D7/57 |
| SF25-85-T-26S-CO | 25uL Smart Syringe with fixed needle for Tool D7/85: Needle length 85mm, PTFE plunger, gauge 26S, scale length 60mm, Point Style conical | D7/85 |
| SFX25-57-T-22S-FL | 25uL Smart Syringe with fixed needle for Tool D7/57: Needle length 57mm, PTFE plunger, gauge 22S, scale length 60mm, Point Style flat, x-Type | D7/57 |
| SF50-57-T-26S-CO | 50uL Smart Syringe with fixed needle for Tool D7/57: Needle length 57mm, PTFE plunger, gauge 26S, scale length 60mm, Point Style conical | D7/57 |
| SF50-85-T-26S-CO | 50uL Smart Syringe with fixed needle for Tool D7/85: Needle length 85mm, PTFE plunger, gauge 26S, scale length 60mm, Point Style conical | D7/85 |
| SF100-57-T-26S-CO | 100uL Smart Syringe with fixed needle for Tool D7/57: Needle length 57mm, PTFE plunger, gauge 26S, scale length 60mm, Point Style conical | D7/57 |
| SF100-57-T-23S-CO | 100uL Smart Syringe with fixed needle for Tool D7/57: Needle length 57mm, PTFE plunger, gauge 23S, scale length 60mm, Point Style conical | D7/57 |
| SF100-85-T-26S-CO | 100uL Smart Syringe with fixed needle for Tool D7/85: Needle length 85mm, PTFE plunger, gauge 26S, scale length 60mm, Point Style conical | D7/85 |
| SFX100-57-T-22-FL | 100uL Smart Syringe with fixed needle for Tool D7/57: Needle length 57mm, PTFE plunger, gauge 22, scale length 60mm, Point Style flat, x-Type | D7/57 |
| SFX100-57-T-22S-FL | 100uL Smart Syringe with fixed needle for Tool D7/57: Needle length 57mm, PTFE plunger, gauge 22S, scale length 60mm, Point Style flat, x-Type | D8/57 |
| SF250-57-T-26-CO | 250uL Smart Syringe with fixed needle for Tool D8/57: Needle length 57mm, PTFE plunger, gauge 26, scale length 60mm, Point Style conical | D8/57 |
| SF250-85-T-26-CO | 250uL Smart Syringe with fixed needle for Tool D8/85: Needle length 85mm, PTFE plunger, gauge 26, scale length 60mm, Point Style conical | D8/85 |
| SF250-57-T-22-FL | 250uL Smart Syringe with fixed needle for Tool D8/57: Needle length 57mm, PTFE plunger, gauge 22, scale length 60mm, Point Style flat | D8/57 |
| SF500-57-T-22-FL | 500uL Smart Syringe with fixed needle for Tool D8/57: Needle length 57mm, PTFE plunger, gauge 22, scale length 60mm, Point Style flat | D8/57 |

CAUTION: SYRINGES FOR LABORATORY USE ONLY - NOT FOR HUMAN IN VIVO USE

Liquid Smart Syringes - Fixed Needle

| Part No. | Description | Tool compatibility |
|--------------------|---|--------------------|
| SF500-57-T-26-CO | 500uL Smart Syringe with fixed needle for Tool D8/57: Needle length 57mm, PTFE plunger, gauge 26, scale length 60mm, Point Style conical | D8/57 |
| SF1000-57-T-22-FL | 1000uL Smart Syringe with fixed needle for Tool D8/57: Needle length 57mm, PTFE plunger, gauge 22, scale length 60mm, Point Style flat | D8/57 |
| SF1000-57-T-23-CO | 1000uL Smart Syringe with fixed needle for Tool D8/57: Needle length 57mm, PTFE plunger, gauge 23, scale length 60mm, Point Style conical | D8/57 |
| SC2500-57-PE-22-FL | 2500uL Smart Syringe with fixed needle for Tool D18/57: Needle length 57mm, Polyethylene plunger, gauge 22, scale length 60mm, Point Style flat | D18/57 |
| SC2500-57-T-19-FL | 2500uL Smart Syringe with fixed needle for Tool D18/57: Needle length 57mm, PTFE plunger, gauge 19, scale length 60mm, Point Style flat | D18/57 |
| SC2500-57-T-22-FL | 2500uL Smart Syringe with fixed needle for Tool D18/57: Needle length 57mm, PTFE plunger, gauge 22, scale length 60mm, Point Style flat | D18/57 |
| SC5000-57-T-19-FL | 5000uL Smart Syringe with fixed needle for Tool D18/57: Needle length 57mm, PTFE plunger, gauge 19, scale length 60mm, Point Style flat | D18/57 |
| SC5000-57-T-22-FL | 5000uL Smart Syringe with fixed needle for Tool D18/57: Needle length 57mm, PTFE plunger, gauge 22, scale length 60mm, Point Style flat | D18/57 |

LCMS Smart Syringes - without Needle

| Part No. | Description | Tool compatibility |
|--------------|--|--------------------|
| SLCMS-100-PE | 100µl Smart LCMS Syringe for LCMS Tool, PE Plunger, Gastight, Scale Length 60 mm, Thread 1/4-28 UNF Removable Needle not included | LCMS |
| SLCMS-100-T | 100µl Smart LCMS Syringe for LCMS Tool, PTFE Plunger, Gastight, Scale Length 60 mm, Thread 1/4-28 UNF Removable Needle not included | LCMS |
| SLCMS-250-PE | 250µl Smart LCMS Syringe for LCMS Tool, PE Plunger, Gastight, Scale Length 60 mm, Thread 1/4-28 UNF Removable Needle not included | LCMS |
| SLCMS-250-T | 250µl Smart LCMS Syringe for LCMS Tool, PTFE Plunger, Gastight, Scale Length 60 mm, Thread 1/4-28 UNF Removable Needle not included | LCMS |

ITEX Smart Syringes - without Needle

| Part No. | Description | Tool compatibility |
|-------------|--|--------------------|
| SITEX1300-T | Smart ITEX Syringe 1300µl for Tool ITEX, ITEX Trap not included. | ITEX |

Liquid Smart Syringes - Exchangeable Needle

| Part No. | Description | Tool compatibility |
|-----------------|--|--------------------|
| SE0-B-57-23P-CO | 0.5µL Smart Syringe body for exchangeable needles for Tool D8/57, Plunger in Needle, Needle length 57mm, gauge 23P, scale length 27mm, Point Style conical, Needle Flange Type NIP_0 | D8/57 |
| SE1-B-57-23P-CO | 1µL Smart Syringe body for exchangeable needles for Tool D8/57, Plunger in Needle, Needle length 57mm, gauge 23P, scale length 27mm, Point Style conical, Needle Flange Type NIP_1 | D8/57 |

Liquid Smart Syringes - Exchangeable Needle

CAUTION: SYRINGES FOR LABORATORY USE ONLY - NOT FOR HUMAN IN VIVO USE

| Part No. | Description | Tool compatibility |
|-----------|---|--------------------|
| SE5-M | 5µL Smart Syringe body for exchangeable needles for Tool D8, Metal Plunger, Needle Flange Type F_5, Needle not included | D8/57 |
| SE10-M | 10µL Smart Syringe body for exchangeable needles for Tool D8, Metal Plunger, Needle Flange Type F_10, Needle not included | D8/57 |
| SE10-T | 10µL Smart Syringe body for exchangeable needles for Tool D8, PTFE Plunger, Needle Flange Type F_10, Needle not included | D8/57 |
| SE25-T | 25µL Smart Syringe body for exchangeable needles for Tool D8, PTFE Plunger, Needle Flange Type F_25, Needle not included | D8/57 |
| SE50-T | 50µL Smart Syringe body for exchangeable needles for Tool D8, PTFE Plunger, Needle Flange Type F_25, Needle not included | D8/57 |
| SE100-T | 100µL Smart Syringe body for exchangeable needles for Tool D8, PTFE Plunger, Needle Flange Type F_25, Needle not included | D8/57 |
| SE250-T | 250µL Smart Syringe body for exchangeable needles for Tool D8, PTFE Plunger, Needle Flange Type F_25, Needle not included | D8/57 |
| SE500-T | 500µL Smart Syringe body for exchangeable needles for Tool D8, PTFE Plunger, Needle Flange Type F_25, Needle not included | D8/57 |
| SE1000-T | 1000µL Smart Syringe body for exchangeable needles for Tool D8, PTFE Plunger, Needle Flange Type F_1000, Needle not included | D8/57 |
| SCE2500-T | 2500µL Smart Syringe body for exchangeable needles for Tool D18, PTFE Plunger, Needle Flange Type F_1000, Needle not included | D18/57 |
| SCE5000-T | 5000µL Smart Syringe body for exchangeable needles for Tool D18, PTFE Plunger, Needle Flange Type F_1000, Needle not included | D18/57 |
| SE10000-T | 10000µL Smart Syringe body for exchangeable needles for Tool D18, PTFE Plunger, Needle Flange Type F_10000, Needle not included | D18/57 |

Headspace Smart Syringes - Fixed Needle - Needle Fixation with Glue

| Part No. | Description | Tool compatibility |
|-------------------|--|--------------------|
| SH1000-65-T-23-SP | 1000µL Smart Headspace Syringe with fixed needle for Tool HS1000, PTFE Plunger, Needle length 65mm, Gauge 23, Point Style Sideport | HS1000 |
| SH2500-65-T-23-SP | 2500µL Smart Headspace Syringe with fixed needle for Tool HS2500, PTFE Plunger, Needle length 65mm, Gauge 23, Point Style Sideport | HS2500 |
| SH5000-65-T-23-SP | 5000µL Smart Headspace Syringe with fixed needle for Tool HS5000, PTFE Plunger, Needle length 65mm, Gauge 23, Point Style Sideport | HS5000 |

Headspace Smart Syringes - Fixed Needle - Glue Free Needle Fixation

| Part No. | Description | Tool compatibility |
|------------------|---|--------------------|
| SHGF2250-T-23-SP | 2250µl Smart Glue Free Headspace Syringe with fixed needle for Tool HS2500, PTFE plunger, needle gauge 23, Point Style Sideport | HS2500 |
| SHGF4500-T-23-SP | 4500µl Smart Glue Free Headspace Syringe with fixed needle for Tool HS5000, PTFE plunger, needle gauge 23, Point Style Sideport | HS5000 |

Headspace Smart Syringes - Exchangeable Needle

| Part No. | Description | Tool compatibility |
|-----------|--|--------------------|
| SHE1000-T | 1000µL Smart Headspace Syringe for exchangeable needles for Tool HS1000, PTFE Plunger, Needle Flange Type HS_1000 - Needle not included! | HS1000 |
| SHE2500-T | 2500µL Smart Headspace Syringe for exchangeable needles for Tool HS2500, PTFE Plunger, Needle Flange Type HS_1000 - Needle not included! | HS2500 |
| SHE5000-T | 5000µL Smart Headspace Syringe for exchangeable needles for Tool HS5000, PTFE Plunger, Needle Flange Type HS_1000 - Needle not included! | HS5000 |

CAUTION: SYRINGES FOR LABORATORY USE ONLY - NOT FOR HUMAN IN VIVO USE

Liquid Smart Syringe Needles - Needle Flange Type F_5

| Part No. | Description | Syringe compatibility |
|---------------|---|-----------------------|
| ND5-57-23S-CO | Exchangeable needle for liquid syringes with Needle Flange type ND5 Body Type 5µL Liquid Syringe, Needle Length 57mm, Gauge 23S, Point Style conical | SE5-M |
| ND5-57-26S-CO | Exchangeable needle for liquid syringes with Needle Flange type ND5 Body Type 5µL Liquid Syringe, Needle Length 57mm, Gauge 26S, Point Style conical | SE5-M |
| ND5-85-23S-CO | Exchangeable needle for liquid syringes with Needle Flange type ND5 Body Type 5µL Liquid Syringe, Needle Length 85mm, Gauge 23S, Point Style conical | SE5-M |
| ND5-85-26S-CO | Exchangeable needle for liquid syringes with Needle Flange type ND5 Body Type 5µL Liquid Syringe, Needle Length 85mm, Gauge 26S, Point Style conical | SE5-M |

Liquid Smart Syringe Needles - Needle Flange Type F_10

| Part No. | Description | Syringe compatibility |
|----------------|---|-----------------------|
| ND10-57-23S-CO | Exchangeable needle for liquid syringes with Needle Flange type ND10 Body Type 10µL Liquid Syringe, Needle Length 57mm, Gauge 23S, Point Style conical | SE10-M, SE10-T |
| ND10-57-26S-CO | Exchangeable needle for liquid syringes with Needle Flange type ND10 Body Type 10µL Liquid Syringe, Needle Length 57mm, Gauge 26S, Point Style conical | SE10-M, SE10-T |
| ND10-85-23S-CO | Exchangeable needle for liquid syringes with Needle Flange type ND10 Body Type 10µL Liquid Syringe, Needle Length 85mm, Gauge 23S, Point Style conical | SE10-M, SE10-T |
| ND10-85-26S-CO | Exchangeable needle for liquid syringes with Needle Flange type ND10 Body Type 10µL Liquid Syringe, Needle Length 85mm, Gauge 26S, Point Style conical | SE10-M, SE10-T |

Liquid Smart Syringe Needles - Needle Flange Type F_25

| Part No. | Description | Syringe compatibility |
|----------------|--|--|
| ND25-57-22S-FL | Exchangeable needle for liquid syringes with Needle Flange type ND25 Body Type 25, 50, 100, 250 or 500µL Liquid Syringe, Needle Length 57mm, Gauge 22S, Point Style flat | SE25-T, SE50-T, SE100-T. SE250-T, SE500-T |
| ND25-57-23-CO | Exchangeable needle for liquid syringes with Needle Flange type ND25 Body Type 25, 50, 100, 250 or 500µL Liquid Syringe, Needle Length 57mm, Gauge 23, Point Style conical | SE25-T, SE50-T, SE100-T. SE250-T, SE500-T |
| ND25-57-23S-CO | Exchangeable needle for liquid syringes with Needle Flange type ND25 Body Type 25, 50, 100, 250 or 500µL Liquid Syringe, Needle Length 57mm, Gauge 23S, Point Style conical | SE25-T, SE50-T, SE100-T. SE250-T, SE500-T |
| ND25-57-26-CO | Exchangeable needle for liquid syringes with Needle Flange type ND25 Body Type 25, 50, 100, 250 or 500µL Liquid Syringe, Needle Length 57mm, Gauge 26, Point Style conical | SE25-T, SE50-T, SE100-T. SE250-T, SE500-T |
| ND25-57-26S-CO | Exchangeable needle for liquid syringes with Needle Flange type ND25 Body Type 25, 50, 100, 250 or 500µL Liquid Syringe, Needle Length 57mm, Gauge 26S, Point Style conical | SE25-T, SE50-T, SE100-T. SE250-T, SE500-T |
| ND25-85-22-FL | Exchangeable needle for liquid syringes with Needle Flange type ND25 Body Type 25, 50, 100, 250 or 500µL Liquid Syringe, Needle Length 85mm, Gauge 22, Point Style flat | SE25-T, SE50-T, SE100-T. SE250-T, SE500-T |
| ND25-85-22S-FL | Exchangeable needle for liquid syringes with Needle Flange type ND25 Body Type 25, 50, 100, 250 or 500µL Liquid Syringe, Needle Length 85mm, Gauge 22S, Point Style flat | SE25-T, SE50-T, SE100-T. SE250-T, SE500-T |

CAUTION: SYRINGES FOR LABORATORY USE ONLY - NOT FOR HUMAN IN VIVO USE

Liquid Smart Syringe Needles - Needle Flange Type F_25

| Part No. | Description | Syringe compatibility |
|----------------|--|--|
| ND25-85-23-CO | Exchangeable needle for liquid syringes with Needle Flange type ND25 Body Type 25, 50, 100, 250 or 500µL Liquid Syringe, Needle Length 85mm, Gauge 23, Point Style conical | SE25-T, SE50-T, SE100-T. SE250-T, SE500-T |
| ND25-85-23S-CO | Exchangeable needle for liquid syringes with Needle Flange type ND25 Body Type 25, 50, 100, 250 or 500µL Liquid Syringe, Needle Length 85mm, Gauge 23S, Point Style conical | SE25-T, SE50-T, SE100-T. SE250-T, SE500-T |
| ND25-85-26-CO | Exchangeable needle for liquid syringes with Needle Flange type ND25 Body Type 25, 50, 100, 250 or 500µL Liquid Syringe, Needle Length 85mm, Gauge 26, Point Style conical | SE25-T, SE50-T, SE100-T. SE250-T, SE500-T |
| ND25-85-26S-CO | Exchangeable needle for liquid syringes with Needle Flange type ND25 Body Type 25, 50, 100, 250 or 500µL Liquid Syringe, Needle Length 85mm, Gauge 26S, Point Style conical | SE25-T, SE50-T, SE100-T. SE250-T, SE500-T |

Liquid Smart Syringe Needles - Needle Flange Type F_1000

| Part No. | Description | Syringe compatibility |
|-----------------|--|-----------------------------------|
| ND1000-57-22-FL | Exchangeable needle for liquid syringes with Needle Flange type ND1000 Body Type 1000, 2500 or 5000µL Liquid Syringe, Needle Length 57mm, Gauge 22, Point Style flat. Caution: not for Headspace syringes | SE1000-T, SCE2500-T, SCE5000-T |
| ND1000-57-22-SP | Exchangeable needle for liquid syringes with Needle Flange type ND1000 Body Type 1000, 2500 or 5000µL Liquid Syringe, Needle Length 57mm, Gauge 22, Point Style side port - Caution: not for Headspace syringes | SE1000-T, SCE2500-T, SCE5000-T |
| ND1000-57-23-CO | Exchangeable needle for liquid syringes with Needle Flange type ND1000 Body Type 1000, 2500 or 5000µL Liquid Syringe, Needle Length 57mm, Gauge 23, Point Style conical - Caution: not for Headspace syringes | SE1000-T, SCE2500-T, SCE5000-T |
| ND1000-57-23-SP | Exchangeable needle for liquid syringes with Needle Flange type F_1000 Body Type 1000, 2500 or 5000µL Liquid Syringe, Needle Length 57mm, Gauge 23, Point Style side port. Caution: not for Headspace syringes | SE1000-T, SCE2500-T, SCE5000-T |
| ND1000-57-26-CO | Exchangeable needle for liquid syringes with Needle Flange type ND1000 Body Type 1000, 2500 or 5000µL Liquid Syringe, Needle Length 57mm, Gauge 26, Point Style conical. Caution: not for Headspace syringes | SE1000-T, SCE2500-T, SCE5000-T |
| ND1000-85-19-SP | Exchangeable needle for liquid syringes with Needle Flange type ND1000 Body Type 1000, 2500 or 5000µL Liquid Syringe, Needle Length 85mm, Gauge 19, Point Style side port. Caution: not for Headspace syringes | SE1000-T, SCE2500-T, SCE5000-T |
| ND1000-85-22-FL | Exchangeable needle for liquid syringes with Needle Flange type ND1000 Body Type 1000, 2500 or 5000µL Liquid Syringe, Needle Length 85mm, Gauge 22, Point Style flat. Caution: not for Headspace syringess | SE1000-T, SCE2500-T, SCE5000-T |
| ND1000-85-22-SP | Exchangeable needle for liquid syringes with Needle Flange type ND1000 Body Type 1000, 2500 or 5000µL Liquid Syringe, Needle Length 85mm, Gauge 22, Point Style side port. Caution: not for Headspace syringes | SE1000-T, SCE2500-T, SCE5000-T |
| ND1000-85-23-CO | Exchangeable needle for liquid syringes with Needle Flange type ND1000 Body Type 1000, 2500 or 5000µL Liquid Syringe, Needle Length 85mm, Gauge 23, Point Style conical - Caution: not for Headspace syringes | SE1000-T, SCE2500-T, SCE5000-T |
| ND1000-85-23-SP | Exchangeable needle for liquid syringes with Needle Flange type ND1000 Body Type 1000, 2500 or 5000µL Liquid Syringe, Needle Length 85mm, Gauge 23, Point Style side port. Caution: not for Headspace syringes | SE1000-T, SCE2500-T, SCE5000-T |
| ND1000-85-26-CO | Exchangeable needle for liquid syringes with Needle Flange type ND1000 Body Type 1000, 2500 or 5000µL Liquid Syringe, Needle Length 85mm, Gauge 26, Point Style conical - Caution: not for Headspace syringes | SE1000-T, SCE2500-T, SCE5000-T |

CAUTION: SYRINGES FOR LABORATORY USE ONLY - NOT FOR HUMAN IN VIVO USE

Liquid Smart Syringe Needles - Needle Flange Type F_ 10000

| Part No. | Description | Syringe compatibility |
|------------------|--|-----------------------|
| ND10000-57-19-SP | Exchangeable needle for liquid syringes with Needle Flange type ND10000 Body Type 10000µL Liquid Syringe, Needle Length 57mm, Gauge 19, Point Style side port | SE10000-T |
| ND10000-57-22-FL | Exchangeable needle for liquid syringes with Needle Flange type ND10000 Body Type 10000µL Liquid Syringe, Needle Length 57mm, Gauge 22, Point Style flat | SE10000-T |
| ND10000-57-22-SP | Exchangeable needle for liquid syringes with Needle Flange type ND10000 Body Type 10000µL Liquid Syringe, Needle Length 57mm, Gauge 22, Point Style side port | SE10000-T |
| ND10000-57-23-CO | Exchangeable needle for liquid syringes with Needle Flange type ND10000 Body Type 10000µL Liquid Syringe, Needle Length 57mm, Gauge 23, Point Style conical | SE10000-T |
| ND10000-57-23-SP | Exchangeable needle for liquid syringes with Needle Flange type ND10000 Body Type 10000µL Liquid Syringe, Needle Length 57mm, Gauge 23, Point Style side port | SE10000-T |

Headspace Smart Syringe Needles - Needle Flange Type HS_ 1000

| Part No. | Description | Syringe compatibility |
|------------------|---|------------------------------------|
| NDH1000-65-22-SP | Exchangeable needle for liquid syringes with Needle Flange type ND1000 Body Type 1000, 2500 or 5000µL Headspace Syringe, Needle Length 65mm, Gauge 22, Point Style side port | SHE1000-T, SHE2500-T, SHE5000-T |
| NDH1000-65-23-CO | Exchangeable needle for liquid syringes with Needle Flange type ND1000 Body Type 1000, 2500 or 5000µL Headspace Syringe, Needle Length 65mm, Gauge 23, Point Style conical | SHE1000-T, SHE2500-T, SHE5000-T |
| NDH1000-65-23-SP | Exchangeable needle for liquid syringes with Needle Flange type ND1000 Body Type 1000, 2500 or 5000µL Headspace Syringe, Needle Length 65mm, Gauge 23, Point Style side port | SHE1000-T, SHE2500-T, SHE5000-T |
| NDH1000-65-26-CO | Exchangeable needle for liquid syringes with Needle Flange type ND1000 Body Type 1000, 2500 or 5000µL Headspace Syringe, Needle Length 65mm, Gauge 26, Point Style conical | SHE1000-T, SHE2500-T, SHE5000-T |



Smart SPME Arrow / Smart SPME Fiber

PAL SPME Kit

| Part No. | Description | Compatible with: | |
|--------------------|---|------------------|-----|
| | | RTC | RSI |
| PAL3-SPME-SFib-Kit | SPME Kit, consisting of: 1 pc SPME Tool 1 pc SPME Fiber Assortment Kit (SFIB-SEL5-S2) 1 pc SPME Performance Evaluation Mix | X | X |

PAL SPME Fibers

| Part No. | Description | Tool compatibility |
|------------------------|---|--------------------|
| SFIB-A-85/10-P3 | Smart SPME fiber Acrylate (Polyacrylate) Fiber thickness 85µm, Fiber length 10mm, color code: gray, pack of 3pcs | SPME Fiber |
| SFIB-C-WR-95/10-P3 | Smart SPME fiber C-WR/PDMS (Carbon Wide Range / Polydimethylsiloxane) Fiber thickness 95µm, Fiber length 10mm, color code: dark blue, pack of 3pcs | SPME Fiber |
| SFIB-DVB-65/10-P3 | Smart SPME fiber DVB/PDMS (Divinylbenzene / Polydimethylsiloxane) Fiber thickness 65µm, Fiber length 10mm color code: violet, pack of 3pc | SPME Fiber |
| SFIB-DVB/C-WR-80/10-P3 | Smart SPME fiber DVB/C-WR/PDMS (Divinylbenzene / Carbon Wide Range / Polydimethylsiloxane) Fiber thickness 50/30µm, Fiber length 10mm color code: dark gray, pack of 3pc | SPME Fiber |
| SFIB-P-100/10-P3 | Smart SPME fiber PDMS (Polydimethylsiloxane) Fiber thickness 100µm, Fiber length 10mm, color code: red, pack of 3pcs | SPME Fiber |
| SFIB-P-30/10-P3 | Smart SPME fiber PDMS (Polydimethylsiloxane) Fiber thickness 30µm, Fiber length 10mm, color code: golden, pack of 3pcs | SPME Fiber |
| SFIB-P-7/10-P3 | Smart SPME fiber PDMS (Polydimethylsiloxane) Fiber thickness 7µm, Fiber length 10mm, color code: green, pack of 3pcs | SPME Fiber |
| SFIB-SEL5-S1 | Collection of five different Smart SPME fibers with fiber length of 10mm and fiber thickness of: PDMS: 7µm, 30µm & 100µm Acrylate: 85µm C-WR: 95µm | SPME Fiber |
| SFIB-SEL5-S2 | Collection of five different Smart SPME fibers with fiber length of 10mm and fiber thickness of: PDMS: 100µm DVB: 65µm C-WR: 95µm Acrylate: 85µm DVB/C-WR/PDMS: 80µm | SPME Fiber |

All PAL Smart SPME Fibers are available as single fiber, pack of 3 fibers and pack of 5 fibers

PAL SPME Arrows

| Part No. | Description | Tool compatibility |
|-------------------------|---|--------------------|
| SARR11-A-100/20-P3 | Smart SPME Arrow 1.10mm: Acrylate (Polyacrylate) Phase thickness 100µm, Phase length 20mm, color code: grey, pack of 3 pcs | SPME Arrow |
| SARR11-C-WR-120/20-P3 | Smart SPME Arrow 1.10mm: C-WR/PDMS (Carbon Wide Range) Phase thickness 120µm, Phase length 20mm, color code: light blue, pack of 3 pcs | SPME Arrow |
| SARR11-DVB-120/20-P3 | Smart SPME Arrow 1.10mm: DVB/PDMS (Divinylbenzene) Phase thickness 120µm, Phase length 20mm, color code: violett, pack of 3 pcs | SPME Arrow |
| SARR11-DVB/CWR120/20-P3 | Smart SPME Arrow 1.10mm: DVB/C-WR/PDMS (Divinylbenzene / Carbon Wide Range / Polydimethylsiloxane), Phase thickness 110µm, Phase length 20mm, color code: dark gray, pack of 3 pc | SPME Arrow |
| SARR11-P-100/20-P3 | Smart SPME Arrow 1.10mm: PDMS (Polydimethylsiloxane) Phase thickness 100µm, Phase length 20mm, color code: red, pack of 3 pcs | SPME Arrow |
| SARR15-C-WR-120/20-P3 | Smart SPME Arrow 1.50mm wide sleeve: C-WR/PDMS (Carbon Wide Range / Polydimethylsiloxane) Phase thickness 120µm, Phase length 20mm, color code: light blue, pack of 3 pcs | SPME Arrow |
| SARR15-DVB-120/20-P3 | Smart SPME Arrow 1.50mm wide sleeve: DVB/PDMS (Divinylbenzene / Polydimethylsiloxane) Phase thickness 120µm, Phase length 20mm, color code: violett, pack of 3 pcs | SPME Arrow |
| SARR15-DVB/CWR120/20-P3 | Smart SPME Arrow 1.50mm wide sleeve: DVB/C-WR/PDMS (Divinylbenzene/Carbon Wide Range/Polydimethylsiloxane) Phase thickness 110µm, Phase length 20mm, color code: dark gray, pack of 3 pc | SPME Arrow |
| SARR15-P-100/20-P3 | Smart SPME Arrow 1.50mm wide sleeve: PDMS (Polydimethylsiloxane) Phase thickness 100µm, Phase length 20mm, color code: red, pack of 3 pcs | SPME Arrow |
| SARR15-P-250/20-P3 | Smart SPME Arrow 1.50mm: PDMS (Polydimethylsiloxane) Phase thickness 250µm, Phase length 20mm, color code: black, pack of 3 pc | SPME Arrow |
| SARR1115-SEL5-S2 | Collection of five different Smart SPME Arrows with phase length of 20mm: 1 pc Smart SPME Arrow 1.10mm Acrylate (Polyacrylate), Phase thickness 100µm 1 pc Smart SPME Arrow 1.10mm C-WR/PDMS (Carbon Wide Range / Polydimethyl- siloxane), Phase thickness 120µm 1 pc Smart SPME Arrow 1.10mm DVB/PDMS (Divinylbenzene / Polydimethylsiloxan), Phase thickness 120µm 1 pc Smart SPME Arrow 1.10mm PDMS (Polydimethylsiloxane) , Phase thickness 100µm 1 pc Smart SPME Arrow 1.10mm DVB/Carbon WR/PDMS (Divinylbenzene / Carbon wide range / Polydimethylsiloxane) , Phase thickness 110µm | SPME Arrow |

All PAL Smart SPME Arrows are available as single Arrow, pack of 3 Arrows and pack of 5 Arrows

PAL SPME Arrow Accessories

| Part No. | Description | Compatible with: | | |
|--------------------------|---|------------------|-----|-----|
| | | RTC | RSI | LSI |
| PAL3-AgiHeatex-Kit | Agitator and Heatex Stirrer Kit for SPME Arrow consisting of: 1 Agitator for the incubation and agitation of up to 6x 20mL vials 1 Heatex Stirrer for intensive heating and stirring of 1x 20mL vial without stir bar Specific Heatex Stirrer insert for 10mL vial available (not included) Specific Agitator insert for 10 mL vials available (not included) requires firmware version 2.3 or higher and adapters for 10mL vials | X | X | |
| PAL3-Heatex-Ins-10mL | Specific insert for the Heatex Stirrer to stir one 10mL vial | X | X | |
| PAL3-Heatex-Ins-2mL | Specific insert for the Heatex Stirrer to stir one 2mL vial | X | X | |
| ARR-Man-Inj-Kit | Manual Injection Kit for SPME Arrow, consisting of: 1x Manual SPME Holder 1x Manual SPME Extraction Guide 1x Manual SPME Injection Guide | | | |
| PAL3-SARR-Start-6890-LCh | Starter Kit Smart SPME Arrow for Agilent GC 6890 with large charcoal filter consisting of: 1 Adaptation for the split/splitless injector with large charcoal filter of AGILENT GC 6890 (ARR-SSL-Inj-GC6890-LCh) 2 SPME Arrow Liner for SSL Injector of Agilent GC 6890 1x PAL Smart SPME Arrow Kit (PAL3-SPME-SArr-Kit) 1 pc SPME Arrow Tool 1 pc Smart SPME Arrow Assortment Kit (SARR1115-SEL5-S2) 1 pc SPME Performance Evaluation Mix | X | X | |
| PAL3-SARR-Start-6890-STD | Starter Kit Smart SPME Arrow for Agilent GC 6890 with standard Assy consisting of: 1 Adaptation for the split/splitless injector with standard Assy of AGILENT GC 6890 (ARR-SSL-Inj-GC6890-STD) 2 SPME Arrow Liner for SSL Injector of Agilent GC 6890 1x PAL Smart SPME Arrow Kit (PAL3-SPME-SArr-Kit) 1 pc SPME Arrow Tool 1 pc Smart SPME Arrow Assortment Kit (SARR1115-SEL5-S2) 1 pc SPME Performance Evaluation Mix | X | X | |
| PAL3-SARR-Start-GC2010 | Starter Kit Smart SPME Arrow for Shimadzu GC-2010 Plus consisting of: 1x Adaptation for the split/splitless injector of SHIMADZU GC-2010 Plus (ARR-SSL-Inj-GC2010) 1 pc Liner Nut 1 pc Screw Cap 2 pcs SPME Arrow Liner 1x PAL Smart SPME Arrow Kit (PAL3-SPME-SArr-Kit) 1 pc SPME Arrow Tool 1 pc Smart SPME Arrow Assortment Kit (SARR1115-SEL5-S2) 1 pc SPME Performance Evaluation Mix (only for SHIMADZU GC-2010 Plus, not suitable for use with GC-2010) | X | X | |
| PAL3-SARR-Start-GC2030 | Starter Kit Smart SPME Arrow for Agilent GC 6890 with standard Assy consisting of: 1 Adaptation for the split/splitless injector with standard Assy of AGILENT GC 6890 (ARR-SSL-Inj-GC6890-STD) 2 SPME Arrow Liner for SSL Injector of Agilent GC 6890 1x PAL Smart SPME Arrow Kit (PAL3-SPME-SArr-Kit) 1 pc SPME Arrow Tool 1 pc Smart SPME Arrow Assortment Kit (SARR1115-SEL5-S2) 1 pc SPME Performance Evaluation Mix | X | X | |

All Starter Kits are compatible with all PAL RTC/PAL RSI models from FW version 2.3. or higher, Requiring firmware version 3.0 or higher and PAL RTC/PAL RSI Series II for Smart functionality

| Part No. | Description | Compatible with: | | |
|-------------------------|---|------------------|-----|-----|
| | | RTC | RSI | LSI |
| PAL3-SARR-Start-GC7890 | Starter Kit Smart SPME Arrow for Agilent GC 7890 consisting of: 1 Adaptation for the split/splitless injector of AGILENT GC 7890 (ARR-SSL-Inj-GC7890) 2 SPME Arrow Liner for SSL Injector of Agilent GC 7890 1x PAL Smart SPME Arrow Kit (PAL3-SPME-SArr-Kit) 1 pc SPME Arrow Tool 1 pc Smart SPME Arrow Assortment Kit (SARR1115-SEL5-S2) 1 pc SPME Performance Evaluation Mix | X | X | |
| PAL3-SARR-Start-GC8890 | Starter Kit Smart SPME Arrow for Agilent GC 8890 / 8860 consisting of: 1 Adaptation for the split/splitless injector of AGILENT GC (ARR-SSL-Inj-GC8890) 2 SPME Arrow Liner for SSL Injector of Agilent GC 1x PAL Smart SPME Arrow kit (PAL3-SPME-SArr-Kit) 1 pc SPME Arrow Tool 1 pc Smart SPME Arrow Assortment Kit (SARR1115-SEL5-S2) 1 pc SPME Performance Evaluation Mix | X | X | |
| PAL3-SARR-Start-Tr1300 | Starter Kit Smart SPME Arrow for Thermo GC Trace1300/1310 consisting of: 1 Adaptation for the split/splitless injector of Thermo GC Trace1300/1310 (ARR-SSL-Inj-Trace1300) 1 SPME Arrow Liner for SSL Injector of Thermo GC Trace1300/1310 Pk of 2 pcs 1x PAL Smart SPME Arrow Kit (PAL3-SPME-SArr-Kit) 1 pc SPME Arrow Tool 1 pc Smart SPME Arrow Assortment Kit (SARR1115-SEL5-S2) 1 pc SPME Performance Evaluation Mix | X | X | |
| PAL3-SARR-Start-TrUltra | Starter Kit Smart SPME Arrow for Thermo GC TraceUltra consisting of: 1 Adaptation for the split/splitless injector of Thermo GC TraceUltra (ARR-SSL-Inj-TraceUltra) 2 SPME Arrow Liner for SSL Injector of Thermo GC TraceUltra 1x PAL Smart SPME Arrow Kit (PAL3-SPME-SArr-Kit) 1 pc SPME Arrow Tool 1 pc Smart SPME Arrow Assortment Kit (SARR1115-SEL5-S2) 1 pc SPME Performance Evaluation Mix | X | X | |
| PAL3-SPME-SArr-Kit | Smart SPME Arrow Kit without GC Adaption Kit consisting of: 1 pc SPME Arrow Tool 1 pc Smart SPME Arrow Assortment Kit (SARR1115-SEL5-S2) 1 pc SPME Performance Evaluation Mix (Requires adaption kit for GC-specific injector) | X | X | |

All Starter Kits are compatible with all PAL RTC/PAL RSI models from FW version 2.3. or higher, Requiring firmware version 3.0 or higher and PAL RTC/PAL RSI Series II for Smart functionality

PAL SPME Arrow Liner

| Part No. | Description |
|-----------------------|--|
| ARR-Liner-CondModule | Liner for SPME Arrow Conditioning Module, Pk of 3 pcs |
| ARRLIN13-GC2010-3 | SPME Arrow Liner ID 1.3mm, for SSL Injector of SHIMADZU GC 2010, Pk of 3 pcs Only for SMPE Arrow OD 1.1mm |
| ARRLIN13-GC6890-3 | SPME Arrow Liner ID 1.3mm, for SSL Injector of AGILENT GC 6890, Pk of 3 pcs Only for SPME Arrow with OD 1.1mm |
| ARRLIN13-GC7890-3 | SPME Arrow Liner ID 1.3mm, for SSL Injector of AGILENT GC 7890, Pk of 3 pcs Only for SMPE Arrow with OD 1.1mm |
| ARRLIN13-Trace1300-3 | SPME Arrow Liner ID 1.3mm, for SSL Injector of Thermo GC Trace1300, Pk of 3 pcs Only for SPME Arrow with OD 1.1mm |
| ARRLIN17-GC2010-3 | SPME Arrow Liner ID 1.7mm, for SSL Injector of SHIMADZU GC 2010, Pk of 3 pcs For any SPME Arrow |
| ARRLIN17-GC6890-3 | SPME Arrow Liner ID 1.7mm, for SSL Injector of AGILENT GC 6890, Pk of 3 pcs For any SPME Arrow |
| ARRLIN17-GC7890-3 | SPME Arrow Liner ID 1.7mm, for SSL Injector of AGILENT GC 7890, Pk of 3 pcs For any SPME Arrow |
| ARRLIN17-Trace1300-3 | SPME Arrow Liner ID 1.7mm, for SSL Injector of Thermo GC Trace1300, Pk of 3 pcs For any SPME Arrow |
| ARRLIN20-TraceUltra-3 | SPME Arrow Liner ID 2.0mm, for SSL Injector of Thermo GC TraceUltra, Pk of 3 pcs For any SPME Arrow |

All PAL SPME Arrow Liners are available as single Liner and as pack of 3 Liners.



μ SPE Options

μSPE Smart Option

| Part No. | Description | Compatible with | | |
|-----------------|--|-----------------|-----|-----|
| | | RTC | RSI | LSI |
| PAL3-Kit-uSPE54 | <p>uSPE hardware kit for PAL RTC, consisting of:</p> <ul style="list-style-type: none"> 1 pc uSPE Liquid Tool D8/57 with Needle Guide for uSPE cartridges 2 pc Liquid smart syringe 1000μL (SF1000-57-T-22-FL), gauge 22, flat tip 1 pc Trayholder 2 pc VT54 rack (PAL3-Rack-VT54) 1 pc uSPE Cartridge Tray, 54 position, aluminum 1 pc uSPE Elution Tray Cover, 54 position aluminum 1 pc Waste Receptacle with drain 4 pc sample cartridges for initial installation, teaching USB stick with installation guide <p>(requires PAL RTC System, 2 mL vials with split septa, uSPE cartridges, Firmware 2.5 or higher)</p> | X | | |

μSPE Options

| Part No. | Description | Compatible with: | | |
|------------------------|---|------------------|-----|-----|
| | | RTC | RSI | LSI |
| Cart-uSPE-LC-QUECHERS | Consumables for LC QuEChERS for 216 Samples (four loads of VT54), incl LC Cartridges 30mg, 1.5 mL Clear Glass Vials 12x32mm and Magnetic Screw Caps | X | | |
| Cart-uSPE-GC-QUECHERS | Consumables for GC QuEChERS for 216 Samples (four loads of VT54), incl GC Cartridges 45mg, 1.5 mL Clear Glass Vials 12x32mm and Magnetic Screw Caps | X | | |
| Cart-uSPE-LC-QUE-0.3mL | Consumables for LC Quechers for 216 Samples (four loads of VT54), incl Cartridges, 0.3mL PP Vials with conical insert, and Caps | X | | |
| Cart-uSPE-GC-QUE-0.3mL | Consumables for GC Quechers for 216 Samples (four loads of VT54), incl Cartridges, 0.3mL PP Vials with conical insert, and Caps | X | | |

Recommended Option for μSPE

| Part No. | Description | |
|---------------|---|---|
| PAL3-WashFast | <p>Fast Wash Module for two solvents (Needle Gauge 19 to 26), consisting of:</p> <ul style="list-style-type: none"> 1 pc Wash Module with two wash positions and one Waste position 1 pc Pump Module with two solvent pumps 2 pc 1 liter solvent bottles with tubings and frits 1 pc PALbus connection cable <p>requires firmware version 2.0.1 or higher</p> |  |



Accessories

LC Valves:

| Part No. | Description | Compatible with: | | |
|------------------------|---|------------------|-----|-----|
| | | RTC | RSI | LSI |
| PAL3-C2V-2006D-CTC-K | 6-port Cheminert Injection Valve 0.40mm bore size incl. Needle Guide, Needle Seal G22, PAL3-Kit LCVlv | X | X | X |
| PAL3-C2V-3606D-CTC-K | 6-port Cheminert Injection Valve, 0.75mm bore size incl. Needle Seal (Gauge 22) and PAL3-Kit LCVlv Prep Application, Flow Rate 5 to 100mL/min | X | X | X |
| PAL3-C5F-1004D-CTC | Selector Valve 8-port 0.25mm, 1000psi | X | X | X |
| PAL3-C72VU-6676D-CTC-K | Cheminert-Valve 6-Port 0.15mm 18'500 psi inkl. Needle Guide, Needle Seal | X | X | X |
| PAL3-C72VX-1670D-CTC-K | Cheminert-Valve 10-Port 0.25mm 15'000 psi | X | X | X |
| PAL3-C72VX-1676D-CTC-K | Cheminert-Valve 6-Port 0.25mm 15'000 psi inkl. Needle Guide, Needle Seal | X | X | X |
| PAL3-C82VX-1676D-CTC-K | 6-port Cheminert Injection Valve 15»000 psi 75 °C 0.25 mm bore size | X | X | X |
| PAL3-ColSel6 | Column Selector Valve for up to 6 columns, consists of 2 pc 6-port Cheminert valves, bore size 0.25mm | X | X | X |
| PAL3-ColumnHolder | Column Holder for Valve Drive with clamps for columns OD 4.6mm | X | X | X |
| PAL3-DCI4WK.5/1-CTC-K | Valco W-Type LC-Valve 4-Port Internal sample loop 0.5/1.0 µL 5000 psi 50 °C, 0.40 mm bore size incl. Needle Guide, Needle Seal and Kit LCVlv | X | X | X |
| PAL3-LCInj-G22 | Standard Needle Guide for Rheodyne/ Vici valves, Gauge 22 | X | X | X |
| PAL3-LCInj-V-G19 | Standard Needle Guide for Vici valves, Gauge 19 | X | X | X |
| PAL3-LCInj20mL-G19 | Needle Guide for Vici valves, Gauge 19 Designed for tools with 20mL magnet | X | X | X |
| PAL3-LCInj20mL-G22 | Needle Guide for Rheodyne/ Vici valves, Gauge 22 Designed for tools with 20mL magnet or foil cutter | X | X | X |
| PAL3-PD718-110-K | 6-port Rheodyne Injection Valve 0.28mm bore size and MBB functionality, stainless steel rotor seal (RPC-13 coated), Pressure up to 1200 bar, incl. Needle Guide, Needle Seal Without Sample Loop | X | X | X |
| PAL3-PD718-202 | 10-port Rheodyne Injection Valve 0.15mm bore size SPC3 coated stator and RPC-13 coated rotor seal Pressure up to 1200 bar | X | X | X |

Mounting Kits:

| Part No. | Description | Compatible with: | | |
|---------------------|---|------------------|-----|-----|
| | | RTC | RSI | LSI |
| PAL3-Kit-5890 | Mounting Kit for Agilent GC model 5890 | X | X | X |
| PAL3-Kit-5975T | Mounting Kit for Agilent GC model 5975T | X | X | X |
| PAL3-Kit-6850 | Mounting Kit for Agilent GC model 6850 | X | X | X |
| PAL3-Kit-7820-D | Mounting Kit for Agilent GC model 7820 (deck mounting) | X | X | X |
| PAL3-Kit-7890-D | Mounting Kit for Agilent GC models 6890 / 7890 (deck mounting) | X | X | X |
| PAL3-Kit-7890-S | Mounting Kit for Agilent GC models 6890 / 7890 (side mounting) | X | X | X |
| PAL3-Kit-88xx-D | Mounting Kit for Agilent GC 6890/7980/8860/8890 with longer legs (deck mounting) | X | X | X |
| PAL3-Kit-88xx-S | Mounting Kit for Agilent GC 6890/7980/8860/8890 with longer legs (side mounting) | X | X | X |
| PAL3-Kit-Bruker-436 | Mounting Kit for Bruker GC model 436 (compatible with Varian GC V3900) | X | X | X |
| PAL3-Kit-Bruker-451 | Mounting Kit for Bruker GC model 451 | X | X | X |
| PAL3-Kit-Bruker-456 | Mounting Kit for Bruker GC model 456 | X | X | X |
| PAL3-Kit-Focus | Mounting kit for Thermo Focus GC | X | X | X |
| PAL3-Kit-Intuvo-L | Mounting Kit for Agilent Intuvo 9000 GC with additional long leg | X | X | X |
| PAL3-Kit-MasterGC | Mounting Kit for DANI Master GC | X | X | X |

| Part No. | Description | Compatible with: | | |
|----------------------|---|------------------|-----|-----|
| | | RTC | RSI | LSI |
| PAL3-Kit-PEC500 | Mounting Kit for PERKIN ELMER Clarus 500 | X | X | X |
| PAL3-Kit-PEC600 | Mounting Kit for PERKIN ELMER Clarus 600 / Clarus 680-iQT | X | X | X |
| PAL3-Kit-S17A | Mounting Kit for Shimadzu GC-S17A model | X | X | X |
| PAL3-Kit-S2010 | Mounting Kit for Shimadzu GC-2010 model | X | X | X |
| PAL3-Kit-S2014 | Mounting Kit for Shimadzu GC-2014 model | X | X | X |
| PAL3-Kit-S2025 | Mounting Kit for Shimadzu GC-2025 model | X | X | X |
| PAL3-Kit-S2030 | Mounting Kit for Shimadzu GC-2030 model | X | X | X |
| PAL3-Kit-Standalone | Standard Support PAL3, 215mm, incl. PAL Interface Cable, 2pcs | X | X | X |
| PAL3-Kit-Trace | Mounting Kit for Thermo Trace GC Ultra | X | X | X |
| PAL3-Kit-Trace1300 | Mounting Kit for Thermo Trace 1300/1310 GC | X | X | X |
| PAL3-Kit-Varian-3800 | Mounting Kit for Varian GC model 3800 | X | X | X |
| PAL3-Kit-YL6500 | Mounting Kit for YoungLin 6500 GC | X | X | X |

Other Accessories

| Part No. | Description | Compatible with: | | |
|----------------------|---|------------------|-----|-----|
| | | RTC | RSI | LSI |
| PAL3-Agi-InsSet-10ml | Set of 6x inserts for the Agitator to place 10mL vials instead of 20mL vials | X | X | X |
| PAL3-Agi-InsSet-2ml | Set of 6x inserts for the Agitator to place 2mL vials instead of 20mL vials | X | X | X |
| PAL3-DilAspKit | Aspiration Tubing Kit, consisting of: 1 pc Tube PFA 1 pc 28mm VitraPOR Glass Frit 40um pore size 1 pc Plug with Feed-through Holes | X | X | |
| PAL3-FoilCutter | Foil cutter to pierce either aluminum or plastic foils prior to the needle penetration into Microtiter or Deep Well plates (Requires PAL3-LCInj-V-G19 or PAL3-LCInj-G22) | X | X | |



Consumables - Vials, Vial Caps and Seals

Vials

| Part No. | Description | Compatible with: | | |
|---------------------|---|------------------|-----|-----|
| | | RTC | RSI | LSI |
| Vial-0.3-ND9-PP-100 | PAL System Micro-Vial 2CV, 0.3mL Clear PP, designed for the PAL Autosampler. 12x32mm, with conical insert, fits ND9 Screw Caps, Pk of 100 Pcs | X | X | X |
| Vial-1.5-ND9-AG-100 | PAL System Vial 2CV, 1.5ml Amber Glass with Label, designed for the PAL Autosampler. 12x32mm, 1st Class Hydrolytic Glass, fits ND9 Screw Caps, Pk of 100 Pcs | X | X | X |
| Vial-1.5-ND9-CG-100 | PAL System Vial 2CV, 1.5ml Clear Glass with Label, designed for the PAL Autosampler. 12x32mm, 1st Class Hydrolytic Glass, fits ND9 Screw Caps, Pk of 100 Pcs | X | X | X |
| Vial-1.5-ND9-PP-100 | PAL System Vial 2CV, 1.5mL Clear PP, designed for the PAL Autosampler. 12x32mm, with filling lines, fits ND9 Screw Caps, Pk of 100 Pcs | X | X | X |
| Vial-10-DC20-CG-100 | PAL System Vial 10CV, 10ml Clear Glass with Label, for PAL Wash Modules. 46x22.5mm, 1st Class Hydrolytic Glass, DIN-Crimp Neck, Pk of 100pcs | X | X | X |
| Vial-10-ND18-AG-100 | PAL System Vial 10CV, 10ml Amber Glass with Label, designed for the PAL Autosampler. 46x22.5mm, 1st Class Hydrolytic Glass, with flat finish for better sealing, fits ND18 Screw Caps, Pk of 100 Pcs | X | X | X |
| Vial-10-ND18-CG-100 | PAL System Vial 10CV, 10ml Clear Glass with Label, designed for the PAL Autosampler. 46x22.5mm, 1st Class Hydrolytic Glass, with flat finish for better sealing, fits ND18 Screw Caps, Pk of 100 Pcs | X | X | X |
| Vial-20-DC20-AG-100 | PAL System Vial 20CV, 20ml Amber Glass for Headspace. 75.5x22.5mm, 1st Class Hydrolytic Glass, DIN-Crimp Neck, Pk of 100pcs | X | X | X |
| Vial-20-DC20-CG-100 | PAL System Vial 20CV, 20ml Clear Glass for Headspace. 75.5x22.5mm, 1st Class Hydrolytic Glass, DIN-Crimp Neck, Pk of 100pcs | X | X | X |
| Vial-20-ND18-AG-100 | PAL System Vial 20CV, 20ml Amber Glass with Label, designed for the PAL Autosampler. 75.5x22.5mm, 1st Class Hydrolytic Glass with flat finish for better sealing, fits ND18 Screw Caps, Pk of 100 Pcs | X | X | X |
| Vial-20-ND18-CG-100 | PAL System Vial 20CV, 20ml Clear Glass with Label, designed for the PAL Autosampler. 75.5x22.5mm, 1st Class Hydrolytic Glass with flat finish for better sealing, fits ND18 Screw Caps, Pk of 100 Pcs | X | X | X |
| Vial-40-ND24-CG-100 | PAL System Vial 40CV, 40ml Clear Glass with Label, designed for the PAL Autosampler. 95x27.5mm, 1st Class Hydrolytic Glass with flat top for optimized sealing. For ND24 Screw Caps, Pk of 100 Pcs | X | X | X |

Vial Sample Set for trial

| Part No. | Description | Compatible with: | | |
|----------------------|---|------------------|-----|-----|
| | | RTC | RSI | LSI |
| PAL3-Vial-Sample-Set | Vial and Cap Sample Set consisting of: 2pc 1.5ml transparent clear glass vial, fits ND9 Screw Cap 2pc 1.5ml amber glass vial, fits ND9 Screw Cap 2pc 10ml transparent clear glass vial, fits ND18 Screw Cap 2pc 10ml amber glass vial, fits ND18 Screw Cap 2pc 20ml transparent clear glass vial, fits ND18 Screw Cap 2pc 20ml amber glass vial, fits ND18 Screw Cap 4pc Screw Cap, ND9, Magnetic, Slitted Septum 1.0mm 4pc Screw Cap, ND18, Magnetic 4pc Screw Re-Cap, ND18, Magnetic | X | X | X |

Vial Caps

| Part No. | Description | RTC | RSI | LSI |
|------------------------|--|-----|-----|-----|
| Cap-ND9-PP-SP10-100 | PAL System Screw Cap 2CV, designed for the PAL Autosampler. ND9, Silicone/PTFE Septa 1.0mm, Pk of 100 pcs | X | X | X |
| Cap-ND9-PP-SP10S-100 | PAL System Screw Cap 2CV, designed for the PAL Autosampler. ND9, Silicone/PTFE/I-Slitted Septa 1.0mm, Pk of 100 pcs | X | X | X |
| Cap-ND9-St-SP10-100 | PAL System Screw Cap 2CV, designed for the PAL Autosampler. ND9, Magnetic, Gold, Silicone/PTFE Septa 1.0mm, Pk of 100 pcs | X | X | X |
| Cap-ND9-St-SP10Sb-100 | PAL System Screw Cap 2CV, designed for the PAL Autosampler. ND9, Magnetic, Gold, Silicone/PTFE/Starburst-Slitted Septa 1.0mm, Pk of 100 pcs | X | X | X |
| Cap-ND18-StC-SP20-100 | PAL System Screw Re-Cap 10CV and 20CV, designed for the PAL Autosampler. ND18 Magnetic, PTFE Metal Coating for optimized screw and sealing performance, Silicone/PTFE Septa 2.0mm, Pk of 100 pcs | X | X | X |
| Cap-ND18-St-SP15-100 | PAL System Screw Cap 10CV and 20CV, designed for the PAL Autosampler. ND18 Magnetic, Silver, Silicone white/PTFE blue Septa 1.5mm, Pk of 100 pcs | X | X | X |
| Cap-ND24-PP-SP32-100 | PAL System Screw Cap 40CV, designed for the PAL Autosampler. ND24, Silicone/PTFE Septa 3.2mm, Pk of 100 pcs | X | X | X |
| Cap-DC20-PE8-SP13Y-100 | "PAL System Clip Cap for PAL Wash Modules. 22mm PE with 8mm center hole, Silicone/PTFE/Y-Slitted Septa 1.3mm, Pk of 100 pcs" | X | X | X |
| Cap-DC20-St-SP15-100 | PAL System Crimp Cap for Headspace, SPME Fiber and SPME Arrow. 20mm with 8mm Center hole, Silikon / PTFE Septa 3.0mm with thinned center 1.5mm. Pack of 100pcs | X | X | X |
| Cap-DC20-St-SP30-100 | PAL System Crimp Cap for Headspace and SPME Arrow. 20mm with 8mm Center hole, Silikon / PTFE Septa 3.0mm. Pack of 100pcs | X | X | X |

Vial and Cap Kit

| Part No. | Description | RTC | RSI | LSI |
|------------------------|---|-----|-----|-----|
| PAL3-Kit-WashVial-Caps | Snap cap (PE) and septum (Silicon/PTFE) for 10mL/100mL vials for Standard Wash Module (set of 10) | X | X | X |
| PAL3-Kit-WashVial-10ml | 10mL solvent vial for Standard Wash Module completed with snap caps and seals (set of 10) | X | X | X |

Needle Seal

| Part No. | Description | RTC | RSI | LSI |
|--------------|--|-----|-----|-----|
| PAL3-NdlSeal | PALseal gauge 22 for Valco and Rheodyne valves (use together with Needle Guide NV-30-25) Pack of 5 pcs | X | X | X |

Spare Parts A-Z



| Part No. | Description | Compatible with: | | |
|------------------------|--|------------------|-----|-----|
| | | RTC | RSI | LSI |
| PAL3-Adap-7890/6890 | Adapter for Agilent GC 7890/6890 | X | X | X |
| PAL3-Adap-GC2010 | Adapter for Shimadzu GC2010 Plus | X | X | X |
| PAL3-Adap-GC2030 | Adapter for Shimadzu GC2030 | X | X | X |
| PAL3-Adap-Trace1300-GC | Adapter for Thermo Trace 1300 GC injectors, suitable for SSL or PTV injector | X | X | X |
| PAL3-Adap-TraceUltra | Adapter for Thermo Trace 1300 Ultra injector | X | X | X |
| PAL3-AdapStdWash-2ml | Adapter for Standard Wash Station to place 2mL vial instead of the 10mL wash vial | X | X | X |
| PAL3-AdapStdWash-Waste | Waste Adapter for Standard Wash Module | X | X | X |
| PAL3-Adapters-TraceGC | Kit of adapters for Thermo Trace GC Ultra injectors: consists of one adapter for SSL/packed/purge packed injector and one adapter for PTV | X | X | X |
| PAL3-Adapters-FocusGC | Adapter kit for Thermo Focus GC injectors: consists of one adapter for SSL/packed/purge packed injector | X | X | X |
| PAL3-Agi-TeachAdap | Agitator Teaching Tool, 2pcs | X | X | X |
| PAL3-HS1000-Oring-kit | Spare for Headspace Tool, consisting of: 4pc O-ring 7x1.5 MVQ for syringe fixation 2pc O-ring 3x1.0 MVQ for gas connection 1pc syringe removing tool 1pc syringe O-ring removing tool | X | X | |
| PAL3-HS2500-Oring-kit | Spare for Headspace Tool, consisting of: 4pc O-ring 9.25x1.78 MVQ for syringe fixation 2pc O-ring 3x1.0 MVQ for gas connection 1pc syringe removing tool 1pc syringe O-ring removing tool | X | X | |
| PAL3-HS5000-Oring-kit | Spare for Headspace Tool, consisting of: 4pc O-ring 13x1.5 MVQ for syringe fixation 2pc O-ring 3x1.0 MVQ for gas connection 1pc syringe removing tool 1pc syringe O-ring removing tool | X | X | |
| PAL3-KIT-LCVlv | Kit Waste-Tubing for VICI and Rheodyne Valves | X | X | X |
| PAL3-LCMSLoop | Kit PAL3 LCMS Holding Loop with Fittings | X | X | |
| PAL3-LSI-PM-Kit | Preventive Maintenance Kit for PAL LSI, consisting of: 1pc Vial Teaching Kit 1pc Test Label Barcode 1pc Lubricant 1pc Brush 3pc Spare Screw for fixation of PALhead 1pc Spare Screw for fixation of Cover 1pc Kit Wash Vials Soft Tissue for outer surface cleaning Rough Tissue for degreasing X-,Y-,Z-Axes or Guide Rollers | | | X |
| PAL3-LubKit | PAL Lubrication Kit including brush | X | X | X |

| Part No. | Description | Compatible with: | | |
|---------------------|---|------------------|-----|-----|
| | | RTC | RSI | LSI |
| PAL3-PM-GC-KIT | Preventive Maintenance Kit for PALSystem GC configuration, consisting of: 1pc Vial Teaching Kit 2pc Driving Belt Agitator 1pc Test Label Barcode 1pc Lubricant 1pc Brush 1pc Spare Screw for fixation of PALhead cover 1pc Spare Screw for fixation of PALhead with O-Ring 1pc Kit Wash Vials Soft Tissue for outer surface cleaning Rough Tissue for degreasing X-,Y-,Z-Axes or Guide Rollers | X | X | |
| PAL3-PM-LC-KIT | Preventive Maintenance Kit for PALSystem LC configuration, consisting of: 1pc Vial Teaching Kit 1pc Test Label Barcode 1pc Lubricant 1pc Brush 5pc PALseal 1pc Spare Screw for fixation of PALhead cover 1pc Spare Screw for fixation of PALhead with O-Ring 1pc Kit Waste-Tubing for VICI and Rheodyne Valves Soft Tissue for outer surface cleaning Rough Tissue for degreasing X-,Y-,Z-Axes or Guide Rollers | X | X | |
| PAL3-POWERSUPPLY | Main Power supply for PALSystem. Input Voltage 100-240V | X | X | X |
| PAL3-TubeDilutor | Transfer tubing kit for Dilutor | X | X | |
| PAL3-TubeFastWash | Tubing Kit for Fast Wash Module | X | X | X |
| PAL3-WashVial-100ml | 100mL solvent vial for Large Wash Module or Solvent Module completed with snap cap (PE) and septum (Silicon/PTFE, one piece) | X | X | X |

PAL SYSTEM

Ingenious sample handling



Distributed by:

For more information on the PAL System visit:

www.palsystem.com



PAL is a registered trademark of CTC Analytics AG, Switzerland