


# EPA Method 8270D Analysis in 10.1 Minutes of an 80 Component Calibration Standard by GC-TOFMS

*Life Science and Chemical Analysis Centre*  
LECO Corporation, Saint Joseph, Michigan USA



**Pegasus<sup>®</sup> HT GC-TOFMS of an 80 component calibration standard by EPA Method**

**Column:**

10 m x 0.18 mm x 0.18  $\mu$ m Rtx-5SilMS

**TOFMS:**

35 to 500 m/z @ 20 spectra/s

For questions on this analysis e-mail us at:  
[life\\_science@leco.com](mailto:life_science@leco.com)

For a complete listing of snapshots and application notes visit us on the web at [www.leco.com](http://www.leco.com)

