



12

Scheren mit Hartmetallschneiden

Scissors with tungsten carbide edges

Tijeras con filos de carburo tungsteno

Forbici con lame in carburo i tungsteno



Kennzeichnung: goldene Ringe
Characteristic: golden rings
Característica: mangos dorados
Caratteristica: anelli dorata





Ausführungen I Specification
Especificación I Specifica




1 Schneide gezahnt
1 serrated blade
1 hoja dentada
1 lama dentata

Schneidegeometrie
Cutting geometry
Geometría del corte
Geometria di taglio



- Aufgeschweißte Hartmetallschneiden, daher keine zur Korrosion neigende Lötnaht
- Dauerhaft scharfer Schnitt
- Gleichmäßiger, weicher Gang
- Welded tungsten carbide cutting edges, therefore no seam tending to corrosion
- Durable sharp cutting edge
- Regular and smooth function
- Filos de carburo de tungsteno soldados, por eso sin soldadura tendiente a la corrosión.
- Agudeza del corte durable
- Corte suave y regular
- Parti taglienti in carburo di tungsteno saldate e di conseguenza nessuna tendenza alla corrosione
- Taglio eccellente di lunga durata
- Funzionamento armonico e dolce






12-000-14
12-000-16
12-000-18



STANDARD Fig. 1

12-000-14	14,0 cm, 5 ½ "
12-000-16	16,0 cm, 6 ¼ "
12-000-18	18,0 cm, 7 "



12-001-14
12-001-16
12-001-18




STANDARD Fig. 1

12-001-14	14,0 cm, 5 ½ "
12-001-16	16,0 cm, 6 ¼ "
12-001-18	18,0 cm, 7 "

STANDARD Fig. 1

12-001-14 W	14,0 cm, 5 ½ "
--------------------	----------------




12-002-11
12-002-14
12-002-16
12-002-18



STANDARD Fig. 2

12-002-11	11,0 cm, 4 ⅜ "
12-002-14	14,0 cm, 5 ½ "
12-002-16	16,0 cm, 6 ¼ "
12-002-18	18,0 cm, 7 "




12-003-14
12-003-16
12-003-18



STANDARD Fig. 2

12-003-14	14,0 cm, 5 ½ "
12-003-16	16,0 cm, 6 ¼ "
12-003-18	18,0 cm, 7 "




12-004-14
12-004-16
12-004-18



STANDARD Fig. 3

12-004-14	14,0 cm, 5 ½"
12-004-16	16,0 cm, 6 ¼"
12-004-18	18,0 cm, 7"



12-005-14
12-005-16
12-005-18



STANDARD Fig. 3

12-005-14	14,0 cm, 5 ½"
12-005-16	16,0 cm, 6 ¼"
12-005-18	18,0 cm, 7"




12-107-15



MIXTER

12-107-15	15,0 cm, 6"
------------------	-------------




12-141-20



WALDMANN

12-141-20	20,0 cm, 8"
------------------	-------------

Damm Schere (Episiotomie)
Perineum scissors (Episiotomy)
Tijeras para perineo (episiotomía)
Forbici per perineo (episiotomia)



12-120-20
12-120-23




12-121-20
12-121-23



SIMS

12-120-20 - 12-121-20 20,0 cm, 8 "

12-120-23 - 12-121-23 23,0 cm, 9 "



12-120-20 W




12-121-20 W
12-121-23 W



SIMS

12-120-20 W - 12-121-20 W 20,0 cm, 8 "

12-121-23 W 23,0 cm, 9 "



12-122-20
12-122-23



12-123-20
12-123-23



SIMS

12-122-20 - 12-123-20 20,0 cm, 8 "

12-122-23 - 12-123-23 23,0 cm, 9 "



12-124-20
12-124-23




12-125-20
12-125-23



SIMS

12-124-20 - 12-125-20 20,0 cm, 8 "

12-124-23 - 12-125-23 23,0 cm, 9 "




12-160-14
12-160-17
12-160-23



MAYO

12-160-14	14,0 cm, 5 ½ "
12-160-17	17,0 cm, 6 ¾ "
12-160-23	23,0 cm, 9 "




12-161-14
12-161-17
12-161-23
12-161-36
12-161-40



MAYO

12-161-14	14,0 cm, 5 ½ "
12-161-17	17,0 cm, 6 ¾ "
12-161-23	23,0 cm, 9 "
12-161-36	36,0 cm, 14 ¼ "
12-161-40	40,0 cm, 16 "



12-161-17 W



MAYO

12-161-17 W	17,0 cm, 6 ¾ "
--------------------	----------------




12-166-22




MAYO-HARRINGTON

12-166-22	22,0 cm, 8 ¾ "
------------------	----------------



12-167-22



12-169-23




12-167-23



12-167-22 22,0 cm, 8 $\frac{3}{4}$ "
12-169-23 23,0 cm, 9 "

WERTHEIM


12-167-23 23,0 cm, 9 "



12-173-25 25,0 cm, 9 $\frac{3}{4}$ "



12-175-25 25,0 cm, 9 $\frac{3}{4}$ "




12-170-15
12-170-17
12-170-21



MAYO-STILLE

12-170-15	15,0 cm, 6 "
12-170-17	17,0 cm, 6 ¾ "
12-170-21	21,0 cm, 8 ¼ "




12-171-15
12-171-17
12-171-21



MAYO-STILLE

12-171-15	15,0 cm, 6 "
12-171-17	17,0 cm, 6 ¾ "
12-171-21	21,0 cm, 8 ¼ "





METZENBAUM

12-200-11	11,0 cm, 4 ⅜ "
12-200-14	14,0 cm, 5 ½ "
12-200-18	18,0 cm, 7 "
12-200-20	20,0 cm, 8 "
12-200-23	23,0 cm, 9 "
12-200-25	25,0 cm, 9 ¾ "
12-200-28	28,0 cm, 11 "
12-200-30	30,0 cm, 12 "





METZENBAUM

12-201-11	11,0 cm, 4 ⅜ "
12-201-14	14,0 cm, 5 ½ "
12-201-18	18,0 cm, 7 "
12-201-20	20,0 cm, 8 "
12-201-23	23,0 cm, 9 "
12-201-25	25,0 cm, 9 ¾ "
12-201-28	28,0 cm, 11 "
12-201-30	30,0 cm, 12 "
12-201-36	36,0 cm, 14 ¼ "
12-201-40	40,0 cm, 16 "





METZENBAUM

12-201-11 W	11,0 cm, 4 $\frac{3}{8}$ "
12-201-14 W	14,0 cm, 5 $\frac{1}{2}$ "
12-201-18 W	18,0 cm, 7 "
12-201-20 W	20,0 cm, 8 "
12-201-23 W	23,0 cm, 9 "



METZENBAUM-S

12-209-20	20,0 cm, 8 "
12-209-23	23,0 cm, 9 "
12-209-25	25,0 cm, 9 $\frac{3}{4}$ "
12-209-28	28,0 cm, 11 "





12-207-23



METZENBAUM

12-207-23 23,0 cm, 9 "





12-207-23 W



METZENBAUM

12-207-23 W 23,0 cm, 9 "



METZENBAUM-FINO

12-210-14	14,0 cm, 5 1/2 "
12-210-18	18,0 cm, 7 "
12-210-20	20,0 cm, 8 "
12-210-23	23,0 cm, 9 "
12-210-25	25,0 cm, 9 3/4 "
12-210-28	28,0 cm, 11 "





METZENBAUM-FINO

12-211-14	14,0 cm, 5 1/2 "
12-211-18	18,0 cm, 7 "
12-211-20	20,0 cm, 8 "
12-211-23	23,0 cm, 9 "
12-211-25	25,0 cm, 9 3/4 "
12-211-28	28,0 cm, 11 "




METZENBAUM-FINO

12-211-14 W	14,0 cm, 5 1/2 "
12-211-18 W	18,0 cm, 7 "
12-211-20 W	20,0 cm, 8 "
12-211-23 W	23,0 cm, 9 "
12-211-25 W	25,0 cm, 9 3/4 "
12-211-28 W	28,0 cm, 11 "





METZENBAUM-FINO

12-214-14	14,0 cm, 5 1/2 "
12-214-18	18,0 cm, 7 "
12-214-20	20,0 cm, 8 "
12-214-23	23,0 cm, 9 "





METZENBAUM-FINO

12-215-11	11,0 cm, 4 3/8 "
12-215-14	14,0 cm, 5 1/2 "
12-215-18	18,0 cm, 7 "
12-215-20	20,0 cm, 8 "
12-215-23	23,0 cm, 9 "





METZENBAUM-FINO-S

12-221-18	18,0 cm, 7 "
12-221-20	20,0 cm, 8 "
12-221-23	23,0 cm, 9 "



METZENBAUM-FINO-S

12-225-18	18,0 cm, 7 "
12-225-20	20,0 cm, 8 "
12-225-23	23,0 cm, 9 "



12-230-16
12-230-21





12-231-16
12-231-21



MAYO-LEXER

12-230-16 - 12-231-16 16,0 cm, 6 ¼"

12-230-21 - 12-231-21 21,0 cm, 8 ¼"





12-233-16 W



LEXER

12-233-16 W 16,0 cm, 6 ¼"



12-232-16




12-233-16



LEXER

12-232-16 - 12-233-16 16,0 cm, 6 ¼"



12-241-20
12-241-23
12-241-28



THOREK

12-241-20 20,0 cm, 8"

12-241-23 23,0 cm, 9"

12-241-28 28,0 cm, 11"



12-243-14



JOSEPH

12-243-14 14,5 cm, 5 $\frac{3}{4}$ "



12-245-12
12-245-15



KILNER

12-245-12 12,0 cm, 4 $\frac{3}{4}$ "
12-245-15 15,0 cm, 5 "




12-247-13
12-247-15



JAMESON

12-247-13 13,0 cm, 5 $\frac{1}{8}$ "
12-247-15 15,0 cm, 5 "




12-248-15



REYNOLDS

12-248-15 15,0 cm, 5 "




12-249-15
12-249-18



REYNOLDS

12-249-15 15,0 cm, 5 "
12-249-18 18,0 cm, 7 "




12-253-14



METZENBAUM-LAHEY

12-253-14 14,0 cm, 5 1/2 "



12-253-20
12-253-23



METZENBAUM-LAHEY

12-253-20 20,0 cm, 8 "
12-253-23 23,0 cm, 9 "



12-253-17




TOENNIS-ADSON

12-253-17 17,5 cm, 7 "



TOENNIS-ADSON

12-253-17W 17,5 cm, 7 "




12-305-17



Tonsillenscheren
Tonsil scissors
Tijera para amígdalas
Forbici per tonsille

DEAN

12-305-17 17,0 cm, 6 $\frac{3}{4}$ "



12-319-13



Nasenscheren
Nasal scissors
Tijera nasales
Forbici nasali

FOMON

12-319-13 13,0 cm, 5 $\frac{1}{8}$ "




12-329-13



Nasenscheren
Nasal scissors
Tijera nasales
Forbici nasali

FOMON

12-329-13 13,0 cm, 5 $\frac{1}{8}$ "





12-331-16



Nasenscheren
Nasal scissors
Tijera nasales
Forbici nasali

COTTLE

12-331-16 16,0 cm, 6 $\frac{1}{4}$ "



12-360-20


12-360-23



GORNEY

12-360-20 20,0 cm, 8 "

12-360-23 23,0 cm, 9 "




12-363-14



GREGORY

12-363-14 14,0 cm, 5 1/2 "



12-365-15

12-365-18



KAYE

12-365-15 14,5 cm, 5 3/4 "

12-365-18 18,0 cm, 7 "




12-367-18



STELLA

12-367-18 18,0 cm, 7 "




12-369-18



GORNEY-FREEMAN

12-369-18 18,0 cm, 7 "




12-370-23



GORNEY-FREEMAN

12-370-23 23,0 cm, 9 "




12-371-18



FREEMAN-KAYE

12-371-18 18,0 cm, 7 "





12-371-23




FREEMAN-GORNEY

12-371-23 23,0 cm, 9 "



12-400-16

12-401-16



12-404-14

12-405-14




KELLY

12-400-16 - 12-401-16 16,0 cm, 6 1/4"

JOSEPH

12-404-14 - 12-405-14 14,0 cm, 5 1/2"



12-416-13


12-417-13



GOLDMANN-FOX

12-416-13 - 12-417-13 13,0 cm, 5 1/8"

Zahnfleischscheren
Gum scissors
Tijeras para encías
Forbici per gengiva




12-424-12

12-425-12




WAGNER


12-424-12 - 12-425-12 12,0 cm, 4 3/4"



12-436-11



12-437-11



12-438-11





12-437-11 W



IRIS

12-436-11 - 12-438-11 11,5 cm, 4 ½ "

12-437-11 W 17,0 cm, 6 ¾ "




12-519-16



Zahnfleischschere
Gum scissors
Tijera para encías
Forbici per gengiva

LOCKLIN

12-519-16 16,0 cm, 6 ¼ "




12-521-16



Zahnfleischschere
Gum scissors
Tijera para encías
Forbici per gengiva

LOCKLIN-S

12-521-16 16,0 cm, 6 ¼ "



DE BAKEY



12-541-16 - 12-545-16 16,0 cm, 6 ¼ "

12-541-23 - 12-545-23 23,0 cm, 9 "



POTTS-SMITH


12-551-19 - 12-555-19 19,0 cm, 7 ½ "



Gefäßscheren
Vascular scissors
Tijeras vasculares
Forbici vascolari



HEGEMANN-DIETRICH

12-571-18 - 12-579-18 18,0 cm, 7 "




UNIVERSAL

12-721-12 12,0 cm, 4 ¾ "



SMITH

12-722-16 - 12-724-16 16,0 cm, 6 ¼ "




12-731-14
12-731-17
12-731-23



MAYO

12-731-14 14,0 cm, 5 ½ "
12-731-17 17,0 cm, 6 ¾ "
12-731-23 23,0 cm, 9 "



12-731-14 W
12-731-17 W
12-731-23 W



MAYO

12-731-14 W 14,0 cm, 5 ½ "
12-731-17 W 17,0 cm, 6 ¾ "
12-731-23 W 23,0 cm, 9 "



12-733-14 W
12-733-18 W
12-733-20 W
12-733-23 W



METZENBAUM

12-733-14 W 14,0 cm, 5 ½ "
12-733-18 W 18,0 cm, 7 "
12-733-20 W 20,0 cm, 8 "
12-733-23 W 23,0 cm, 9 "



12-733-26 W



METZENBAUM

12-733-26 W 26,0 cm, 10 ¼ "

