


Sol-Care™ Safety Needle

Canadian Medical Surgical Supplies
Tel: 647-514-6290, F 905-890-4912



Canadian Medical Surgical Supplies
Tel: 647-514-6290, F 905-890-4912

Sol-Care™ Safety Needle

Sol-Care™ Safety Needles are designed to help prevent needlestick injuries. Every needle is equipped with a pivoting safety shield mechanism that locks in place to lower the risk of accidental injury after medication administration. There are two methods of safety shield activation that help to keep hands and fingers behind and away from the needle tip for protection from a contaminated needle.

Features:


- Two options for protective arm activation: hard surface activation and thumb activation
- Simple one-handed safety shield mechanism
- Audible click indicates that the safety mechanism is secured and locked in place
- Hubs, protective shields, and packaging are color-coded for facilitating the identification of needles gauges (ISO 6009)
- Protective arm designed not to obstruct visualization of injection site
- Bevel-up needle allows for low-angle injections
- Compatible with all standard Luer Lock and Luer Slip Tip syringes (ISO 80369-7)
- Designed according to ISO 7864 and 23908
- Not manufactured with Latex, PVC and DEHP



Color-coded
for facilitating
the identification
needles gauges



Locking
mechanism
designed to
secure
the needle



Bevel-up
needle allows
low-angle
injections

CE 2797

Sol-Care™ Safety Needle



NEEDLE

Ref. Number	Product Description	Gauge	Length (inches)	Length (mm)	Box Qty	Case Qty
SN1810	Sol-Care Safety Needle 18G x 1"	18G	1"	25 mm	100	800
SN1815	Sol-Care Safety Needle 18G x 1 1/2"	18G	1 1/2"	40 mm	100	800
SN1820	Sol-Care Safety Needle 18G x 2"	18G	2"	50 mm	100	800
SN1910	Sol-Care Safety Needle 19G x 1"	19 G	1"	25 mm	100	800
SN1915	Sol-Care Safety Needle 19G x 1 1/2"	19 G	1 1/2"	40 mm	100	800
SN1920	Sol-Care Safety Needle 19G x 2"	19 G	2"	50 mm	100	800
SN2010	Sol-Care Safety Needle 20G x 1"	20G	1"	25 mm	100	800
SN2015	Sol-Care Safety Needle 20G x 1 1/2"	20G	1 1/2"	40 mm	100	800
SN2110	Sol-Care Safety Needle 21G x 1"	21G	1"	25 mm	100	800
SN21125	Sol-Care Safety Needle 21G x 1 1/4"	21G	1 1/4"	32 mm	100	800
SN2115	Sol-Care Safety Needle 21G x 1 1/2"	21G	1 1/2"	40 mm	100	800
SN2105	Sol-Care Safety Needle 21G x 2"	21G	2"	50 mm	100	800
SN2210	Sol-Care Safety Needle 22G x 1"	22G	1"	25 mm	100	800
SN22125	Sol-Care Safety Needle 22G x 1 1/4"	22G	1 1/4"	32 mm	100	800
SN2215	Sol-Care Safety Needle 22G x 1 1/2"	22G	1 1/2"	40 mm	100	800
SN2205	Sol-Care Safety Needle 22G x 2"	22G	2"	50 mm	100	800
SN2310	Sol-Care Safety Needle 23G x 1"	23G	1"	25 mm	100	800
SN23125	Sol-Care Safety Needle 23G x 1 1/4"	23G	1 1/4"	32 mm	100	800
SN2315	Sol-Care Safety Needle 23G x 1 1/2"	23G	1 1/2"	40 mm	100	800
SN2410	Sol-Care Safety Needle 24G x 1"	24G	1"	25 mm	100	800
SN2558	Sol-Care Safety Needle 25G x 5/8"	25G	5/8"	16 mm	100	800
SN2510	Sol-Care Safety Needle 25G x 1"	25G	1"	25 mm	100	800
SN2515	Sol-Care Safety Needle 25G x 1 1/2"	25G	1 1/2"	40 mm	100	800
SN2612	Sol-Care Safety Needle 26G x 1/2"	26G	1/2"	12.5 mm	100	800
SN2712	Sol-Care Safety Needle 27G x 1/2"	27G	1/2"	12.5 mm	100	800
SN2758	Sol-Care Safety Needle 27G x 5/8"	27G	5/8"	16 mm	100	800
SN2734	Sol-Care Safety Needle 27G x 3/4"	27G	3/4"	20 mm	100	800
SN27125	Sol-Care Safety Needle 27G x 1 1/4"	27G	1 1/4"	32 mm	100	800
SN2715	Sol-Care Safety Needle 27G x 1 1/2"	27G	1 1/2"	40 mm	100	800
SN3012	Sol-Care Safety Needle 30G x 1/2"	30G	1/2"	12.5 mm	100	800

Please contact your local SOL-MILLENNIUM® representative for additional sizes of the Sol-Care™ Safety Needle.




NEEDLE



Two options for protective arm activation:



Hard Surface Activation



Thumb Activation



Intended use: Sol-Care™ Safety Needle is used in conjunction with a standard syringe. This device is used for aspiration and administration of medication.

Intended user: Licensed healthcare professionals (HCP).

Cautions:

- Single use device. Re-use or use if the package is damaged may lead to infection or other illness/ injury.
- Sharp object. Risk of sharp injuries.

The content of this material is directed at and restricted to a professional audience only (Healthcare Professionals). This material is not medical advice, guideline or a recommendation. Before using any particular product always read the package insert, product label and/or instructions for use.

This material is distributed by Sol-Millennium Europe Sp z o.o. with a registered address at Twarda 18, 00-105 Warsaw, Poland, Registry Court: District Court for Warsaw, XIIth Economic Department KRS 0000434783, EU VAT: PL5252535720, Share capital: 500.000 PLN.

Sol-Millennium Europe Sp. z o.o. All rights reserved.

For more information, please contact:

EMA-MD-488-0623



Sol-Millennium Europe Sp. z o.o.
Warsaw HUB, Rondo Daszyńskiego 2B
00-843 Warsaw, Poland



info.eu@sol-m.com



Sol M

Building
a Healthier
Tomorrow