


Passive mattress PrimaCare 10, 200 x 86 x 14

Article number: **4PPD21100AS**

Description:

- Standard mattress suitable for use in healthcare and nursing care.
- PUR, Monoblock
- The core is cut transversely and longitudinally in the upper part.
- Possibility of overturning: NO
- Option to rotate: YES
- Cutouts on the underside: NO
- Contact layer profiling: YES
- Waterproof, vapor permeable cover
- Cover binding technology: stitching
- Zipper protection: YES



Properties

BD:T2	External Dimensions: 200 x 86 x 14
CC:T1	Corner Cuts: Square
DE:T1	Product Option: Mattress
DP[3]:03	Type of foam: Other foam
FRS:T0	Firesock: None
FS:00	Fixation, Tight Downs: Without fixation, tight downs
HA:00	Handles: None
LA:0A	User Manual Language: CZ,EN,DE,ES
MO:TD	Basic model: PrimaCare 10
PG:SP	Packing: Standard
SWLn:200	SWL: 200 kg
Wn:10	Weight: 10 kg
ZIP:01	Zip: Zip covered 180°

Compatibility

E1	Eleganza 1
EZ1	EZ1
EZ2	Essenza 200
EZ3	Essenza 300
IM3	Image 3
LA/LT	Latera Acute / Latera Thema
P1/P2	Praktika 1 / Praktika 2