

Instrument Trolley Gracie

INFORMATION

Article number **1IT-349**
Configuration state Valid article

DESCRIPTION

The instrument trolley Gracie saves space in the office and provides the doctor with maximum comfort during examinations.

TECHNICAL SPECIFICATION

Ancho	520 mm
Profundidad	600 mm
Altura	1125 mm
Peso	47 kg
Volume of front bin - right	13 l
Volume of front bin - left	18 l
Volume of the rear basket	18 l
Recommended volume of liquid in the black bin	9 l
Loading capacity of the middle shelf	15 kg
Loading capacity of the shelf for storing mirrors	3 kg

ILLUSTRATIVE IMAGE



PICTURE	DESCRIPTION	CODE
	Dieños por región: Mexico <i>Modification according to the local legislation and specific product demands.</i>	RU:MX
	Diseño básico: Mesa sin calefactar	MO:BN
	Color: Amarillo maiz F2/6586 <i>Svitap</i>	SK:M7
	Lado: Plástico	ACC:D4
	Cable y enchufe: Sin componentes eléctricos	PL:00
	Conductor design: Sin componentes eléctricos	PL[1]:00
	Color: Sin componentes eléctricos	PL[2]:00
	Voltaje: Sin componentes eléctricos	VO:00
	Idioma del manual de uso: Español	LA:05
	Accessories: Sin opción	HO:00
	Petición especial: Sin petición especial del cliente	AW:00
	Embalaje: Caja de cartón en el palet	PG:B1