

Table S3. Primer pairs used for ChIP analysis

Gene	Sequence
<i>ANGPTL2</i> E2-4	Forward: 5'-AGTGCCCAACGAGGCATCC-3' Reverse: 5'-CCAAAGCTCAGGCAGCTC-3'
<i>ANGPTL2</i> E2/3	Forward: 5'-AATGAGGCCTCCACGAAG-3' Reverse: 5'-GGGTTTAGATCCGCTCCAG-3'
<i>ANGPTL2</i> E4	Forward: 5'-ACAGCTGGAGCGGATCTAAAC-3' Reverse: 5'-GGCAGGCCCCCTCCTTCCTT-3'
<i>PER2</i> E-box	Forward: 5'-ATGTATGCAGATGAGACGGAGTCG-3' Reverse: 5'-ACAGCTGCACGTATCCCCTCAG-3'
<i>GAPDH</i>	Forward: 5'-TCGAACAGGAGGAGCAGAGAGCGA-3' Reverse: 5'-TACTAGCGGTTTTACGGGCG-3'