

Table S4. Primer sequences used for real-time PCR analysis.

Transcript ID	Forward Primer	Reverse Primer
Cr_TC35622	TGGAGTACTGGAAATGGCAAAGA	CCACGGTAGCAGGGAAATCA
Cr_TC01142	GCAAGAGCGAGCTAGGGAAGA	TTAAGCGTTCCAATTCAGGTTTGT
Cr_TC57924	GTCGAAGAATGACCAAACAAAGAA	GCTAGCGTGCGAGATGCTATT
Cr_TC33880	CCAATTTGTCATTCCTGGGATAA	CCCGAACCACTTGACGACAT
Cr_TC02224	GAAGCAAACCAGCACCAGAGT	TGATTCAAATACCTGGGCTTCCT
Cr_TC07184	ACGAAATTTGGGCAGCAGTT	CCATTTAAAAACATCCCCATTGTT
Cr_TC51001	TTCTGTTGCTGTTGCTTTTGGT	TCCGTACAATCTTGGCGAAAC
Cr_TC07063	CCAACTCAGGCCATACATTCTCA	AGACGATGACCAGTGCATTGAA
Cr_TC33950	GATCTTTCTCCTACTACTCCTCGTCAA	GTTAGCCCAAGAATTCGATTTCA
Cr_TC57680	TTGGTACAGCACGCAAAGGAT	ATTTGGATATGGAGGAGAAGTTTGG
Cr_TC35401 (Actin)	CCATCTCCAGAGTCCAGAACAAT	CCATTCAAGCTGTTTTATCCCTCTA