

# TRK Total Revision Knee System

## Product Catalog



Biotech GmbH  
Hauptstraße 113.  
56598 Rheinbrohl Germany  
Tel : +49 2635 92221-0  
Email : office-de@biotech-medical.net  
www.biotech-medical.com



# BIOTECH

REF  
VB-004-PROSP-PC-06-EN-00

Revision  
2  
Publication date :  
01.04.2022.


© 1011

## "Movement is Life"

The **Biotech TRK Revision Knee System** is a completely modular system, compatible with all Biotech Primary Knee Systems. The system provides variety of solutions covering a wide range of revision indications.

**Main features:**

- Total modularity of parts
- **Revision Femoral (A)** components, CoCr, Anatomical design – 7 sizes
- **Revision Tibial Tray (B)** components, TiAl6V4 – 7 sizes
- **Revision inserts**, UHMWPE – 7 sizes, 6 heights
- **Tibial and Femoral Augmentation Blocks**, TiAl6V4 – with screw fixation
- **Universal stems (C)** for Tibial and Femoral components, TiAl6V4
- **Adapters (D)** for stems and offset adapters, TiAl6V4 – 5°, 7°, 9°
- **360° rotatable offset adapters**, TiAl6V4 – provide **4 mm offset**. Short offset (**E**) adapter (44mm) for tibial and femoral side and long offset (**F**) adapter (75mm) for tibial side.



\* distance between the two parallel axis: implant axis and stem axis ► 4mm offset



**TRK Cemented Revision  
femoral components  
CoCr DIN EN ISO 5832-4**

**Right**

XXS	054-1311-0558
XS	054-1311-1562
S	054-1311-2566
M	054-1311-3570
L	054-1311-4574
XL	054-1311-5578
XXL	054-1311-6582

**Left**

XXS	054-1312-0558
XS	054-1312-1562
S	054-1312-2566
M	054-1312-3570
L	054-1312-4574
XL	054-1312-5578
XXL	054-1312-6582

**TRK Revision  
Universal Cemented Stems  
TiAl6V4 DIN EN ISO 5832-3**

10x50 *	054-5410-0510
10x80 *	054-5410-0810
10x110 *	054-5410-1110
10x140 *	054-5410-1410
10x170 *	054-5410-1710
12x50 *	054-5410-0512
12x80	054-5410-0812
12x110	054-5410-1112
12x140	054-5410-1412
12x170 *	054-5410-1712
14x50 *	054-5410-0514
14x80	054-5410-0814
14x110	054-5410-1114
14x140	054-5410-1414
14x170 *	054-5410-1714
16x50 *	054-5410-0516
16x80	054-5410-0816
16x110	054-5410-1116
16x140	054-5410-1416
16x170 *	054-5410-1716
18x50 *	054-5410-0518
18x80 *	054-5410-0818
18x110 *	054-5410-1118
18x140 *	054-5410-1418
18x170 *	054-5410-1718

**Stem adapters  
TiAl6V4 DIN EN ISO 5832-3**

Adapter for Revision Universal Stem 5°, 9°	054-9415-0001
Adapter for Revision Universal Stem 7°	054-9415-0002

**Stem offset adapters  
TiAl6V4 DIN EN ISO 5832-3**

Short offset adapter for Femoral and Tibial Stem	054-6415-0500
Long offset adapter for Tibial Stem	054-5415-0500

**TRK Revision Tibial Tray  
TiAl6V4 DIN EN ISO 5832-3**

XXS	054-3410-0500
XS	054-3410-1500
S	054-3410-2500
M	054-3410-3500
L	054-3410-4500
XL	054-3410-5500
XXL	054-3410-6500

**TRK Revision Tibial Inserts  
UHMWPE DIN EN ISO 5834-2**

XXS 8mm	054-4500-0508
XXS 10mm	054-4500-0510
XXS 12mm	054-4500-0512
XXS 15mm	054-4500-0515
XXS 18mm	054-4500-0518
XXS 20mm	054-4500-0520
XS 8mm	054-4500-1508
XS 10mm	054-4500-1510
XS 12mm	054-4500-1512
XS 15mm	054-4500-1515
XS 18mm	054-4500-1518
XS 20mm	054-4500-1520
S 8mm	054-4500-2508
S 10mm	054-4500-2510
S 12mm	054-4500-2512
S 15mm	054-4500-2515
S 18mm	054-4500-2518
S 20mm	054-4500-2520
M 8mm	054-4500-3508
M 10mm	054-4500-3510
M 12mm	054-4500-3512
M 15mm	054-4500-3515
M 18mm	054-4500-3518
M 20mm	054-4500-3520
L 8mm	054-4500-4508
L 10mm	054-4500-4510
L 12mm	054-4500-4512
L 15mm	054-4500-4515
L 18mm	054-4500-4518
L 20mm	054-4500-4520
XL 8mm	054-4500-5508
XL 10mm	054-4500-5510
XL 12mm	054-4500-5512
XL 15mm	054-4500-5515
XL 18mm	054-4500-5518
XL 20mm	054-4500-5520
XXL 8mm	054-4500-6508
XXL 10mm	054-4500-6510
XXL 12mm	054-4500-6512
XXL 15mm	054-4500-6515
XXL 18mm	054-4500-6518
XXL 20mm	054-4500-6520



**“Movement is Life”**



**TRK TAB**  
TiAl6V4 DIN EN ISO 5832-3

XXS 6mm 054-8400-0006  
XS 6mm 054-8400-1006  
S 6mm 054-8400-2006  
M 6mm 054-8400-3006  
L 6mm 054-8400-4006  
XL 6mm 054-8400-5006  
XXL 6mm 054-8400-6006

XXS 9mm 054-8400-0009  
XS 9mm 054-8400-1009  
S 9mm 054-8400-2009  
M 9mm 054-8400-3009  
L 9mm 054-8400-4009  
XL 9mm 054-8400-5009  
XXL 9mm 054-8400-6009

XXS 12mm 054-8400-0012  
XS 12mm 054-8400-1012  
S 12mm 054-8400-2012  
M 12mm 054-8400-3012  
L 12mm 054-8400-4012  
XL 12mm 054-8400-5012  
XXL 12mm 054-8400-6012

XXS 15mm 054-8400-0015  
XS 15mm 054-8400-1015  
S 15mm 054-8400-2015  
M 15mm 054-8400-3015  
L 15mm 054-8400-4015  
XL 15mm 054-8400-5015  
XXL 15mm 054-8400-6015

XXS 18mm 054-8400-0018  
XS 18mm 054-8400-1018  
S 18mm 054-8400-2018  
M 18mm 054-8400-3018  
L 18mm 054-8400-4018  
XL 18mm 054-8400-5018  
XXL 18mm 054-8400-6018



**TRK FAB**  
TiAl6V4 DIN EN ISO 5832-3

**4mm**

XXS 4mm LM/RL DISTAL 054-0403-0004  
XXS 4mm LL/RM DISTAL 054-0400-0204  
XXS 4mm POSTERIOR 054-0404-1104  
XS 4mm LM/RL DISTAL 054-0403-1004  
XS 4mm LL/RM DISTAL 054-0400-1204  
XS 4mm POSTERIOR 054-0404-2104  
S 4mm LM/RL DISTAL 054-0403-2004  
S 4mm LL/RM DISTAL 054-0400-2204  
S 4mm POSTERIOR 054-0404-3104  
M 4mm LM/RL DISTAL 054-0403-3004  
M 4mm LL/RM DISTAL 054-0400-3204  
M 4mm POSTERIOR 054-0404-4104  
L 4mm LM/RL DISTAL 054-0403-4004  
L 4mm LL/RM DISTAL 054-0404-4204  
L 4mm POSTERIOR 054-0404-5104  
XL 4mm LM/RL DISTAL 054-0403-5004  
XL 4mm LL/RM DISTAL 054-0400-5204  
XL 4mm POSTERIOR 054-0404-6104  
XXL 4mm LM/RL DISTAL 054-0403-6004  
XXL 4mm LL/RM DISTAL 054-0400-6204  
XXL 4mm POSTERIOR 054-0404-7104

**8mm**

XXS 8mm LM/RL DISTAL 054-0403-0008  
XXS 8mm LL/RM DISTAL 054-0400-0208  
XXS 8mm POSTERIOR 054-0404-1108  
XS 8mm LM/RL DISTAL 054-0403-1008  
XS 8mm LL/RM DISTAL 054-0400-1208  
XS 8mm POSTERIOR 054-0404-2108  
S 8mm LM/RL DISTAL 054-0403-2008  
S 8mm LL/RM DISTAL 054-0400-2208  
S 8mm POSTERIOR 054-0404-3108  
M 8mm LM/RL DISTAL 054-0403-3008  
M 8mm LL/RM DISTAL 054-0400-3208  
M 8mm POSTERIOR 054-0404-4108  
L 8mm LM/RL DISTAL 054-0403-4008  
L 8mm LL/RM DISTAL 054-0404-4208  
L 8mm POSTERIOR 054-0404-5108  
XL 8mm LM/RL DISTAL 054-0403-5008  
XL 8mm LL/RM DISTAL 054-0400-5208  
XL 8mm POSTERIOR 054-0404-6108  
XXL 8mm LM/RL DISTAL 054-0403-6008  
XXL 8mm LL/RM DISTAL 054-0400-6208  
XXL 8mm POSTERIOR 054-0404-7108

**12mm**

XXS 12mm LM/RL DISTAL 054-0403-0012  
XXS 12mm LL/RM DISTAL 054-0400-0212  
XXS 12mm POSTERIOR 054-0404-1112  
XS 12mm LM/RL DISTAL 054-0403-1012  
XS 12mm LL/RM DISTAL 054-0400-1212  
XS 12mm POSTERIOR 054-0404-2112  
S 12mm LM/RL DISTAL 054-0403-2012  
S 12mm LL/RM DISTAL 054-0400-2212  
S 12mm POSTERIOR 054-0404-3112  
M 12mm LM/RL DISTAL 054-0403-3012  
M 12mm LL/RM DISTAL 054-0400-3212  
M 12mm POSTERIOR 054-0404-4112  
L 12mm LM/RL DISTAL 054-0403-4012  
L 12mm LL/RM DISTAL 054-0400-4212  
L 12mm POSTERIOR 054-0404-5112  
XL 12mm LM/RL DISTAL 054-0403-5012  
XL 12mm LL/RM DISTAL 054-0400-5212  
XL 12mm POSTERIOR 054-0404-6112  
XXL 12mm LM/RL DISTAL 054-0403-6012  
XXL 12mm LL/RM DISTAL 054-0400-6212  
XXL 12mm POSTERIOR 054-0404-7112



\* - optional

**“Movement is Life”**



REF  
VB-004-PROSP-PC-06-EN-00

Revision  
2

Publication date :  
01.04.2022.

©1011



Biotech GmbH  
Hauptstraße 113.  
56598 Rheinbrohl Germany  
Tel : +49 2635 92221-0  
Email : office-de@biotech-medical.net  
www.biotech-medical.com