

CA

.....

Bl

60

RÖS'S


Binar



Aesthetic technology division

Advanced professional technology combining diathermy and electro-stimulation. Beauty and wellness.

6



**Restores cell pulse Recovers** natural beauty



## Regenerates with non-invasive energy

Binary PREMIUM TISSUE REGENERATION AND MUSCLE TONE

# The new multifunctional device that interacts with professionals

## RÖS'S Estética introduces its new BINARY PREMIUM model, the advanced multifunctional diathermy device combined with electro-stimulation.

**Binary Premium** is an innovative technological model that provides a full response to the most common beauty concerns in the cabin. With an elegant design and extended benefits, Binary strengthens its position as an absolute reference for the treatments of tissue regeneration, muscle toning, body sculpting and recovery of wellness.

RÖSA



**ANTI-FALL CLICK-MAGNETIC SYSTEM** in all its electrodes. Accident prevention is a priority.





**ADDITION OF NEW** ACCESSORIES to extend the application options. Indirect application wristband. Optional range of concave and endocavitary applicators.

RÖS'S CAPACITIVE **TOUCH SCREEN 10.1'** Fast and sensitive to touch.

BI

**EVERYTHING IN ITS PLACE** Designed thinking about

the location of everything you need in your daily life.



**MOBILITY AND STABILITY** Because we know how dynamic cabin life is, the design adapts to you.

## No one offers you so much

**RÖS'S** keeps the **Binary** characteristics that have made it a successful system and extends its benefits improving the response time of the software, improving the design and increasing the equipment.



PULSED-ENERGY SYSTEM

delivery in real time based

on the response of each

client's or patient's tissues.

| ÆJ |  |  |
|----|--|--|

NON

FREE NON-STOP KIT

with spare parts in handles It ensures consistent energy and return plate wire. It will protect you against any incident.



We meet your needs. If you have several cabins, the all-in-one model meets your needs. If instead you have no space, we propose to you the desktop model.

We reinvent multifunctional technology for beauty treatments. We search for excellence in outcomes and integral wellness.




## ERGO-RETURN PLATE

New return plate. Customisable and skin-friendly.



## WE EXTEND THE SOFTWARE CAPACITY 30 free channels to record your programs and more programming options for experts and physical therapists.



## **Action from** the heart of the cell

## **DIATHERMY 448. Science applied to regeneration** and recovery.

Any unsightly imperfection or disorder results in a change in the cell electrical potential. The consequences triggered by this affect our morphology in a highly diverse manner, so we need a precise energy that activates or slows down cell pulse. This rebalancing and regenerating global action recovers the damaged body tissues progressively.

**Binary Premium emits an adequate high-frequency** signal of 448 KHz, scientifically supported to achieve intercellular ion mobilisation and action from the heart of the cell. Its modes of application, capacity monopolar and resistive monopolar, allow to reach all tissues, from the dermis to the bone tissue.

## 448 KHz frecuency

RÖS'S

## **ELECTRO-STIMULATION.** The most efficient electrotherapy.

**Binary** treatment is complemented with electro-stimulation to achieve an optimum result, recovering muscle tone and reducing fatty tissue.

## FACIAL ELECTROLIFTING

- Activation of microcirculation and general metabolism with a revitalising effect for tissues.
- Tone normalisation. Improves muscle quality affected by the gravitational effect that occasionally causes more marked fine lines.



Pointer electrodes facilitate precision and control in each application manoeuvre according to the muscle being worked on.

RÖS'S

## ELECTRO-LIPOLYSIS

The lipolytic current designed for **Binary** is intended to stimulate the local level sympathetic system and trigger a triglyceride degradation process, while it activates the circulatory system to facilitate the removal of fat content. This is possible thanks to the hands-free electro-stimulation bands that synchronize two processes: the degradation of fat and its drainage in an automatic way.

## ELECTRO-FIRMING

The **Binary** excitomotor current has a firming function. Its design stimulates sequential muscle contractions, which makes them ideal in toning and muscle recovery treatments.



The electro-stimulation bands distribute the current evenly on the whole muscle group. For specific work on specific muscles the connections system also allows us to use individual electrodes.

## **Binary** PREMIUM **TISSUE REGENERATION** AND MUSCLE TONE

## **5 techniques in the same** device that together complement and strengthen each other

The medical-aesthetic professional community has studied and reached the conclusion that to deal effectively with the most sought after facial and body treatments in booths various agents need to be used, which work together well and enhance their effects.

The in sequence application of the capacitive and resistive diathermy and the electro-stimulation is a methodological commitment able to globally treat all causes which create problems such as ageing skin, dermal laxity, loss of muscle tone, cellulite processes or excess of volume.

## **Binary Method**

| 1 | Capacitive<br>biostimulation | More specific action in dermis and muscles.                            |
|---|------------------------------|--|
|   |                              | Active hyperthermia increases the production of collagen and elastine. |
| 2 | Resistive biostimulation     | More specific action on fat and dense tissues.                         |
|   |                              | Recovery treatments and reducing fatty tissue.                         |
|   |                              | ••••••   |
| 3 | Electrolifting               | Specific action on elongated or<br>retracted face muscles.             |
|   |                              | Recovery of muscle tone.   |
|   |                              |  |
|   | Electro-lipolysis            | Lipolytic action.  |
|   | Electro-<br>stimulation      | Firming action.  |
|   |                              |  |

## PROGRAMS

Binary Premium introduces three programming modes. The pre-established programs are supported by the protocol guide for enhancing each stage step by step. In addition, the device has a free programming mode to draw up customized programs and free mode for health and wellness professionals.

| Predefined<br>programs    |   |  |
|---------------------------|---|--|
|                           |   |  |
|                           | FACIAL PROGRAMS                                       | BODY PROGRA  |
|                           | 1. Firming<br>2. Revitalising                         | 1. Anti-cellul<br>2. Anti-cellul                                 |
|                           | 3. Upper intensive                                    | 3. Anti-cellul   |
|                           | 4. Middle intensive                                   | 4. Reduction.  |
|                           | 5. Lower intensive                                    | 5. Firming   |
|                           | Guide of protocols and tutorials                      | available.   |
|                           |   |  |
|                           |   |  |
| Free .                    | EXPERT MODE WITH 30 CHANN                             | ELS AVAILABLE FOR  |
| programming               |   |  |
|                           |   |  |
|                           | FACE  | BODY   |
|                           | Under-eye bags<br>Double chin<br>Couperose<br>Rosacea | Hair treatmer<br>Breast firmin<br>Ongoing stret<br>Recovery of s |
|                           | Oil regulation  | Marked laxity<br>Keloids<br>Recovery of v                        |
|                           |   |  |
| Individual<br>techniques* | MEDICAL AND PHYSIO. MODE                              | SPECIFIC TECHNOL   |
|                           |   |  |
|                           |   |  |
|                           | AESTHETIC MEDICINE                                    | GYNAECO-AES  |
|                           | Anti-aging  | Recovery of p  |
|                           | Healing   | Vaginal dysfu  |
|                           | Recovery from skin injuries                           |  |
|                           | Pre and post-operative periods                        | 5  |
|                           | Treatments for pain                                   |  |
|                           |   |  |
|                           | *Binary Premium device allows                         |  |
|                           | professionals. If your company h                      | ias a license ior med  |

ness opportunities in collaboration with other health and wellness nedical applications, this is your opportunity.



### AMS

lite. Soft cellulite lite. Hard cellulite lite. Oedematous cellulit n. Localised adiposity



ROS

OR RECORDING YOUR OWN PROGRAMMES.

etch marks scars v (abdomen, arms)

wellness



### STHETICS

pelvic floor unctions

### PHYSICALTHERAPY

Rehabilitation Sports physical therapy Recovery of pelvic floor Physio-aesthetics

## Protocols tested under medical supervision

The preliminary case study, performed under medical supervision, evidences the efficiency of the **Binary system** in the classic cabin treatments and specific treatments for recovery from skin injuries.

## **Body treatments**

## Facial Treatments

- Mobilisation of localised fat and body sculpting.
- Correction of flaccidity and irregularities in the skin surface.
- Very significant improvement of muscle tone.
- Acceleration of skin recovery in damaged tissues (keloids, scars, burns or ongoing stretch marks).
- Breast firming.

## Revitalization of skin and a

- significant improvement in tone and luminosity.
- Improvement of elasticity and firmness in the skin around the eyes, neckline and double chin.
- There is a relaxation of muscle contractures and an improvement in abnormally relaxed muscle tone.
- Improve of rosacea and couperose.


Skin firming and muscle toning in the abdomen.



3

Pro-aging revitalisation with specific work on fine lines.

Remarkable case histories





Recovery from skin injuries. Burn from chemical peeling.

5

**RÖS'S** 

BINARY

Binar

Recovery from skin injuries. Surgical scar.

Before treatment
























After treatment

# Recognition and commercial success

**RÖS'S** offers you a full commercial support pack to promote your treatments. Welcome to the BINARY Community, our team is at your disposal.



## Technical specifications and range of accessories

| Code    | Accessory                        |
|---------|----------------------------------|
| 04BI000 | ELECTROLIFTING APPLICATOR        |
| 33B0000 | ELECTROLIFTING CABLE             |
| 04BI060 | STIMULATION BAND 58CM            |
| 04BI085 | STIMULATION BAND 74CM            |
| 33R001A | STIMULATION CABLE A              |
| 33R001B | STIMULATION CABLE B              |
| 33R001C | STIMULATION CABLEC               |
| 04BI140 | ELASTIC BAND 40CM                |
| 04BI160 | ELASTIC BAND 60CM                |
| 04BI100 | ELASTIC BAND 100CM               |
| 04BPRF2 | CAPACITIVE MAGNETIZED ANGLED HA  |
| 04BPRF4 | RESISTIVE MAGNETIZED ANGLED HAN  |
| 04BPC30 | CAPACITIVE MAGNETIZED ELECTRODE  |
| 04BPC50 | CAPACITIVE MAGNETIZED ELECTRODE  |
| 04BPC70 | CAPACITIVE MAGNETIZED ELECTRODE  |
| 04BPR30 | RESISTIVE MAGNETIZED ELECTRODE   |
| 04BPR50 | RESISTIVE MAGNETIZED ELECTRODE   |
| 04BPR70 | RESISTIVE MAGNETIZED ELECTRODE   |
| 04B0000 | RETURN PLATE                     |
| 33BP030 | RETURN PLATE CABLE               |
| 04R0010 | CAPACITIVE SCREEN PEN            |
| 04BIRP0 | INDIRECT APPLICATION WRISTBAND   |
| 04BP000 | ELASTIC BAND FOR WRISTBAND       |
| 33BP040 | WRISTBAND CABLE                  |
| 99B0000 | GLASS BOWL 0.07L 7.5X3.2CM.      |
| 0401330 | SELF-ADHESIVE ELECTRODES         |
| 6600200 | BIENERGY EMULSION 1.000ML        |
| 6600210 | BIENERGY ELECTRO-STIMULATION GEL |
|         |                                  |

## **KIT NON-STOP**

| 3BP030 | RETURN PLATE CABLE              |
|--------|---------------------------------|
| 4BPRF2 | CAPACITIVE MAGNETIZED ANGLED HA |
| 4BPRF4 | RESISTIVE MAGNETIZED ANGLED HAI |

