Integral treatment system Do not rule out any treatment options

Advantages of using two techniques simultaneously

- I Electro-stimulation is a very pleasant sensation
- I Shorter application times
- I Increased electrical conductivity from sweating
- I Benefits of both static and dynamic heat

Benefits derived from TEI System

- I Lipolytic effect from penetration of infrared heat from 1 cm to 3.5 cm
- I Stimulation of blood and lymphatic circulation
- I Increased metabolic activity
- I Improved vascular system and venous return
- I Increased cell tropism, healing capacity and connective tissue elasticity
- I Analgesia and muscle relaxation
- I Muscle firming with all treatments



TEI SYSTEM . REF.	0080000	
DIMENSIONS	Width Height Depth	490 mm 1.000 mm 490 mm
WEIGHT POWER CONSUMPTION		34 kg 290 W

- 2 580×190 mm 35 W thermal bands COD. 24TE235 4 240×260 mm 25 W thermal bands COD. 24TE125
- 4 240×170 mm 15 W thermal bands COD. 24TE115
- 2 $180 \times 67 \text{ mm 4 W thermal bands COD. } 24\text{TE}104$
- 10 Band cables labelled from A to J. COD. 33TE114 4 100 cm elastic bands COD. 0400800
- 6 60 cm elastic bands COD, 0401000
- 6 40 cm elastic bands COD. 0401100
- 1 Band hanger COD. 30T0110
- 1 Firming ElectroGel (1000ml) 6600071 1 Draining ElectroGel (1000ml) 6600070

- 1 580×190 mm 35 W thermal band COD. 24TE235
- 1 240×260 mm 25 W thermal band COD. 24TE125
- 1 240×170 mm 15 W thermal band COD. 24TE115
- 1 180 × 67mm 4 W thermal band COD, 24TE104
- 1 Band cable 33TE114

Tei System belongs to SUMMA range of RÖS'S Profitable treatments for the professional beauticians





TEI System Integral thermo-stimulation

Infrared band and variable current equipment for joint heat and muscle stimulation treatment.

RÖS'S pioneered the development of thermostimulation and its team of engineers, in collaboration with beauty professionals, has designed the most sophisticated application bands to provide customers with the best possible results and maximum safety.



INTEGRAL TREATMENT FOR THE ABDOMEN, THIGHS, CALF MUSCLES AND FEET

Integral treatments



OVERALL ACTION FOR A WELL-BALANCED PROFESSIONAL RESULT

The pleasure of technology Reduction of volume in just four sessions

How does TEI System work?

Two basic professional beauty-treatment agents, HEAT and STIMULATION, are the key factors in TEI System. Thanks to deep infrared heat, the adipose cells in fat tissue are heated, aiding their liquefaction by electric current. The controlled emission of heat combined with different types of electric current guarantee the results.

The benefits of thermo-stimulation have been medically proven by Dr Juan Ramón Zaragoza.

All-round TEI treatments. Your appointments diary is full!

18 electro-currants combined in 21 body programmes Soft cellulite

Hard cellulite

Oedematous cellulite

Firming

Relaxation

Intense firming

I TENS

Post-pregnancy Excess weight

| Lipolytic Passive

Localised excess

gymnastics

Calming passive

weight

gymnastics Sequential

Circulatory

Draining

boosting | Lipo-Draining

Reduction Sculpting

Upper trunk

Toning



The easily applied thermo-stimulation strips are the key to **TEI treatment.** Flexible, adaptable, tough and safe, their laminar circuit provides completely uniform infrared emission, while their flexible cover adapts to the body to ensure full contact.

Supplied with 12 long-lasting strips which adapt to the body. Rapid application. Easy to clean. Four different sizes. There is also a Courtesy Kit to cope with any incident.

