

# Skin Analysis






- 1.Camera
- 2.General power switch
- 3.Power supply line socket
- 4.USB socket
- 5.Moisture test pen Start button
- 6.Camera flash



Moisture test



visible light Sketch



Ultraviolet Sketch



Clear the screen view



Hide and show test value



Open images



Save images



Show outline of face



Shows the balance of face by contour



Eidt



Move to result page



Full screen with zoom function



Show 3D view of skin



Compare two faces by overwrapping



Print out images



Simulation



Delete the diagnosis date list