



Kerrison Catalog



Table of Contents

Kerrison Portfolio Overview	3-4
Symmetry Sharp Kerrisons®	5-7
• Reusable Handles	7
• Single-Use & Reusable Tips	7
RapidClean® Convertible Kerrisons	8-10
• Ceramic Coated	9
• Stainless Steel	10
RapidClean® Detachable Kerrisons	11-12
• Ceramic Coated	12
One-Piece Kerrisons	13-22
• Ceramic Coated	14-15
• Ceramic Coated - Opti-Length®	16
• Classic® Stainless Steel	17
• Stainless Steel	18-19
• Stainless Steel - Extended Length	20
• Stainless Steel - Spurling Handle	21
• Ultra Coat - Titanium Nitride Coating	22
Specialty Kerrisons	23-24
• Bayonet-Shaped	23
• Curved-Up	23
• Laminectomy	24
• Schlesinger Cervical Laminectomy	24
Customer Support	25

Kerrison Portfolio

Overview


Symmetry Surgical has a rich heritage in Kerrison design dating back to the early 1900s. Collaborations between Thomas Codman and Harvey Cushing M.D. yielded unique insight into the needs of surgeons. Dr. Cushing was one of the first surgeons to specialize in neurosurgery and his initial input and designs are evident in modern day Kerrisons.

Our Kerrison designs, materials and performance have resulted in products that have stood the test of time. Our Kerrison portfolio has a pedigree of breadth and quality.



Kerrison Portfolio

Kerrison Anatomy



Symmetry Sharp Kerrisons®

What If...


Your biggest Kerrison challenges could be solved with one product?



Single-Use or Reusable Tip



- You never had to use a dull kerrison again?
- You could have a productive bite every time?
- You could minimize the risk of dural tears?




- You could reduce hand and arm fatigue?
- You could easily open the kerrison for cleaning?
- You no longer had to open multiple trays?

Introducing the
RapidClean
Symmetry Sharp Kerrison®

Symmetry Sharp Kerrisons®

Features & Benefits



“I choose to use the Symmetry Sharp Kerrisons for my ACD cases because they reduce the decompression time and improve the efficiency of the procedure.”

Neurosurgeon, Louisiana

“If one has to tug and twist the Kerrison, the risk of a dural tear is much higher.”

Neurosurgeon, Virginia

Symmetry Sharp Kerrisons®

Reusable Handles

Ceramic Coated



DESCRIPTION

RAPIDCLEAN®

ONE-PIECE

5" (125 mm) Handle for 7" (178 mm) Kerrison	53-1674RC	53-1674
6" (152 mm) Handle for 8" (203 mm) Kerrison	53-1675RC	53-1675
7" (178 mm) Handle for 9" (229 mm) Kerrison	53-1676RC	53-1676

Single-Use & Reusable Tips

40° Up, 9 mm Opening, 2" (50 mm) Tip Length

Single-Use packaged sterile, 3 per box

Reusable packaged non-sterile, 2 per pack



DESCRIPTION

SINGLE-USE THIN FOOTPLATE


SINGLE-USE STANDARD FOOTPLATE

REUSABLE STANDARD FOOTPLATE

1 mm Tip	53-1661T	53-1661	
2 mm Tip	53-1662T	53-1662	53-1662R
3 mm Tip	53-1663T	53-1663	53-1663R
4 mm Tip	53-1664T	53-1664	53-1664R

RapidClean® Convertible Kerrisons

Our RapidClean® Convertible Kerrisons are designed to allow for thorough cleaning and decontamination. These Kerrisons offer a convenient hinge design that opens easily for cleaning, inspection and lubrication while eliminating the risk of misplacing pieces during processing.



RapidClean[®] Convertible Kerrisons

Ceramic Coated

- Smooth, wear-resistant surface
- Reduces glare
- Optimal gliding
- Stays sharp longer



59-3200J-BD

		Bite (mm)	Opening (mm)	7" (178 mm)	8" (203 mm)	9" (229 mm)	
		Thin Footplate					
40° Up	Without Ejector	1	9	58-1180-BD	58-1200-BD	58-1230-BD	
		2	9	58-2180-BD	58-2200-BD	58-2230-BD	
		3	9	58-3180-BD	58-3200-BD	58-3230-BD	
		4	9	58-4180-BD	58-4200-BD	58-4230-BD	
		5	9	58-5180-BD	58-5200-BD	58-5230-BD	
		6	9			58-6230-BD	
		Standard Footplate					
		1	9	59-1180-BD	59-1200-BD	59-1230-BD	
		2	9	59-2180-BD	59-2200-BD	59-2230-BD	
		3	9	59-3180-BD	59-3200-BD	59-3230-BD	
		4	9	59-4180-BD	59-4200-BD	59-4230-BD	
		5	9	59-5180-BD	59-5200-BD	59-5230-BD	
	2	15		58-2262-BD			
	3	15		58-2263-BD			
	4	15		58-2264-BD			
	5	15		58-2265-BD			
	With Ejector	Thin Footplate					
		2	9		58-23200J-BD		
		3	9	58-3180J-BD	58-3200J-BD	58-3230J-BD	
4		9	58-4180J-BD	58-4200J-BD	58-4230J-BD		
5		9	58-5180J-BD	58-5200J-BD	58-5230J-BD		
Standard Footplate							
3		9		59-3200J-BD			
4		9		59-4200J-BD			
5		9		59-5200J-BD			

RapidClean® Convertible Kerrisons

Stainless Steel

- Durable
- Economical



58-1180

		Bite (mm)	Opening (mm)	7" (178 mm)	8" (203 mm)	9" (229 mm)	
		Thin Footplate					
40° Up	Without Ejector	1	9	58-1180	58-1200		
		2	9	58-2180	58-2200		
		3	9	58-3180	58-3200		
		4	9		58-4200		
		Standard Footplate					
		1	9	59-1180	59-1200	59-1230	
		2	9	59-2180	59-2200	59-2230	
		3	9	59-3180	59-3200	59-3230	
		4	9		59-4200	59-4230	
		5	9		59-5200	59-5230	
	2	15		58-2262			
	3	15		58-2263			
	4	15		58-2264			
	5	15		58-2265			
	With Ejector	Thin Footplate					
3		9		58-3200J			
4		9	58-4180J	58-4200J			
5		9		58-5200J			

RapidClean® Detachable Kerrisons

Symmetry's RapidClean® Detachable Kerrisons are easily disassembled into two pieces for cleaning and inspection. They feature a frictionless channel, secondary safety locks to ensure secure latching and unique markings on the top and bottom sections to ensure straight-forward reassembly.



Secondary safety lock to ensure latch is secure

Unique markings on top and bottom segments to ensure proper match for reassembly

RapidClean® Detachable Kerrisons

Ceramic Coated

- Smooth, wear-resistant surface
- Reduces glare
- Optimal gliding
- Stays sharp longer



53-1572CD

40° Up	Without Ejector	Bite (mm)	Opening (mm)	7" (178 mm)	8" (203 mm)	9" (229 mm)
		Thin Footplate				
		1	10	53-1540CD	53-1570CD	53-1548CD
		2	10	53-1542CD	53-1572CD	53-1550CD
		3	10	53-1544CD	53-1574CD	
		4	10	53-1546CD	53-1576CD	
		Standard Footplate				
	With Ejector	2	10		53-1624LD	53-1642LD
		2	12	53-1612LD		
		3	10	53-1615LED	53-1627LED	53-1645LED
		4	10	53-1618LED	53-1630LED	53-1648LED
		5	10	53-1621LED	53-1633LED	53-1651LED
		6	10		53-1636LED	53-1654LED

One-Piece Kerrisons

Symmetry's one-piece Kerrisons have the standard Kerrison design that you have known and trusted for years. We offer a variety of bite sizes, lengths and finishes for any preference or procedure.



One-Piece Kerrisons

Ceramic Coated

- Smooth, wear-resistant surface
- Reduces glare
- Optimal gliding
- Stays sharp longer



57-0761-BD

		Bite (mm)	Opening (mm)	7" (178 mm)	8" (203 mm)	9" (229 mm)	
		Thin Footplate					
40° Up	Without Ejector	1	9	57-0721-BD	57-0761-BD		
		2	9	57-0722-BD	57-0762-BD		
		3	9	57-0723-BD	57-0763-BD		
		4	9	57-0724-BD	57-0764-BD		
		5	9	57-0725-BD	57-0765-BD		
		6	9		57-0766-BD		
		Standard Footplate					
		1	9			57-0891-BD	
		2	9			57-0892-BD	
		3	9			57-0893-BD	
		4	9			57-0894-BD	
		5	9			57-0895-BD	
	6	9			57-0896-BD		
	2	15		57-0022-BD			
	3	15		57-0023-BD			
	4	15		57-0024-BD			
	5	15		57-0025-BD			
	With Ejector	Standard Footplate					
2		9	57-0392-BD	57-0692-BD			
3		9	57-0393-BD	57-0693-BD	57-0893J-BD		
4		9	57-0394-BD	57-0694-BD	57-0894J-BD		
5		9	57-0395-BD	57-0695-BD			
6		9	57-0396-BD	57-0696-BD			
40° Down	Without Ejector	Thin Footplate					
		1	9		57-0801-BD		
		2	9		57-0802-BD		
		3	9		57-0803-BD		
		4	9		57-0804-BD		

One-Piece Kerrisons

Ceramic Coated

- Smooth, wear-resistant surface
- Reduces glare
- Optimal gliding
- Stays sharp longer



57-0740-BD

		Bite (mm)	Opening (mm)	7" (178 mm)	8" (203 mm)	9" (229 mm)
		Thin Footplate				
90° Up	Without Ejector	1	9	57-0701-BD	57-0739-BD	
		2	9	57-0702-BD	57-0740-BD	
		3	9	57-0703-BD		
		4	9	57-0704-BD		
		5	9	57-0705-BD		
Thin Footplate						
90° Down	Without Ejector	1	9		57-0751-BD	
		2	9		57-0750-BD	
		3	9		57-0752-BD	
		4	9		57-0754-BD	

One-Piece Kerrisons

Ceramic Coated - Opti-Length®

- Smooth, wear-resistant surface
- Reduces glare
- Optimal gliding
- Stays sharp longer
- Extended Length



67-7267

		Bite (mm)	Opening (mm)	10" (254 mm)	12" (305 mm)	15" (381 mm)
		40° Up	Without Ejector	Standard Footplate*		
1	10					67-7325
2	10			67-7267	67-7282	67-7326
3	10			67-7270	67-7285	67-7327
4	10			67-7273	67-7288	67-7328
5	10			67-7276		67-7329
90° Up	Without Ejector	Standard Footplate*				
		1	10			67-7330
		2	10			67-7331
		3	10			67-7332
		4	10			67-7333
		5	10			67-7334

*Includes drainage holes for flushing

Classic® One-Piece Kerrisons

Classic® Stainless Steel Thin Footplate

- Gold screw for easy identification of thin footplate designed to easily fit under the cervical lamina



53-1506

		Bite (mm)	Opening (mm)	6" (152 mm)	7" (178 mm)	8" (203 mm)
40° Up	Without Ejector	Cervical Footplate*				
		0.6	9		53-1536	53-1538
		1	9		53-1537	53-1539
		2	9	53-1500C	53-1512	53-1524
		3	9	53-1503	53-1515	53-1527
		4	9	53-1506	53-1518	53-1530
		5	9	53-1509	53-1521	53-1533
90° Up	Without Ejector	Cervical Footplate*				
		2	9	53-1501	53-1513	53-1525
		3	9	53-1504	53-1516	53-1528
		4	9	53-1507	53-1519	53-1531
		5	9	53-1510	53-1522	53-1534
90° Down	Without Ejector	Cervical Footplate*				
		2	9	53-1502	53-1514	53-1526
		3	9	53-1505	53-1517	53-1529
		4	9	53-1508	53-1520	53-1532
		5	9	53-1511	53-1523	53-1535

*Includes drainage holes for flushing

One-Piece Kerrisons

Stainless Steel

- Durable
- Economical



57-0892

		Bite (mm)	Opening (mm)	7" (178 mm)	8" (203 mm)	9" (229 mm)	11" (280 mm)
		Thin Footplate					
40° Up	Without Ejector	1	9	57-0721	57-0761	57-0891	
		2	9	57-0722	57-0762	57-0892	
		3	9	57-0723	57-0763	57-0893	
		4	9	57-0724	57-0764	57-0894	
		5	9	57-0725	57-0765	57-0895	
		6	9	57-0726	57-0766	57-0896	
	Standard Footplate*						
		2	10	53-1400	53-1480	53-1482	
		3	10	53-1402	53-1414	53-1484	53-1426
		4	10	53-1405	53-1417	53-1485	53-1427
		5	10	53-1408	53-1420C	53-1486	53-1428
		6	10	53-1411	53-1423	53-1487	53-1429
With Ejector	Standard Footplate						
	2	9	57-0392	57-0692			
	3	9	57-0393	57-0693			
	4	9	57-0394	57-0694			
	5	9	57-0395	57-0695			

*Includes drainage holes for flushing

One-Piece Kerrisons

Stainless Steel

- Durable
- Economical



57-0701

		Bite (mm)	Opening (mm)	7" (178 mm)	8" (203 mm)	9" (229 mm)	
		Thin Footplate					
90° Up	Without Ejector	1	9	57-0701	57-0739	57-0901	
		2	9	57-0702	57-0740	57-0902	
		3	9	57-0703	57-0743	57-0903	
		4	9	57-0704	57-0744	57-0904	
		5	9	57-0705	57-0745	57-0905	
		Standard Footplate*					
		2	10		53-1481		
		3	10	53-1403	53-1415		
		4	10	53-1406	53-1418		
		5	10	53-1409	53-1421		
6	10		53-1424				
90° Down	Without Ejector	Thin Footplate					
		2	9	57-0712	57-0750		
		3	9	57-0713	57-0752		
		4	9		57-0754		
		5	9		57-0755		
		6	9		57-0756		
		Standard Footplate*					
		3	10	53-1404	53-1416		
		4	10	53-1407	53-1419		
		5	10	53-1410	53-1422		
6	10		53-1425				

*Includes drainage holes for flushing

One-Piece Kerrisons

Stainless Steel - Extended Length

- Durable
- Economical



57-0991

		Bite (mm)	Opening (mm)	10" (254 mm)	11" (280 mm)	12" (305 mm)	
		Thin Footplate					
40° Up	Without Ejector	1	9	57-0101	57-0991		
		2	9	57-0102	57-0992		
		3	9	57-0103	57-0993		
		4	9	57-0104	57-0994		
		5	9	57-0105	57-0995		
		6	9		57-0996		
		Standard Footplate					
		2	9			57-0112	
		3	9			57-0113	
		4	9			57-0114	
		5	9			57-0115	
90° Up	Without Ejector	Thin Footplate					
		1	9	57-0911			
		2	9	57-0912			
		3	9	57-0913			

One-Piece Kerrisons

Stainless Steel - Spurling Handle

- Durable
- Economical



57-0582

		Bite (mm)	Opening (mm)	7" (178 mm)	8" (203 mm)	9" (229 mm)	11" (280 mm)
40° Up	Without Ejektor	Thin Footplate					
		2	15	57-0552	57-0582	57-0782	57-0592
		3	15	57-0553	57-0583	57-0783	57-0593
		4	15	57-0554	57-0584	57-0784	57-0594
		5	15	57-0555	57-0585	57-0785	57-0595
		6	15	57-0556	57-0586	57-0786	57-0596
		90° Up	Without Ejektor	Thin Footplate			
3	15			57-0533	57-0563		
4	15			57-0534	57-0564		
5	15			57-0535	57-0565		
6	15			57-0536	57-0566		
90° Down	Without Ejektor	Thin Footplate					
		2	15		57-0572		
		3	15	57-0543	57-0573		
		4	15	57-0544	57-0574		
		5	15	57-0545	57-0575		

One-Piece Kerrisons

Ultra Coat - Titanium Nitride Coating

- Increases strength of metal
- Stays sharp longer than stainless steel



57-0761UC

40° Up	Without Ejector	Bite (mm)	Opening (mm)	8" (203 mm)	9" (229 mm)	11" (280 mm)
		Thin Footplate				
		1	9	57-0761UC	57-0891UC	
		2	9	57-0762UC	57-0892UC	57-0992UC
		3	9	57-0763UC	57-0893UC	57-0993UC
		4	9	57-0764UC	57-0894UC	57-0994UC
		5	9	57-0765UC	57-0895UC	57-0995UC

Specialty Kerrisons

Bayonet

- One-piece kerrisons
- Ceramic coated
- Tips are offset from the handle's axis



56-1002-BD

40° Up	Without Ejector	Bite (mm)	Opening (mm)	Working Length 6.25" (160 mm)
		Standard Footplate		
1		9	56-1001-BD	
2		9	56-1002-BD	
3		9	56-1003-BD	
4		9	56-1004-BD	
5		9	56-1005-BD	
90° Up	Standard Footplate			
	1	9	56-1201-BD	
	2	9	56-1202-BD	
	3	9	56-1203-BD	
	4	9	56-1204-BD	
	5	9	56-1205-BD	

Curved-Up

- One-piece kerrisons
- Stainless steel



57-1100

70° Up	Without Ejector	Bite (mm)	Opening (mm)	8" (203 mm)
		Thin Footplate		
2		9	57-1100	
3		9	57-1102	
4		9	57-1104	
5	9	57-1105		

Specialty Kerrisons

Laminectomy

- One-piece kerrisons
- Stainless steel



53-1385

40° Up	Without Ejector	Bite (mm)	Opening (mm)	6" (152 mm)	7.5" (192 mm)	
		Standard Footplate				
3		12.5	53-1385	53-1387		
5		12.5	53-1386	53-1388		
90° Up		Standard Footplate				
		3	12.5	53-1370	53-1372	
5		12.5	53-1371	53-1373		
90° Down		Standard Footplate				
		3	12.5	53-1374	53-1376	
	5	12.5	53-1375	53-1377		

Schlesinger Cervical Laminectomy

- One-piece kerrisons
- Stainless steel



53-1380

90° Up	With Ejector	Bite (mm)	Opening (mm)	6" (152 mm)	
		Standard Footplate			
3		12.5	53-1380		
5		12.5	53-1381		
90° Down		Standard Footplate			
		3	12.5	53-1382	
5	12.5	53-1383			

Customer Support

Our highly-trained customer support team is available to take your calls from 7:00 am to 5:30 pm (CST), Monday - Friday. This group is dedicated to providing a full range of support services, including:

- Order Placement
- Order Inquiries
- Delivery Information
- Product Information
- Pricing Inquiries

If we can be of assistance, please contact us using the methods below:

PHONE

US & Canada: +1 800-251-3000

International: +1 615-964-5532

EMAIL

customerservice@symmetrysurgical.com

orders@symmetrysurgical.com

WEB

www.symmetrysurgical.com

-Live Chat on symmetrysurgical.com:

Click the “Can We Help You” window in the bottom corner of any page on symmetrysurgical.com to be connected directly to a member of our sales or customer service team. We can chat in real time and answer your questions about our company and surgical instrument portfolio.

(Available Monday through Friday, 7:30 am - 5:30 pm CST)

MAILING ADDRESS

Symmetry Surgical Inc.

Attn: Customer Service

3034 Owen Drive

Antioch, TN 37013



3034 Owen Drive | Antioch, TN 37013 | 800.251.3000 (US) | +41 (0) 52 557 59 00 | customerservice@symmetrysurgical.com | symmetrysurgical.com

©2018 Symmetry Surgical Inc. All rights reserved. INTERNATIONAL KERRISON CAT 09.30.2020.V1