


SpectralWorkbench.org

Instructions and usage at
PublicLab.org/wiki/paper-spec

This spectrometer can be mounted on
 a camera phone, laptop, or webcam.

Line up carefully so that the rainbow
 is in the middle of the image, and
 tape down firmly so that the
 spectrometer stays rigid.



Open hardware by Public Lab contributors;
 Copy, share, & distribute with attribution.

Papercraft Spectrometer



Point at a compact
 fluorescent to calibrate

Fold flap back and tape flat

peel apart



Cut DVD, trim and tape
 onto reverse side
 here

