

## Product Specifications

## RHINOstic™ Nasal Swab Collection Device

**Description:**

Innovative and proprietary, fully molded flexible sterile nasal swab with integral de-capper cap. Increase testing rate exponentially over typical manual process flow. Patents pending.

**Physical Properties:**

Material: Polypropylene

Color: Clear

Overall length: 1.92 inches (48.7mm)

Diameter of collection zone: .12 inches (3mm)

Length of collection zone: .61inches (15.4mm)

Cap diameter: .34inches (8.6mm)

**Qualities:**

Medical grade quality

Highly flexible comfortable tip

High article collection capacity

Fits standard 1 ml cryotube

Made in the USA

**Applications:**

Anterior nares nasal swab collection

**Contact Information:**

**Rhinostics Inc.**  
179 Concord Road  
Wayland, MA 01778  
(877) 746-6789  
(877) RINOSTX



**Rhinostics Inc.**  
quality@rhinostics.com  
www.rhinostics.com



RH-100001/LV-2D00001/LV-2D1D001/LV-2D1D002 RHINOstic™ Nasal Swab Collection Device



RH-100001/LV-2D00001/LV-2D1D001/LV-2D1D002 RHINOstic™ Nasal Swab Collection Device