

Authorized CAD/CAM devices for GC materials

Manufacturer	Device	CERASMART270	Initial LRF	Initial LiSi Block
Dentsply Sirona	CEREC MC	update 4.6.1	update 4.6.1	update 4.6.1
Dentsply Sirona	CEREC MC X	update 4.6.1	update 4.6.1	update 4.6.1
Dentsply Sirona	CEREC MC XL	update 4.6.1	update 4.6.1	update 4.6.1
Dentsply Sirona	inLab MC L/XL	update 16.1	update 16.1	update 16.1
Dentsply Sirona	inLab MC X5	update 16.1	update 16.1	update 16.1
Dentsply Sirona	PrimeMill	update 5.1.3	update 5.1.3	update 5.1.3

Please note Initial LiSi Block can be fired in the SpeedFire CEREC SW device and the process lasts 12 minutes.