


**When ordering Replacement Bulbs for your Illuminator Cabinet -- Please refer to directions below:**

1. Locate your Illuminator Cabinet Label on the Right Side (as you face the cabinet) - Above your Power Switch
2. Determine your cabinet Model or Cat Number which will begin with 2425. See Diagram A Below for examples. If your cabinet label is not present in this list, please [review this message](#) before continuing to step #3.
3. After you have located your Cat 2425 Number, please proceed to Diagram B for ordering instructions.

**DIAGRAM A**  
Cabinet Label Samples



**DIAGRAM B**

Choosing the Right Bulb  
Each Cabinet requires (2) bulbs.

- Step #1: Match your model and voltage on the left below**  
**Step #2: Choose Pre-Burned or No Pre-Burn**  
**Step #3: Order (2) of the bulb ID# found to the right of your model and voltage.**

<u>Model:</u>	<u>Voltage:</u>	<u>Pre-Burned 96 Hours</u>	<u>No Pre-Burn</u>
2425	120	Bulb2427D-BI-120V1	Bulb2427D-120V1
2425E	230	Bulb2427D-BI-230V1	Bulb2427D-230V1
2425-v2	120	Bulb2427D-BI-120V2	Bulb2427D-120V2
2425E-v2	230	Bulb2427D-BI-230V2	Bulb2427D-230V2
2425Ev3	99-264	Bulb2427D-BI-V3	Bulb2427D-V3

**Still not sure? That is okay!**

Take a photo of your cabinet label and email it to us at: [info@precision-vision.com](mailto:info@precision-vision.com)  
 or Call us at 815-223-2022. We are here to serve!