

ACTUAL ACUITY
 LETTER SIZE $\frac{2}{\text{METER}}$ $\frac{6.5}{\text{FOOT}}$

LOGARITHMIC VISUAL ACUITY CHART "ETDRS"
 WITH NOTATIONS FOR TESTING AT 2 METERS (6.5 FEET)
 CHART "2"

EQUIVALENT VALUES
 DECI-MAL $\frac{6}{\text{METER}}$ $\frac{20}{\text{FOOT}}$ Log MAR

20M $\frac{2}{20}$ $\frac{6.5}{65}$

Z R K D C

0.10 $\frac{6}{60}$ $\frac{20}{200}$ 1.0

16M $\frac{2}{16}$ $\frac{6.5}{52}$

D N C H V

0.12 $\frac{6}{48}$ $\frac{20}{160}$ 0.9

12.5M $\frac{2}{12.5}$ $\frac{6.5}{41}$

C D H N R

0.16 $\frac{6}{38}$ $\frac{20}{125}$ 0.8

10M $\frac{2}{10}$ $\frac{6.5}{32.5}$

— R V Z O S —

0.20 $\frac{6}{30}$ $\frac{20}{100}$ 0.7

8M $\frac{2}{8}$ $\frac{6.5}{25.8}$

O S D V Z

0.25 $\frac{6}{24}$ $\frac{20}{80}$ 0.6

6.3 $\frac{2}{6.3}$ $\frac{6.5}{20.5}$

N O Z C D

0.32 $\frac{6}{19}$ $\frac{20}{63}$ 0.5

5M $\frac{2}{5}$ $\frac{6.5}{16}$

— R D N S K —

0.40 $\frac{6}{15}$ $\frac{20}{50}$ 0.4

4M $\frac{2}{4}$ $\frac{6.5}{13}$

O K S V Z

0.50 $\frac{6}{12}$ $\frac{20}{40}$ 0.3

3M $\frac{2}{3}$ $\frac{6.5}{10}$

K S N H O

0.63 $\frac{6}{9.5}$ $\frac{20}{32}$ 0.2

2.5M $\frac{2}{2.5}$ $\frac{6.5}{8}$

H O V S N

0.80 $\frac{6}{7.5}$ $\frac{20}{25}$ 0.1

2M $\frac{2}{2}$ $\frac{6.5}{6.5}$

— V C S Z H —

1.00 $\frac{6}{6}$ $\frac{20}{20}$ 0.0

1.5M $\frac{2}{1.5}$ $\frac{6.5}{5}$

C Z D R V

1.25 $\frac{6}{4.8}$ $\frac{20}{16}$ -0.1

1.2M $\frac{2}{1.2}$ $\frac{6.5}{4}$

— S H R Z C —

1.60 $\frac{6}{3.8}$ $\frac{20}{12.5}$ -0.2

1M $\frac{2}{1}$ $\frac{6.5}{3.2}$

D N O K R

2.00 $\frac{6}{3}$ $\frac{20}{10}$ -0.3