

When ordering Replacement Bulbs for your Illuminator Cabinet -- Please refer to directions below:

- 1. Locate your Illuminator Cabinet Label on the Right Side (as you face the cabinet) - Above your Power Switch
- 2. Determine your cabinet Model or Cat Number which will begin with 2425. See Diagram A Below for examples. If your cabinet label is not present in this list, please [review this message](#) before continuing to step #3.
- 3. After you have located your Cat 2425 Number, please proceed to Diagram B for ordering instructions.

DIAGRAM A
Cabinet Label Samples




DIAGRAM B
Choosing the Right Bulb
Each Cabinet requires (2) bulbs.

- Step #1: Match your model and voltage on the left below**
Step #2: Choose Pre-Burned or No Pre-Burn
Step #3: Order (2) of the bulb ID# found to the right of your model and voltage.

<u>Model:</u>	<u>Voltage:</u>	<u>Pre-Burned 96 Hours</u>	<u>No Pre-Burn</u>
2425	120	2427D-BI-120V1	2427D-120V1
2425E	230	2427D-BI-230V1	2427D-230V1
2425-v2	120	2427D-BI-120V2	2427D-120V2
2425E-v2	230	2427D-BI-230V2	2427D-230V2
2425Ev3	99-264	2427D-BI-V3	2427D-V3

Still not sure? That is okay!

Take a photo of your cabinet label and email it to us at: info@precision-vision.com
or Call us at 815-223-2022. We are here to serve!