



## SURESTOP EMERGENCY TRAUMA BANDAGE

### ABSORBABLE HEMOSTATIC TRAUMA BANDAGE

The Surestop Trauma Dressing is changing performance at the point of care and at the secondary level of trauma care. It has many advantages over traditional gauze agents.

When Surestop TD comes into contact with blood and exudates, it quickly transforms into a viscous gel. It fills wound voids and blocks capillary ends.

No mechanical removal is necessary, Surestop TD leaves no secondary wounds.



Controls the following bleeds:

- Femoral artery
- Carotid artery
- Class V liver bleeds

FDA 510(k)  
CE MARK Class IIB  
Utilizes Oxidized Regenerated Cellulose (ORC)



#### FEATURES:

- Neutral pH of 6.4
- Bacteriostatic
- Works under anti-coagulation treatments
- Moist wound environment
- No secondary wounds
- Water soluble
- Bio-absorbable
- Sterile
- No known side effects
- Treats major bleeding situations

PRODUCT CODE	DESCRIPTION	UNIT	UNIT/BOX	BOX/CTN
SSSBG-8	10cm x 100cm / 4"x 40"	10 PCS	1	40



## SURESTOP HEMOSTATIC DRESSING STERILE SOLUBLE HEMOSTATIC DRESSING

This dressing comes in seven sizes ranging from 1.2 x 1.2" in. up to 4 x 8" in. It has been clinically proven to stop bleeding nearly 50% faster than its competitors. It has a near neutral Ph, there are zero side effects, it's plant based and it turns to a bio absorbable gel. This product works with anti-coagulants and is bacteriostatic. When First Responders utilize SureStop to stop bleeding, upon arrival at the ER it doesn't need to be removed like traditional gauze which interrupts the clotting and exposes the wound to bacteria. It is a water soluble gel and is absorbable in the body, so there is no need to pull it away and irritate the wound.



- Innovative sterile soft white knitted fabric
- Widely used in first aid and surgery hemostats
- Quickly turns into a viscous gel on contact with blood or exudates
- Actively promotes the clotting system to help stop the bleeding
- Wide range of uses

PRODUCT CODE	DESCRIPTION	UNIT	UNIT/BOX	BOX/CTN
SSSBG-1	3cm x 3cm / 1.2"x 1.2"	12 PCS	10	6
SSSBG-2	5cm x 3cm / 2" x 1.2"	12 PCS	10	6
SSSBG-5	5cm x 5cm / 2" x 2"	12 PCS	10	6
SSSBG-3	7cm x 5cm / 2.8" x 2"	12 PCS	10	6
SSSBG-4	10cm x 8cm / 4"x 3.2"	12 PCS	10	6
SSSBG-6	10cm x 10cm / 4" x 4"	10 PCS	10	6
SSSBG-7	10cm x 20cm / 4" x 8"	5 PCS	10	6
SSSBG-D	1.9cm x 1.9cm / 0.76" x 0.76"	20 PCS	1	200





## SURESTOP ISLAND DRESSING CONFORMABLE ISLAND DRESSING

This item has a non-woven fabric backing with a large strip of absorbable material imbedded with SureStop hemostatic dressing. These are essential when there is a focus on stopping the bleeding while at the same time protecting large wounds or incisions from bacteria, friction and contact with foreign items. These are ideal for Operating Rooms and Post Op Care.



- Conforms comfortably to most body surfaces.
- Provides a sterile barrier over the wound.
- Decreases the need for multiple layers of gauze, abdominal/surgical pads and tapes.
- Allows for effective drainage of exudates.
- Non adherent pad will not stick to wound.
- Wide range of uses.

PRODUCT CODE	DESCRIPTION	UNIT	UNIT/BOX	BOX/CTN
SSCG-2510	25cm x 10cm / 10"x 4"	12 PCS	1	20
SSCG-74	7cm x 4cm / 2.8"x1.6"	20 PCS	1	32



## SURESTOP NASAL DRESSING STERILE NASAL DRESSING

SureStop nasal dressing with string is an absorbable gauze that has a full sleeve of SureStop encasing the gauze. When placed in the nasal cavity for bleeding noses and nasal injuries, the gauze absorbs and expands thereby pushing the SureStop dressing onto the bleeding capillaries and stops the bleeding. Once the strip is ready to be removed, there is no irritation of the clotting as the gel that has transformed provides a lubricated surface and remains in the nasal cavity to continue to control the bleeding. Great for football, soccer and other contact sports.



- Nasal dressing sponge expands on contact with blood or exudates.
- The expanded sponge provides pressure to help stop bleeding.
- Easy removal.
- Recommended for short term use.



PRODUCT CODE	DESCRIPTION	UNIT	UNIT/BOX	BOX/CTN
SSNDSNG-102015	10cm x 2.0cm x 1.5cm / 4"x0.8"x0.6"	10 PCS	10	4
SSNDSNG-082015	8cm x 2.0cm x 1.5cm / 3.2"x0.8"x0.6"	10 PCS	10	4
SSNDSNG-552015	5.5cm x 2.0cm x 1.5cm / 2.2"x0.8"x0.6"	10 PCS	10	4



## SURESTOP NASAL DRESSING WITH GAUZE

### STERILE HEMOSTATIC NASAL DRESSING

SureStop nasal dressing with string is an absorbable gauze that has a full sleeve of SureStop encasing the gauze. When placed in the nasal cavity for bleeding noses and nasal injuries, the gauze absorbs and expands thereby pushing the SureStop dressing onto the bleeding capillaries and stops the bleeding. Once the strip is ready to be removed, there is no irritation of the clotting as the gel that has transformed provides a lubricated surface and remains in the nasal cavity to continue to control the bleeding. Great for football, soccer and other contact sports.

WE'VE GOT BLEEDING  
COVERED



- Expansive gauze covered sponge Nasal Dressing
- Expands on contact with blood or exudates
- Provides pressure once expanded to stop bleeding
- Easy removal

PRODUCT CODE	DESCRIPTION	UNIT	UNIT/BOX	BOX/CTN
SSNDGS-102015	10cm x 2.0cm x 1.5cm / 4"x0.8"x0.6"	10 PCS	10	4
SSNDGS-082015	8cm x 2.0cm x 1.5cm / 3.2"x0.8"x0.6"	10 PCS	10	4
SSNDGS-552015	5.5cm x 2.0cm x 1.5cm / 2.2"x0.8"x0.6"	10 PCS	10	4
SSANDS-80215	8cm x 2.0cm x 1.5cm / 3.2"x0.8"x0.6"	10 PCS	10	4





## SURESTOP FIRST AID STRIP & DOT DRESSING STERILE FIRST AID STRIP DRESSING

This item was initially designed for dialysis applications. It has the appearance of a traditional band-aid but it is wrapped with SureStop dressing. Many nurses prefer this product for post catheter care/post IV care. The first aid strip is simply placed over the injection site and stops bleeding within seconds rather than minutes. This eliminates medical clinics having to burn soft dollars on staff taking more time with patients.



- Non woven strip first aid dressing
- Quickly helps to stop the bleeding
- Non adherent pad will not stick to the wound
- Conforms well to most body extremities

PRODUCT CODE	DESCRIPTION	UNIT	UNIT/BOX	BOX/CTN
SSCG-10018	10cm x 1.8cm / 4"x .72"	20 PCS	1	30

### DOT DRESSING

PRODUCT CODE	DESCRIPTION	UNIT	UNIT/BOX	BOX/CTN
SSDD-74	7cm x 4cm / 2.8"x 1.6"	20 PCS	1	32



## SURESTOP HYDROGEL ADHESIVE DRESSING

### STERILE ADHESIVE HYDROGEL WOUND DRESSING

SureStop Hydrogel adhesive is the same structure as the Hydrogel Dressing but has an additional adhesive film to affix the dressing to a wound. This provides the patient with a waterproof environment that allows a patient to bath. They can even use the Hydrogel adhesive dressing for coverage during longer periods that require freedom of movement



Suitable for use on:  
Partial thickness  
burns, skin donor sites,  
wound and epidermis  
defects, leg ulcers,  
pressure sores, and  
diabetic ulcers

- Breathable sterile island film dressing for wounds and burns.
- Absorbs water and exudates
- Hydrophillic 3D structure provides optimal moist wound healing environment
- Semi transparent for wound observation
- Protects against abrasion, friction, desiccation, and contamination
- Will not affect general or local drug administration if used on infected wounds

PRODUCT CODE	DESCRIPTION	UNIT	UNIT/BOX	BOX/CTN
SSAH-7545	7.5cm x 4.5cm / 3"x 1.8"	12 PCS	1	80



## RESPONDER GUN SHOT WOUND KIT HEMOSTATIC TAMPONADE & CHEST SEAL

Designed for the First Responder or active military member that is in the line of fire. These kits are designed to be placed in the agents pocket for easy access. Its compact, lightweight nature makes it easy to carry and practical for everyday packing. The vacuum sealed package can be quickly torn open for access to two life saving devices: Surestop Hemostatic Tamponade and an adhesive chest seal.

The tamponade is simply inserted into the wound and begins working immediately. This item is embedded with the same life saving agent as the Trauma bandage. Once inserted, the chest seal can be placed over the wound to help control blood loss and protect the wound.

This life saving combination is effective and practical for agents around the globe.

WE'VE GOT BLEEDING  
**COVERED**

### Tamponade

- Expands on contact with blood or exudates.
- Includes easy removal pull cord.

### Chest Seal

- Uses an advanced hydrogel adhesive that sticks over hair and debris.
- Transparent and suitable for burns.



PRODUCT CODE	DESCRIPTION	UNIT	UNIT/BOX	BOX/CTN
4441001	4.8" x3.2" Dressing & Tamponade	1 PCS	5	25
4441002	Two 6" Chest Seals & Two Tamponades	1 PCS	5	25



## SURESTOP HYDROGEL DRESSING STERILE HYDROGEL WOUND DRESSING



SureStop hydrogel is a hydrophilic 3D structure designed for burns and any other wound requiring hemostatis and a moist wound healing environment. Once opened and applied, it has a cooling agent to ease discomfort. It is non adherent and will not stick to irritated tissue. It is 85% water so it releases moisture into the wound while at the same time absorbing puss and other exudates from the wound site. It is semi-transparent so medical staff can monitor the wound without removing the dressing



*Suitable for use on partial thickness burns, skin donor sites, wound and epidermis defects, leg ulcers, pressure sores and diabetic ulcers.*

- Sterile wound and burn dressing.
- Free of antimicrobial agents.
- Hydrophilic 3D structure provides optimal moist wound healing environment.
- Protects against abrasion, friction, desiccation, and contamination.

PRODUCT CODE	DESCRIPTION	UNIT	UNIT/BOX	BOX/CTN
SSHD-9055	9cm x 5.5cm / 3.6"x 2.2"	12 PCS	1	40
SSHD-12080	12cm x 8cm / 4.8"x 3.2"	12 PCS	1	40

## SURESTOP HYDROGEL WOUND DRESSING STERILE HYDROGEL LIQUID



This is the liquid form of SureStop dressing. It can be applied to wounds that cannot be dressed with a traditional hemostatic dressing due to an awkward location. For instance, it works great in clinics dealing with ulcers of the feet and between toes. It is embedded with the SureStop matrix and will provide a moist wound healing environment while absorbing exudates at the same time.



- Sterile transparent wound dressing gel
- Free of antimicrobial agents
- Absorbs water and exudates
- Gently rehydrates wounds
- Protects against abrasion, friction, desiccation and contamination
- Suitable for use on partial thickness burns, skin donor sites, wound and epidermis defects, leg ulcers, pressure sores and diabetic ulcers
- Will not affect general or local drug administration if used on infected wounds

PRODUCT CODE	DESCRIPTION	UNIT	UNIT/BOX	BOX/CTN
SSAHD25	25g	1 PCS	20	24



## RESPONDER RAZOR

### PORTABLE SINGLE USE UTILITY RAZOR


This item compliments so many life saving devices responders are using in the field. Ideal for quickly preparing the patients skin in the area needing attention. The Responder Razor clears unwanted hair from sites of the hairy patients, allowing the medic to apply a chest seal, a wound dressing, or find a vein. This item uses a precision stainless steel blade and comb combination for fast, smooth shaving. The Razor is compact, quick and cost effective.

#### BENEFITS:

- Contoured round handle for comfortable grip.
- Medical single blade body hair razor.
- Individually packaged and compact.
- Easy to use.
- Safe for the hairiest of patients.

#### INCLUDES:

1. Razor with Comb
2. Stainless single blade
3. Round handle with traction coating



PRODUCT CODE	DESCRIPTION	UNIT	UNIT/BOX	BOX/CTN
4441003	Single blade with comb compact razor	1	50	30



**SURESTOP HEMOSTATIC POWDER**  
HEMOSTATIC WOUND POWDER  
*FOR ANIMAL USE ONLY*

SureStop Hemostatic powder is allowing veterinary clinics to perform procedures more efficiently and cost effectively. It has been actively tested on animals with groundbreaking results. In recent swine tests it stopped a liver bleed for 2 hours while also passing arterial and capillary tests.



- Hemostatic Powder for all types of wounds
- Quickly turns into a viscous gel on contact with blood or exudates
- Actively promotes the clotting system, to help stop the bleeding
- Wide range of uses for many types of wounds
- Puff pack for easy administration
- More than 50 doses per pack
- **FOR ANIMAL USE ONLY**

PRODUCT CODE	DESCRIPTION	UNIT	UNIT/BOX	BOX/CTN
SSHP5	5g Sterile Pack	10 PCS	1	8
SSHP150	150g Plastic Bottle	1 PCS	1	50