

# PressureCareReliefProducts®



A – Regular Bed Pad



B – Regular Bed Pad with Wings



C – Full Mattress Pad Anchor with 4 Corners Elastics




D – Full Mattress Pad Contour/Fitted



E – Regular Pillow Case  
F – Pillow Case Underlined with Polyvinyl

## PressureCareReliefProducts® – THE LATEST & EASIEST SOLUTION TO ELIMINATE THE RISK OF PRESSURE SORES.

The new, simple and cost effective way to prevent pressure sore, offering 24-hour protection from bed to chair using the same **PCRP®** pad. The thousands of polyester fibers in each layer combine to conform to the natural body contours in both supine and seated positions thereby offering pressure relief. Furthermore, air permeability tests conclude that an average of 170 cubic centimeters of air can flow through every square centimeter of surface area proving this simple product is high effective in keeping the user cool thus the findings by Dr. Steve Young in the Clinical Evaluation stating all patients noted they were much cooler and slept better on a **PressureCareReliefProducts®** over and above a product which uses pumped air.



### BENEFITS:

- Assisted in the healing Stage 2,3 & 4 Open Wounds within 4 weeks during several clinical trials.
- A recent clinical study on 6 cases (*Clinical Evaluation available upon request*) showed improved blood flow maintaining healthy tissue.
- Proven to reduce users temperature and has a cooling effect over and above other surfaces.
- No pumps, no electricity, no power needed, no adjustments required.
- Breathable, lightweight and hardwearing, *patend pending fabric density* made in Quebec, Canada. The internal Polyester fibers work in such a way to achieve minimal shearing. It can be frequently washed at 60C degrees and bleached to avoid cross infection. The pad is easily transferred from bed to seat offering 24-hour protection.
- Pad resistant to abrasion.
- Easily transportable (e.g. 2-layer pad foldable to 40cm x 25cm x 10cm, weighing 900g) so can fit in any suitcase.

### WHAT KIND OF COMPLICATIONS CAN BE CAUSED BY PRESSURE SORES?

- Can be life threatening.
- Infection can spread to the blood, heart, bone.
- Amputations.
- Prolonged bed rest.
- Autonomic dysreflexia

**DISCLAIMER** considering individual circumstances may be appropriate for the use of these products standard medical care must be followed at all times. These products are not intended to replace any procedure or guidelines in place for the prevention of pressure wounds. For additional protection whilst seated we recommend the pad is folded and always placed on a foam based chair (recliner, sofa, etc...).

Due to personal hygiene reasons and to avoid cross infection being passed from one person to another we regret we are not able to accept returns of this product once it has been opened and used.

**WARRANTY** this product is guaranteed to be free from defects in materials and workmanship under normal use for a period of 2 years.

Our PressureCareReliefProducts®			
	Code	Description & Size	
A	PCRP3038R	Regular Bed Pad 30 x 38 in (76 x 96 cm)	
B	PCRP3038W	Regulas Bed Pad with Wings 30 x 38 in (76 x 96 cm)	
C	PCRP3672A	Full Mattress Pad Anchor with 4 Corner Elastics 36 x 72 in (91 x 183 cm)	
D	PCRP3976C	Full Mattress Pad Contour/Fitted 39 x 76 in (99 x 193 cm)	
E	PCRP2127R	Regular Pillow Case	
F	PCRP2127L	Pillow Case Underlined with Polyvinyl	

Any feedbacks on our products are welcome at [casacsmedicare@gmail.com](mailto:casacsmedicare@gmail.com).

---



---



---



---



---



---

Simple & Easy to Prevent  
or Treat Pressure Sore

### CASA C.S. Inc.

9850 Ray Lawson Blvd  
Montreal (Quebec) H1J 1L1  
Tel: (514) 815-6877  
Toll Free: 1-844-900-7277  
Fax: (514) 352-3563  
[casacsmedicare@gmail.com](mailto:casacsmedicare@gmail.com)  
[www.pcrproducts.ca](http://www.pcrproducts.ca)

