

# Plumbing and Dialysis Wall Station Chase

## SPECIALTY WATER TECHNOLOGIES

### *Ultra-Pure Technologies*

**Function:** Houses wall stations and provides enclosure that conceals process plumbing and electrical connections. No degradation compared to traditional wooden chase work. Easily wipes clean. Lids provide ease of access to valves, regulators and ports. Dialysis Wall stations sold separately.

**Materials:** High Density Polyethylene Seaboard



**Mechanics:** Installed on wall to conceal process plumbing and electrical connections.

**User Quality Checks:** Follow all facility and regulatory policies and protocols regarding this enclosure. Allow only adequately trained and qualified personnel to service or repair enclosure. Routinely inspect integrity of enclosure.

**Factors Effecting Operation:** Ensure adequate drain slopes on installation.



**Specialty Water Technologies, Inc. (SWT)** provides high quality, innovative water purification systems, components and supplies primarily in the medical/dialysis market. A FDA 510k company SWT is dedicated to high quality equipment manufacturing with exemplary customer service and satisfaction. The owners and employees of SWT have a diverse background of medical, dialysis, manufacturing and business management expertise.



# Wall Station

# A-SWT-5450-XXX

## SPECIALTY WATER TECHNOLOGIES

### *Ultra-Pure Technologies*

**Function:** Provides Water, Concentrate and Bicarbonate connection and drain for hemodialysis machine.

**Materials:** Plastic Box with Color Coded Concentrate and Bicarbonate Ports, Quick Connect Water and Drain connections. Wall Stations ports are customized to facility needs. Options include pressure regulators on concentrate and bicarbonate lines, connections to Teflon, PE or PVC loops. Wall Station Labeling can be customized with company logo.



**Mechanics:** Dialysis purified water and dialysate is plumbed to the wall station with tubing and fittings.

Customized per facility need. Can be fitted into existing chase or SWT installed chase.

User connects lines from hemodialysis machine into coordinating ports.

**User Quality Checks:** : Routine inspection is recommended. Follow all facility and regulatory policies and protocols regarding this piece of equipment. Allow only adequately trained and qualified personnel to service this equipment. Inspect for leaks or broken components.

**Factors Effecting Operation:** Ensure machine connections are seated completely in ports.

Specialty Water Technologies, Inc. (SWT) provides high quality, innovative water purification systems, components and supplies primarily in the medical/dialysis market. A FDA 510k company SWT is dedicated to high quality equipment manufacturing with exemplary customer service and satisfaction. The owners and employees of SWT have a diverse background of medical, dialysis, manufacturing and business management expertise.

## SPECIALTY WATER TECHNOLOGIES

Phone: 615-654-4441

Fax: 615-654-4449

Email: sales@specialtyh2o.com