LIFT MORE THAN SCOOTER



FREERIDER INTRO



It is not the weight of owning a mobility device; it is the way of lifting it. Also, you will lift more by ours.

FREERIDER CORP

WIRING FREE

No more wiring


DRILLING FREE


No more modified your vehicle

DOCKING FREE

No more paying









INSTALL & LIFT SAFE STEP

No Tools and Helmet Required.



Attaching the hitch

Inserting linear bar

Adjusting length & height of arm

,



PLUG CONTROLLER

one goes to controller

PLUG POWER

one goes to the scooter





Connecting power before folding.





Matching the plug quickly





Tightening the remote





HOOK REAR ROD

Using Y-Strap to hook the rear rod first

PASS-THROUGH BUCKLE

Pass the buckle through the top of the front bumper.









Let's we name it Y-Strap.

The tiller can be placed in between the seat and the seat back.





LIFTING

Now, you can lift it.

ROTATING

Lifting to the proper height and rotate into the trunk.





UNPLUG CORD


Unpluging from the controller and secure to the metal ring

PULL OUT

Pull out the motor pole from the linear bar, which can stay on the hitch.








LIFTING OUR LUGGIE

Llift various Luggie and various positions.



LIFITNG YOUR STUFF

You can lift anything you feel heavy.









mower 52 lb BBQ grill
84 lb

portable freezer 75 lb

HEAVY STUFF



DELIVERING

CAMPING

FISHING

GOOD FOR EVERY ACTIVITIES



PORTABLE POWER SUPPLY KITS

Docking Station

Luggie does not need to be standing by.

SuperCharger

Juice up the battery quickly

24V li-ion battery

Lightweight and powerful li-ion battery

CONTACT INFORMATION

Office Address

8696 Utica Ave Rancho Cucamonga CA 91730

Email Address

sales@freeriderusa.com

Phone Number

909-466-7996