



OR-ganizer™

Instructions for Use

Intended Use: Laparoscopic Instruments organization for quick exchanges during surgery


Instructions:

- 1- Decide where to place the OR-ganizer, pick a spot that allows for ergonomic access to the laparoscopic instruments and in a way where the tip of the instrument would not hit the patient's body. Keep in mind that it may be placed on the skin, the drape or on both simultaneously.
- 2- Peel Liner
- 3- Adhere to skin or drape

Handling, sterilization and storage:

- Single use only
- Do not re-sterilize
- Do not use if packaging is damaged
- Please examine expiration date and do not use if expired Sterilized with radiation
- Latex Free
- Prescription only
- Store in a cool and dry place

Entabi Innovation, Inc.
1070 N. Cherry St.
Tulare, CA 93291
Tel/Fax: (833) 3Entabi
www.entabiinnovation.com



Patent Pending

Proudly made in the US
Jan,2019