



Krishna's
Educational Publishers
Since 1942



ADVANCED

Practical Physical CHEMISTRY

J.B. Yadav

BRIEF CONTENTS

Dedication.....	(iii)
Preface.....	(iv)
To The Students.....	(v)
Brief Contents.....	(vi)
Detailed Contents.....	(vii-xxiv)

Chapter 1: General Information.....	(01-16)
Chapter 2: Errors in Measurements and Treatment of Experimental Data.....	(17-26)
Chapter 3: Molecular Weights of Volatile Substances.....	(27-38)
Chapter 4: Determination of Molecular Weights in Solution.....	(39-66)
Chapter 5: Density of a Liquid.....	(67-76)
Chapter 6: Viscosity.....	(77-90)
Chapter 7: Surface Tension.....	(91-116)
Chapter 8: Vapour Pressure of Liquids.....	(117-122)
Chapter 9: Solubility.....	(123-134)
Chapter 10: Phase Equilibria.....	(135-168)
Chapter 11: Transition Temperature.....	(169-176)
Chapter 12: Thermochemistry.....	(177-200)
Chapter 13: Optical Measurements in Chemistry.....	(201-264)
Chapter 14: Electric Conductance of Electrolytes.....	(265-298)
Chapter 15: Transport Numbers.....	(299-308)
Chapter 16: Chemical Kinetics.....	(309-364)
Chapter 17: The Colloidal State.....	(365-376)
Chapter 18: Absorption.....	(377-382)
Chapter 19: Chromatographic Separations.....	(383-416)
Chapter 20: Electromotive Force of Cells.....	(417-470)
Chapter 21: Polarography.....	(471-490)
Chapter 22: Measurement of Electric and Magnetic Properties.....	(491-504)
Chapter 23: Equilibrium and Dissociation Constants.....	(505-516)
Chapter 24: Blood Examination.....	(517-528)
Chapter 25: Urine Examination.....	(529-534)
Chapter 26: Estimation of Casein in Milk.....	(535-536)
Chapter 27: Oils and Lubricants.....	(537-546)
Chapter 28: Water Analysis.....	(547-568)
Appendix.....	(01-22)
Log Tables.....	(23-26)

DETAILED CONTENTS

Chapter 1: General Information(01–16)

1.1	Calibration of Weights	01
1.1.1	Method using Standardised Weights	01
1.1.2	Kohlrausch's Method	01
1.2	Cleaning of Volumetric Glassware	03
1.3	Calibration of Volumetric Glassware	03
1.3.1	Calibration of Volumetric Flask	03
1.3.2	Calibration of Pipettes	04
1.3.3	Calibration of Burettes	04
1.4	Concentration of Solution	05
1.4.1	Expressing Concentration in Physical Units	05
1.4.2	Expressing Concentration in Chemical Units	05
1.5	Preparation of Standard Solutions	06
1.5.1	Calculation of Composition of a Mixture in Terms of Mole Fraction	11
1.6	Preparation of Solutions of Indicators	11
1.6.1	Drying Agents	13
1.6.2	Purified Water	13
1.7	Temperature Control Thermostat	14

Chapter 2: Errors in Measurements and Treatment of Experimental Data.....(17–26)

2.1	Errors in Experimental Data	17
2.2	Accuracy and Precision of Measurements	18
2.3	Statistical Treatment of Random Errors	18
2.3.1	Accuracy of Results	21
2.3.2	Calculations	23
2.3.3	Significant Figures	24

Chapter 3: Molecular Weights of Volatile Substances.....(27–38)

3.1	Vapour Density	27
Exp. 3.1:	Determine the molecular weight of a given volatile substance	29
Exp. 3.2:	Determine the molecular weight of a volatile compound, say chloroform, by Duma's method.	35
3.1.1	Observations and Calculations	36
Exp. 3.3:	Determine the composition of a given binary solution made up of volatile liquids by vapour density measurements.	37
Exp. 3.4:	Determine the solubility of carbon disulphide in methyl alcohol at room temperature by vapour density method.	37

Chapter 4: Determination of Molecular Weights in Solution (Cryoscopic and Ebullioscopic Methods).....(39–66)

4.1	Introduction	39
4.2	Measurement of Depression of Freezing Point (Cryoscopy)	39

7.3.3	Drop Fall Method	
Exp. 7.4:	Determine the surface tension of nitrobenzene at 30°C by the drop weight method.	100
Exp. 7.5:	Determine the surface tension of a liquid, say carbon tetrachloride at 30°C by stalagometer method.	101
7.3.4	du Noy's Tensiometric or Torsion Balance Method	
Exp. 7.6:	Study the variation of surface tension of a mixture of ethanol and water with composition by torsion balance method.	103 106
Exp. 7.7:	Determine molecular surface energy and the association factor for a liquid say ethyl alcohol.	107
Exp. 7.8:	Determine the percentage composition of a given mixture of two liquids, say benzene and nitrobenzene, by surface tension measurements.	108
Exp. 7.9:	Study the variation of surface tension of solutions of <i>n</i> -propyl alcohol with concentration and hence determine the limiting cross-sectional area of alcohol molecule.	109
Exp. 7.10:	Determine the critical micelle concentration of a soap (sodium laurate, sodium palmitate etc.) by surface tension measurements.	109
Exp. 7.11:	Determine the surface tension of a mixture of two miscible liquids and hence the parachor of the mixture.	110
Exp. 7.12:	Determine the parachor of a solid organic substance soluble in an organic solvent, say <i>p</i> -dichlorobenzene, <i>p</i> -dibromobenzene, benzoic acid, naphthalene etc. in benzene.	111
Exp. 7.13:	Determine the interfacial tension between two immiscible liquids, say benzene and water, at room temperature.	111 112
Exp. 7.14:	Determine the surface tension of the following liquids : methyl acetate, ethyl acetate, <i>n</i> -heptane, chloroform and ethylene, bromide and hence calculate the atomic parachors of C, H, Cl and Br.	113
Exp. 7.15:	Compare cleansing powers of two samples of a detergent.	114
Exp. 7.16:	Study the variation of surface tension of water with concentration of a soap or detergent and determine the critical micelle concentration of the detergent.	115

Chapter 8: Vapour Pressure of Liquids.....(117-122)

8.1	Introduction	117
8.2	Measurement of Vapour Pressure	118
Exp. 8.1:	Determine the vapour pressure of a pure liquid (benzene) at a series of temperatures, and hence determine the heat of vaporisation of the liquid.	118
Exp. 8.2:	Study the variation of vapour pressure of a pure liquid (carbon tetrachloride) using an isoteniscope.	119
Exp. 8.3:	Determine the molar latent heat of vaporisation of ethyl acetate by measuring its vapour pressure at a series of temperatures by the method of the Ramsay and Young.	121

Chapter 9: Solubility(123-134)

9.1	Solubility	123
9.1.1	Solubility of Solids in Liquids	123
9.1.2	Determination of Solubility	124
Exp. 9.1:	Determine the solubility of a given inorganic salt, e.g., KCl, NaCl, KNO ₃ , NaNO ₃ , K ₂ SO ₄ , etc. at different temperatures and obtain the solubility curve	124
Exp. 9.2:	Determine the solubility of sodium bromide dihydrate (NaBr.2H ₂ O) at a series of temperatures and hence the transition temperature of the salt.	126

Exp. 9.3:	Determine the heat of solution of a given substance (e.g., oxalic acid, benzoic acid, etc.) by solubility method.	126
Exp. 9.4:	Study the effect of addition of an electrolyte such as NaCl, KCl, Na ₂ SO ₄ , K ₂ SO ₄ etc. on the solubility of an organic acid (benzoic acid or salicylic acid) at room temperature.	128
Exp. 9.5:	Study the variation of solubility of Ca(OH) ₂ in NaOH solution and hence determine its solubility product.	128
Exp. 9.6:	Study the variation of the solubility of AgBrO ₃ in KBrO ₃ solutions and hence determine the solubility product of AgBrO ₃ .	130
Exp. 9.7:	Study the influence of the ionic strength on the solubility of CaSO ₄ and hence determine its thermodynamic solubility product and mean ionic activity.	131
Exp. 9.8:	Study the variation of solubility of potassium hydrogen tartrate with ionic strength using a salt having common ion and hence determine the mean ionic activity coefficients.	133

Chapter 10: Phase Equilibria.....(135–168)

10.1	Liquid-Vapour Equilibria	135
10.1.1	Completely Miscible Liquids	135
Exp. 10.1:	Determine the boiling point-composition curve for the following binary liquid mixtures of two miscible liquids :	136
	(a) A – Benzene (80·2°) B – Toluene (110·6°)	
	(b) A – Benzene (80·2°) B – Hexane (69·0°)	
	(c) A – Benzene (80·2°) B – Carbon disulphide (46·3°)	
	(d) A – Acetone (56·3°) B – Water (100°)	
	(e) A – Benzene (80·2°) B – Ether (34·6°)	
	(f) A – Chlorobenzene (152°) B – Bromobenzene (156·2°)	
Exp. 10.2:	Determine the boiling point-composition curve for the following binary mixtures of miscible liquids.	138
	(a) A – Ethanol (78·3°) B – Benzene (80·2°)	
	(b) A – Methanol (64·7°) B – Benzene (80·2°)	
	(c) A – Iso propyl alcohol B – Benzene (80·2°)	
	(d) A – Carbon tetrachloride (76·8°) B – Methanol (64·7°)	
	(e) A – Chloroform (61·2°) B – Methanol (64·7°)	
	(f) A – Acetone (56·3°) B – Carbon tetrachloride (76·8°)	
Exp. 10.3:	Construct the boiling point-composition curves for completely miscible binary liquid systems forming maximum boiling point azeotropes.	138
	(a) A — Acetone (56·3°) B — Chloroform (61·2°)	
	(b) A — Water (100°) B — Formic acid (99·2°)	
10.1.2	Completely Immiscible Liquids (Steam Distillation)	138
Exp. 10.4:	Determine the vapour pressure of chlorobenzene by steam distillation.	139
Exp. 10.5:	Determine the molecular weight of benzene, nitrobenzene, toluene, etc. by steam distillation.	141
10.2	Liquid-Liquid System	141
10.2.1	Partial Miscibility of Liquids (Mutual Solubility)	141
Exp. 10.6:	Determine the mutual solubility curve of phenol and water and hence the consolute point.	142
Exp. 10.7:	Determine the critical solution temperature and the composition of the solution at CST for the following pairs of partially miscible liquids :	144
	(a) Aniline-hexane, (b) Aniline-water,	
	(c) Methanol-cyclohexane and (d) Carbon disulphide-methanol	

Exp. 11.1:	Determine the transition temperature of a given substance (hydrated sodium bromide) by thermometric method.	170
11.2.2	Dilatometric Method	172
Exp. 11.2:	Determine the transition temperature of sodium sulphate decahydrate ($\text{Na}_2\text{SO}_4 \cdot 10\text{H}_2\text{O}$) by dilatometric method.	172
Exp. 11.3:	Determine the transition temperature of hydrated sodium bromide by solubility method.	174
Exp. 11.4:	Determine the transition temperatures of calcium chloride by solubility method in the temperature range 15-60°C.	175

Chapter 12: Thermochemistry.....(177-200)

12.1	Introduction	177
12.2	Measurement of Heat Change	180
Exp. 12.1:	Determine the heat of neutralization of hydrochloric acid using a classical calorimeter.	181
Exp. 12.2:	Determine the heat of neutralization of sulphuric acid using Dewar's vacuum flask as the calorimeter.	185
Exp. 12.3:	Determine the enthalpy of neutralization of a weak acid, say acetic acid and hence calculate its enthalpy of ionisation.	187
Exp. 12.4:	Determine the enthalpy of ionisation of a weak base, NH_4OH , calorimetrically.	188
Exp. 12.5:	Determine enthalpy of neutralization of a weak acid, say CH_3COOH , with a weak base, say NH_4OH .	188
Exp. 12.6:	Determine the enthalpies of neutralization of two acids, e.g., HCl and CH_3COOH and hence their relative strength.	189
Exp. 12.7:	Determine the enthalpy changes for the different stages of neutralization of a polybasic acid (mol. mass 126) and hence find the basicity of the acid.	189
Exp. 12.8:	Determine the integral heat of solution of a salt such as KNO_3 , KCl etc. in water using Dewar's vacuum bottle as calorimeter.	190
Exp. 12.9:	Determine the enthalpy of hydration of CuSO_4 , i.e., the heat of crystallisation of $\text{CuSO}_4 \cdot 5\text{H}_2\text{O}$.	192
Exp. 12.10:	Determine the integral heat of solution of a salt at two concentrations and hence the integral heat of dilution.	192
Exp. 12.11:	Determine the integral heats of dilution of sulphuric acid starting with solutions of different concentrations.	192
Exp. 12.12:	Determine the enthalpies of hydrolysis of (i) HC_2O_4^- ion and (ii) CO_3^{2-} ion.	193
Exp. 12.13:	Determine enthalpy of precipitation of a salt, say BaSO_4 .	194
Exp. 12.14:	Determine the enthalpy of transition of $\text{Na}_2\text{SO}_4 \cdot 10\text{H}_2\text{O}$ by calorimetry.	195
Exp. 12.15:	Determine the enthalpy of displacement of copper or silver from aqueous solution of their salts by zinc, iron or magnesium.	195
Exp. 12.16:	Determine the enthalpy of formation of (a) zinc oxide and (b) magnesium oxide calorimetrically.	196
Exp. 12.17:	Determine the calorific value (fuel value) of a given fuel using a bomb calorimeter.	197

Chapter 13: Optical Measurements in Chemistry.....(201-264)

13.1	Introduction	201
13.2	Refractometry	201
Exp. 13.1:	Determine the refractive index of a given liquid by Abbe refractometer and find the specific and molar refractions.	207

Exp. 14.4:	Determine the relative strength of chloroacetic acid and acetic acid by conductance measurements.	280
Exp. 14.5:	Determine the solubility of lead sulphate in water at 25°C.	281
Exp. 14.6:	Determine the hydrolysis constant of aniline hydrochloride at 25°C.	283
Exp. 14.7:	Determine the degree of hydrolysis and hydrolysis constant of (i) CH_3COONa and (ii) NH_4Cl .	284
Exp. 14.8:	Determine the basicity of an organic acid, say tartaric acid, by conductometric method.	284
Exp. 14.9:	Determine the end point of the following titrations by the conductometric method.	287
Exp. 14.10:	Titrate a moderately strong acid like salicylic acid or malonic acid conductometrically.	287
Exp. 14.11:	Estimate conductometrically the concentration of a dilute solution (0.01 to 0.001 N) of a weak acid.	289
14.4	Titration of Very Weak Acids with a Strong Base	290
14.5	Titration of Very Weak bases with Strong Acid	291
Exp. 14.13:	Titrate a very weak base, like <i>p</i> -toluidine with HCl conductometrically.	291
Exp. 14.14:	Determine the composition of a mixture of acetic acid and hydrochloric acid by conductometric titration.	291
Exp. 14.15:	Titrate the following acids solution, with NaOH, conductometrically.	292
Exp. 14.16:	Titrate a mixture of hydrochloric acid and oxalic acid conductometrically.	292
Exp. 14.17:	Titrate a solution containing acetic acid and oxalic acid.	293
Exp. 14.18:	Estimate normalities of nitric acid and sulphuric acid in a mixture of the two acids.	293
14.6	Precipitation Titrations	293
Exp. 14.19:	Carry out the following precipitation titrations conductometrically :	294
	(a) 50 cm^3 0.02 N AgNO_3 with <i>N</i> -HCl	
	(b) 50 cm^3 0.02 N AgNO_3 with <i>N</i> -KCl	
	(c) 50 cm^3 0.004 N MgSO_4 with 0.1 N-Ba(OH) ₂	
	(d) 50 cm^3 0.002 N BaCl_2 with <i>N</i> - Li_2SO_4	
	(e) 50 cm^3 0.02 N BaCl_2 with <i>N</i> - K_2SO_4	
Exp. 14.20:	Determine conductometrically the concentrations of KCl and KI in a given solution.	295
Exp. 14.21:	Titrate 50 cm^3 of 0.01 M KI solution with 0.5 M HgCl_2 solution.	295
14.7	Miscellaneous Titrations	295
Exp. 14.22:	Determine the concentration of a salt of a weak acid and a strong base like sodium acetate (or sodium oxalate) solution conductometrically.	295
Exp. 14.23:	Conductometrically estimate the concentration (0.05 M) of NH_4Cl (a salt of weak base and a strong acid).	296
Exp. 14.24:	Determine the concentration of each component in the following solutions :	296
	(a) Hydrochloric acid and ammonium chloride in a solution.	
	(b) Ammonium hydroxide and ammonium chloride in a solution.	
Exp. 14.25:	Estimate the concentration of H_2SO_4 , CH_3COOH and $\text{CuSO}_4 \cdot 5\text{H}_2\text{O}$ in a given solution by carrying out conductometric titration with NaOH solution.	297
Exp. 14.26:	Titrate a mixture of HCl, CH_3COOH and CuSO_4 with an alkali conductometrically.	297
Exp. 14.27:	Titrate conductometrically the mixture of silver nitrate and nitric acid with sodium hydroxide solution.	297
Exp. 14.28:	Estimate the concentrations of HCl, HClO_4 and H_2SO_4 in a given aqueous solution by conductometric titration.	298

Chapters 15: Transport Numbers.....	(299–308)
15.1 Introduction	299
15.2 Determination of Transport Number	300
Exp. 15.1: Determine the transference numbers of silver and nitrate ions in a solution of silver nitrate by Hittorf's method.	301
Exp. 15.2: Determine the transport numbers of copper and sulphate ions in a copper sulphate solution by Hittorf method.	305
Exp. 15.3: Determine the transport number of H^+ and Cl^- ions in HCl by moving boundary method.	306
Chapter 16 : Chemical Kinetics.....	(309–364)
16.1 Introduction	309
16.1.1 Reaction Rates	309
16.1.2 Reaction Molecularities	310
16.1.3 Order of a Reaction	310
16.1.4 First order Reaction	310
16.1.5 Second order Reaction	311
16.1.6 Third Order Reaction	313
16.2 Determination of the Order of a Reaction	314
16.3 Drawing Tangents to a Curve	317
16.4 The Temperature Coefficient (Energy of Activation)	318
Exp. 16.1: Determine the rate constant of hydrolysis of an ester such as methyl acetate catalysed by an acid, say 0.5M HCl. Determine also the energy of activation of the reaction.	319
Exp. 16.2: Determine the relative strength of two acids, say hydrochloric acid and sulphuric acid, studying the hydrolysis of an ester (say methyl acetate).	321
Exp. 16.3: Study the hydrolysis of methyl acetate catalysed by HCl solution and equinormal solution of urea hydrochloride, and hence determine the degree of hydrolysis of the salt.	322
Exp. 16.4: Investigate the inversion of cane sugar in presence of an acid, say molar HCl at 30°C (or room temperature). Determine also the energy of activation of the reaction.	323
Exp. 16.5: Study inversion of cane sugar in presence of hydrochloric acid and sulphuric acid, and hence determine the relative strength of the acids.	325
Exp. 16.6: Determine the activation energy of enzyme (invertase) catalysed inversion of sucrose.	325
Exp. 16.7: Investigate the mutarotation of glucose catalysed by an acid or base.	327
Exp. 16.8: Determine the velocity constant of hydrolysis of ethyl acetate by sodium hydroxide (saponification of an ester).	328
Exp. 16.9: Determine the order of hydrolysis of ethyl acetate by sodium hydroxide. Also study the salt effect on the rate constant of the reaction.	331
Exp. 16.10: Study the kinetics of saponification of ethyl acetate by sodium hydroxide at two temperatures by conductance measurements, and hence determine the energy of activation of the reaction.	332
Exp. 16.11: Study the kinetics of hydrolysis of tertiary amyl iodide and determine the order and energy of activation of the reaction.	334
Exp. 16.12: Investigate the reaction between hydrogen peroxide and hydrogen iodide (I^- ion in acidic medium) following the kinetics as an iodide clock reaction.	336
Exp. 16.13: Study the kinetics of decomposition of benzene diazonium chloride at a series of temperatures (between 30° to 50°C).	338
Exp. 16.14: Study the kinetics of decomposition of H_2O_2 catalysed by iodide ion.	339
Exp. 16.15: Study the kinetics of depolymerisation of diacetone alcohol catalysed by sodium hydroxide using a dilatometer.	341
Exp. 16.16: Investigate the kinetics of hydrolysis of acetal using a dilatometer.	344

Exp. 16.17:	Study the oxidation of iodide ion by persulphate ion following the kinetics as iodine clock reaction and find the order with respect to each reacting species by isolation method measuring the initial rate.	345
Exp. 16.18:	Study the kinetics of the reaction between potassium persulphate and potassium iodide. (a) Determine the rate constant and order of the reaction. (b) Study the influence of ionic strength on the rate constant.	347
Exp. 16.19:	Investigate the reaction between acetone and iodine.	349
Exp. 16.20:	Investigate the autocatalytic reaction between potassium permanganate and oxalic acid.	352
Exp. 16.21:	Study the anionotropic rearrangement of phenyl allyl alcohol to cinnamyl alcohol, and determine the rate constant of the reaction.	353
Exp. 16.22:	Study the kinetics of reaction between bromic and hydrobromic acids, and determine the order of the reaction.	355
Exp. 16.23:	Study the kinetics of decomposition of the complex formed between sodium sulphide and sodium nitroprusside spectrophotometrically, and find the order and rate constant of the reaction.	356
Exp. 16.24:	Determine colorimetrically the order and the energy of activation for the decomposition of violet coloured complex of ceric ions and N-phenylanthranilic acid.	357
Exp. 16.25:	Study the kinetics of oxidation of ethanol by potassium dichromate in acidic medium and determine (i) rate constant and (ii) energy of activation of the reaction.	358
Exp. 16.26:	Study the kinetics of the reaction between ammonia and monochloroacetic ion.	360
16.5 Photochemical Reactions		362
Exp. 16.27:	Study the decompositions of oxalic acid in solution photosensitized by uranyl sulphate.	363

Chapter 17: The Colloidal State.....(365-376)

17.1	Introduction	365
17.2	Dialysis of Colloidal Solution	366
17.3	Precipitation of Colloids by Electrolytes	366
17.4	Protection of Colloids	367
Exp. 17.1:	Prepare colloidal solutions of (i) arsenious sulphide, (ii) antimony sulphide (iii) cadmium sulphide, (iv) sulphur, (v) ferric hydroxide, (vi) silver, (vii) molybdenum blue and (viii) gelatin.	367
Exp. 17.2:	Study the coagulation of arsenious sulphide sol. with NaCl, BaCl ₂ , and AlCl ₃ and hence find their coagulating values.	369
Exp. 17.3:	Determine the precipitating value of KCl, K ₂ SO ₄ , 0.0005 M K ₃ Fe(CN) ₆ for ferric hydroxide sol.	370
Exp. 17.4:	Study the protective action of a hydrophilic colloid (such as gelatin, albumin, starch, etc.) on the precipitation of lyophobic sols.	370
Exp. 17.5:	Study the mutual coagulation of As ₂ S ₃ sol. and Fe(OH) ₃ sol. and determine the optimum ratio for precipitation.	371
Exp. 17.6:	Determine the nature of the charge carried by particles in a given sol., and measure their electrophoretic mobility and zeta potential of the sol.	371
17.5	Colloidal Electrolytes or Association Colloids	374
Exp. 17.7:	Determine the critical concentration of micelle formation of a soap by the study of change of spectral behaviour of pinacyanol chloride.	374
Exp. 17.8:	Study the salt effect on the critical concentration of micelle formation in aqueous solution of colloidal electrolytes by change of spectral behaviour of a dye.	375

Chapter 18: Adsorption.....	(377-382)
18.1 Introduction	377
18.2 Adsorption Isotherms	377
Exp. 18.1: Investigate the adsorption of oxalic acid from aqueous solutions by activated charcoal, and examine the validity of classical and Langmuir's adsorption isotherms.	379
Exp. 18.2: Determine the adsorption isotherm of acetic acid from aqueous solution by charcoal.	380
Exp. 18.3: Study the adsorption of iodine from alcoholic solution by charcoal.	381
Chapter 19 : Chromatographic Separations.....	(383-416)
19.1 Chromatography	383
19.2 Column Chromatography	384
19.3 Adsorption Chromatography	386
Exp. 19.1: Separate a mixture of methylene blue and fluorescein (sodium salt) on an alumina column.	387
Exp. 19.2: Separate a mixture of <i>o</i> - and <i>p</i> -nitroanilines on an alumina column.	388
Exp. 19.3: Separate a mixture of 2:4-dinitrophenyl hydrazones by adsorption chromatographic technique.	388
Exp. 19.4: Separate leaf pigments (chlorophylls and carotenoids) from spinach by column chromatography.	389
19.4 Paper Chromatography	391
19.4.1 One Dimensional Paper Chromatography	391
19.4.2 Two Dimensional Paper Chromatography	394
Exp. 19.5: Study the separation of organic acids by one dimensional paper chromatography.	397
Exp. 19.6: Study the separation of amino acids by one and two dimensional paper chromatography.	398
Exp. 19.7: Study the separation of sugars by one dimensional descending paper chromatography.	399
Exp. 19.8: Study the presence of lactose in milk by descending paper chromatography.	400
Exp. 19.9: Demonstrate the separation of inorganic ions by paper chromatography.	401
19.5 Thin-Layer Chromatography	403
Exp. 19.10: Demonstrate the separation of dyes (azobenzenes) by thin layer chromatography.	405
Exp. 19.11: Determine the R_f values of the given amino acids by thin layer chromatography and identify the amino acids present in a given mixture.	406
Exp. 19.12: Separate the plant pigments from the extract of spinach by thin layer chromatography.	407
Exp. 19.13: Study the separation of 2:4 dinitrophenyl hydrazones by thin layer chromatography.	408
19.6 Ion Exchange or Displacement Chromatography	408
19.7 Ion Exchange Capacity	409
Exp. 19.14: Prepare free acid or base from the salt of an organic acid (say sodium citrate) or base (say aniline hydrochloride) using cation and anion exchange resins.	410
Exp. 19.15: Determine the concentration of a salt solution by ion exchange chromatography.	411
Exp. 19.16: Determine the composition of a solution containing an acid and its salt (sodium acetate and acetic acid).	412
19.8 Paper Electrophoresis or Zone Electrophoresis	412
Exp. 19.17: Study the separation of inorganic cations by paper electrophoresis.	414

- Exp. 19.18: Study the separation of amino acids in a mixture by electrophoresis. 415
 Exp. 19.19: Determine the isoelectric point of glutamic acid by paper electrophoresis. 415

Chapter 20: Electromotive Force of Cells(417-470)

20.1	Galvanic Cells	417
20.2	Electrode Potential	417
20.3	Hydrogen Scale	418
20.4	Liquid Junction Potential	418
20.5	Types of Electrodes	419
20.6	Indicating Electrodes	420
20.7	Measurement of e.m.f. of a Cell	428
20.8	Applications of e.m.f. Measurements	431
Exp. 20.1:	Determine the valency of mercurous ions.	431
Exp. 20.2:	(a) Determine the mean ionic activity coefficient of hydrochloric acid at different concentrations using a concentration cell without transference.	432
Exp. 20.3:	Determine the activity coefficient of silver ions using a concentration cell with a salt bridge.	434
Exp. 20.4:	Study the effect of ionic strength on the activity coefficient of silver ions in 0.01 M silver nitrate solution.	436
Exp. 20.5:	Determine the heat of reaction and equilibrium constant of the reaction between metallic zinc and lead nitrate solution.	436
Exp. 20.6:	Determine the heat of reaction, entropy change and equilibrium constant of the reaction between metallic zinc and Cu^{2+} ions in solution.	438
Exp. 20.7:	Determine the transport number of Ag^+ and NO_3^- in solution using 0.1 and 0.01 M AgNO_3 solutions (Mean ionic activity coefficients of AgNO_3 in 0.01 M and 0.1 M solutions are 0.89 and 0.73).	438
Exp. 20.8:	Determine the solubility of silver chloride in water electrometrically at 25°C.	440
Exp. 20.9:	Determine the degree of hydrolysis of aniline hydrochloride in 0.001 M solutions at room temperature, and hence calculate the hydrolysis constant of the salt and dissociation constant of the base.	440
Exp. 20.10:	Determine the equilibrium constant of the reaction :	442
	$\text{Ag}(\text{NH}_3)_2^+ \rightleftharpoons \text{Ag}^+ + 2\text{NH}_3$	
Exp. 20.11:	Determine the ionic product of water by e.m.f. method at 25°C.	443
20.9	Determination of pH Values	445
Exp. 20.12:	Determine the pH of a number of buffer solutions using hydrogen electrode.	446
Exp. 20.13:	Determine the pH of various mixtures of sodium acetate and acetic acid in aqueous solutions, and hence determine the dissociation constant of the acid.	446
Exp. 20.14:	Determine the pH of a buffer solution (with pH less than 8) using a quinhydrone electrode.	447
Exp. 20.15:	Construct the calibration curve (pH versus E_{obs}) for quinhydrone electrode, and hence determine the standard oxidation potential of the quinhydrone electrode.	448
20.10	pH Standards for Glass Electrodes at 25°C	449
Exp. 20.16:	Determine the pH of a given solution using glass electrode.	449
20.11	Potentiometric Titrations	450
Exp. 20.17:	Titrate 0.1 M HCl solution with 0.1 M NaOH using hydrogen electrode.	452
Exp. 20.18:	Titrate a given solution of HCl with NaOH using quinhydrone as the indicator electrode.	454
Exp. 20.19:	Determine the strength of unknown solution of HCl titrating it with NaOH solution using pH-meter.	455

20.12	Dissociation Constants of Weak Acids and Bases	455
Exp. 20.20:	Titrate a weak acid (e.g., acetic acid, benzoic acid, lactic acid; dibasic acid with groups separated—malonic acid or oxalic acid) potentiometrically. From the results, obtain the dissociation constant of the acid.	457
Exp. 20.21:	Determine the acid and basic dissociation constants of an amino acid, and hence the isoelectric point of the acid.	458
Exp. 20.22:	Titrate a mixture of HCl and CH ₃ COOH, potentiometrically.	460
Exp. 20.23:	Titrate 0.1 M solution of ortho phosphoric acid solution with 0.5 M NaOH solution potentiometrically.	461
Exp. 20.24:	Titrate ferrous ammonium sulphate against potassium dichromate (or KMnO ₄) potentiometrically, and hence determine the formal redox potential of ferrous-ferric system.	462
Exp. 20.25:	Titrate ferrous ammonia sulphate against ceric sulphate potentiometrically. Hence, determine the formal redox potential of ferrous-ferric and cerous-ceric system.	463
Exp. 20.26:	Titrate pure solutions of KCl, KBr, KI and their mixture with standard AgNO ₃ solution potentiometrically.	463
Exp. 20.27:	Determine the solubility product of a silver halide (say AgCl) by potentiometric titration.	464
Exp. 20.28:	Titrate cobaltous solutions with potassium ferricyanide, and hence determine the concentration of cobalt (II).	465
Exp. 20.29:	Titrate potentiometrically Zn (II) by K ₄ Fe(CN) ₆ , and verify composition of the complex (precipitate) : K ₂ Zn ₃ [Fe(CN) ₆] ₂	466
20.13	Titrations without reference Electrodes	467
Exp. 20.30:	Titrate a silver nitrate solution with potassium chloride by the different potentiometric technique.	467
Exp. 20.31:	Titrate iodine solution with sodium thiosulphate by the dead stop end point or polarisation method.	468
Exp. 20.32:	Titrate ferrous ammonium sulphate solution with potassium dichromate solution potentiometrically using a bimetallic electrode pair.	470

Chapter 21: Polarography.....(471–490)

21.1	Introduction	471
21.2	Manual Polarograph	473
21.3	Commercial Polarograph	473
21.4	Decomposition Potential and Over Voltage	474
21.5	Polarographic Maxima	477
Exp. 21.1:	Plot current-voltage curves for 0.05M sulphuric acid, 0.1 M sodium sulphate, copper sulphate, zinc sulphate using platinized platinum electrodes.	477
Exp. 21.2:	Study the polarographic waves produced by dissolved oxygen.	477
Exp. 21.3:	Determine the half wave potentials of Cd ²⁺ and Zn ²⁺ ions in 0.1 M KCl solutions.	478
Exp. 21.4:	Study the variation of diffusion current with concentration, and hence construct a wave height-concentration curve for cadmium ion.	480
Exp. 21.5:	Plot a polarogram for a mixture of Cd ²⁺ , Zn ²⁺ and Mn ²⁺ ion.	480
21.6	Amperometric Titrations	481
Exp. 21.6:	Determine the concentration of lead in solution by carrying out amperometric titration with K ₂ Cr ₂ O ₇ solution.	483
Exp. 21.7:	Titrate potassium sulphate solution amperometrically with lead nitrate solution.	484
Exp. 21.8:	Determine nickel in solution by amperometric titration with dimethyl glyoxime.	484
Exp. 21.9:	Titrate amperometrically bismuth, lead and calcium in solution with EDTA.	485
21.7	Formula and Stability Constant of the Complex Ion	486

- Exp. 21.10: Determine the formula and the stability constant of a metal ion complex (lead oxalate complex). 487
- Exp. 21.11: Determine the formation (stability) constant of copper glycinate polarographically. 488

Chapter 22: Measurement of Electric and Magnetic Properties.....(491-504)

- 22.1 Capacitance and Dielectric Constant 491
- 22.1.1 Dipole Moment of Polar Molecules 491
- 22.1.2 Electric or Distortion Polarisation 492
- Exp. 22.1: Determine the dipole moment of a liquid such as chlorobenzene, chloroform, nitrobenzene etc. 494
- 22.2 Magnetic Susceptibility 496
- 22.2.1 Diamagnetism 497
- 22.2.2 Paramagnetism 497
- 22.3 Magnetic Susceptibility Measurements 500
- Exp. 22.2: Determine the magnetic susceptibility of Mohr's Salt - $(\text{NH}_4)_2\text{SO}_4 \cdot \text{FeSO}_4 \cdot 6\text{H}_2\text{O}$ at room temperature, and hence the magnetic moment. 502

Chapter 23: Equilibrium and Dissociation Constants.....(505-516)

- 23.1 Reversible Reactions 505
- 23.2 Equilibrium Constant and Standard Free Energy Change 506
- 23.3 Equilibrium Constant and Enthalpy (ΔH) of the Reaction 506
- 23.4 Determination of Equilibrium Constants 507
- Exp. 23.1: Determine the equilibrium constant of an esterification reaction between acetic acid and ethanol. 507
- Exp. 23.2: Determine the equilibrium constant of the following reversible reaction :

$$2\text{Ag}^+ + \text{CaSO}_4(s) \rightleftharpoons \text{Ag}_2\text{SO}_4(s) + \text{Ca}^{2+}$$
 509
- Exp. 23.3: Determine the equilibrium constant of the keto-enol tautomerism of ethyl acetoacetate. 511
- Exp. 23.4: Determine the dissociation constant of picric acid by studying its distribution between benzene and water. 513
- Exp. 10.17: Determination of equilibrium constant of the reaction :

$$\text{KI} + \text{I}_2 \rightleftharpoons \text{KI}_3$$
 [Partition method] 515
- Exp. 13.22: Determination of stability constant or formation constant of the complex between Fe^{3+} ion and salicylic acid. 516
- Exp. 13.24: Dissociation constant of methyl red by spectrophotometric method. 516
- Exp. 13.25: Dissociation constant of bromophenol blue. 516
- Exp. 14.1: Dissociation constant of a weak acid. 516
- Exp. 14.6: Hydrolysis constant of aniline hydrochloride. 516
- Exp. 14.7: Hydrolysis constant of sodium acetate and ammonium chloride. 516
- Exp. 20.10: Equilibrium constant of the reaction

$$[\text{Ag}(\text{NH}_3)_2]^+ \rightleftharpoons \text{Ag}^+ + 2\text{NH}_3$$
 516
- Exp. 20.20: Dissociation constant of a weak acid. 516
- Exp. 20.21: Dissociation constant (K_a and K_b) of an amino acid. 516
- Exp. 20.22: Dissociation constant of phosphoric acid $K_{a(1)}$, $K_{a(2)}$ and $K_{a(3)}$. 516

Chapter 24: Blood Examination.....(517-528)

- 24.1 Cell Counts 517
- Exp. 24.1: Determine the White Blood Cell count in a given sample of blood. 522
- Exp. 24.2: Determine the Red Blood Cell (RBC) count in a given sample of blood.

Exp. 24.3 :	Determine the platelet count in a given sample of blood.	524
Exp. 24.4 :	Assess absolute eosinophil count (AEC) in a blood sample.	525
24.2	Estimation of Glucose in Blood (Diabetes Blood Test)	526
Exp. 24.5 :	Estimate glucose concentration in blood (plasma/serum and CSF) calorimetrically.	527
Chapter 25: Urine Examination.....		(529– 534)
25.1	Introduction	529
25.2	Collection of Urine Specimen	529
Exp. 25.1 :	Detect and estimate the ketone bodies in a urine sample.	530
25.3	Estimation of Ketone Bodies (Acetoacetic Acid) (Dipstick Method)	531
25.4	Reducing Substances in Urine	532
Exp. 25.2.	(a) Detect and estimate reducing substance in a urine specimen.	532
	(b) Estimate glucose concentration in a given urine sample.	
Chapter 26: Estimation of Casein in Milk.....		(535–536)
26.1	Introduction	535
Exp. 26.1:	Isolate pure casein from skim milk.	535
Chapter 27: Oils and Lubricants.....		(537–546)
27.1	Introduction	537
Exp. 27.1:	Determine the iodine value/number of a fat or oil sample.	537
Exp. 27.2:	Determine the saponification value of an oil or fat.	538
Exp. 27.3:	Determine the acid value of an oil.	540
Exp. 27.4:	Determine the flash point and fire point of a given sample of oil by Pensky-Marten's apparatus.	541
Exp. 27.5:	Determine the aniline point of an oil.	543
Exp. 27.5:	Determine the viscosity (in Redwood seconds) of a liquid hydrocarbon (oil) and the effect of temperature on it.	544
Chapter 28: Water Analysis.....		(547–568)
28.1	Introduction	547
28.2	Hardness of Water	548
Exp. 28.1:	Determine the (i) total hardness (ii) permanent hardness and (iii) temporary hardness of a given sample of water by EDTA method.	549
Exp. 28.2:	Determine calcium and magnesium hardness in a given sample of water separately.	552
Exp. 28.3:	Determine chloride ion concentration in a given sample of water by Mohr's method.	553
Exp. 28.4:	Determine the residual chlorine in a potable water sample.	555
Exp. 28.5:	Determine the percentage of available chlorine in a given sample of bleaching powder.	556
Exp. 28.6:	Determine dissolved oxygen present in a given water sample by Winkler's method.	558
Exp. 28.7:	Determine Biological Oxygen Demand (B.O.D.) of a given waste water sample.	560
Exp. 28.8:	Determine the Chemical Oxygen Demand (C.O.D) in a waste water sample.	562
Exp. 28.9 :	Determine the alkalinity of a given water sample.	563
Exp. 28.10:	Determine the acidity of a water sample.	565
Exp. 28.11:	Determine free carbon dioxide in a given sample of water.	566
Exp. 28.12:	Determine iron in a given water sample spectrophotometrically by 1, 10-phenanthroline method.	568