

SALTER LABS®

Salter-Style® Demand Cannula

A Salter-Style cannula specifically designed for use with Dual Port Oxygen Conservation Delivery systems. Adult, Small Adult/Pediatric sizes with 4', 5', 7' supply tubes.

Now...available in an Adult **TLCannula™** configuration.

Soft foam, closed cell cushions are pre-installed on the head piece tubing. This provides extra comfort to the sensitive over-the-ear tissue. The TLC model is offered with 4' or 7' connecting tubes.



Features

- Dual Tubing Set
- Soft, curved Nasal Prongs
- Anatomically curved face piece to fit upper lip
- Over-the-ear style
- NEW TLCannula** model has soft foam cushions on head tubes
- Available in Small Adult/Pediatric and Adult sizes with 4', 5', 7' tube lengths

Benefits

- Provides separate independent tubes for oxygen delivery and sensing inspiratory/expiratory effort
- Improved comfort and anatomical fit
- Provides a "natural" fit, lowers patient anxiety
- A secure fit, easy to put on or take off, even for a reclining patient
- The foam cushions provide extra comfort for the over-the-ear area and aid in holding the cannula securely in place
- Select the appropriate model for your stationary or ambulatory patients

Variety of uses. . . Hospitals, Skilled Nursing Facility, Clinics, Homecare, Assisted Living, Extended Care




slc 08

Salter-Style® Demand Cannula*


A special cannula designed to give both adult and pediatric users of dual port oxygen conservation delivery systems all the comfort and conveniences of a Salter-Style® cannula. Unique dual tubing delivers oxygen on demand through one tube while inspiratory / expiratory effort is sensed through the other tube. Meets dual port device manufacturers' specifications.

Oxygen Delivery



Patient Sensing

Permanent Barrier



Oxygen Delivery

#4807 (pictured)

Now Available...

The NEW Salter-Style DEMAND TLCannula™. Setting a NEW STANDARD in patient comfort and convenience for long term oxygen therapy (LTOT) patients.

The Demand TLCannula is assembled with gray colored, closed cell, soft foam cushions placed over the cannula head tubes. The cushions provide a new dimension in comfort for the sensitive over-the-ear tissue.

The closed-cell foam cushions are also available as separate items for self installation on nasal cannulas. (See SLC-71 Salter Accessories E-Z Wrap™ P/N 1016)



Gray colored, closed cell, soft foam cushions are pre assembled over the head set tubing.

TLCannula P/N 4804TLC

Ordering Information

Salter-Style® Demand Cannula*

Designed for dual port oxygen conservation delivery systems.
For flow rates up to 6 LPM

	Units (per case)	Part Number
Adult , Demand cannula with 4' supply tubes	25	4804
Adult , Demand cannula with 5' supply tubes	25	4805
Adult , Demand cannula with 7' supply tubes	25	4807
Small Adult/Pediatric Demand cannula with 4' supply tubes	25	4824
Small Adult/Pediatric Demand cannula with 7' supply tubes	25	4827
NEW Adult Salter-Style Demand TLCannula with foam cushions and 4' supply tubes	10 or 25	4804TLC
NEW Adult Salter-Style Demand TLCannula with foam cushions and 7' supply tubes	10 or 25	4807TLC

SALTER LABS®

100 W. Sycamore Road - Arvin, CA 93203
Phone: 661-854-3166
1-800-421-0024
1-800-235-4203
Fax: 661-854-3850
Toll Free Fax: 1-800-628-4690
E-mail: salterlabs@us.salterlabs.com
www.salterlabs.com

ISO 9001 & ISO 13485 0482