Hill-Rom VersaCare

Hospital Bed



For safety, quality and comfort. Intelligent, well designed, clinically-effective and easy to use.

OUR PROCESS INCLUDES

- Professionally refurbished medical beds to manufacturer's specifications
- > JCAHO required documentation
- > New paint & labels
- > Fully function tested
- > Reconditioned PC boards and side rail switches
- Reconditioned side rails
- Reconditioned head/foot boards
- > New led plus trend indicat or
- > Fresh hydraulic fluid
- > Rebuilt casters with new components and wheels
- > 1 year warranty (Domestic only)

OPTIONS

- > A wide variety of room accessories
- > IntelliDrive® power-assist transport system
- > Digital head of bed angle indicator
- > Advanced 3 mode bed exit
- > Service Contracts for preventative maintenance
- > 220-Volt transformers

SPECIFICATIONS

- > "Brake not set" alarm
- Frame and mattress integration for superior therapy/safety performance
- ➤ FlexAfoot™ mechanism
- > Enhanced turn assist (on A.I.R. surfaces only)
- > Point of Care® controls Large, easy-to-read design
- ➤ FlexAfoot[™] mechanism
- > Low position surface edge to floor 18" with optional
- > IntelliDrive® power transport 22"
- > Maintains low height when placed in chair
- > High position surface center to floor 37"
- Headboard/push handles maintain same height regardless of patient surface height



983-2022-01-14

