

Tympanic Probes

FEATURES

Shore Medical's Tympanic Probe is designed to be a less intrusive way to accurately measure body temperature. By measuring body temperature from the tympanic membrane - through the ear canal - the patient is not forced to experience the discomfort of placing an esophageal or rectal probe.

- Reliable way to measure body temperature
- Convenient, less-invasive alternative to esophageal or rectal probes
- Soft memory-cell foam conforms to the ear canal, reducing risk of injury to the tympanic membrane while sealing out any ambient air which could affect the temperature reading
- Disposable single-use design reduces the risk of cross-infection and ensures a quality product each time a new unit is used
- Available with various temperature sensors including but not limited to 400, 700 and 10K sensors
- Latex-free
- All sensor wires come coiled and taped with Latex free Co-Adhesive tape*
- Private labeling available



TP 400

*Co-adhesive tape not shown



TP 700

*Co-adhesive tape not shown

ORDERING INFORMATION AND TECHNICAL DATA

Item No.	Description	Case Qty.
TP 400	Probe with 400 series sensor	50
TP 700	Probe with 700 series sensor	50


400 Series Sensor also available with 48" tail.
Temperature interface cables available for most monitors.

Note To Patients:

All the devices shown here are prescription products and must be obtained from a licensed homecare provider or physician. Patients cannot purchase directly from Shore Medical.



400 series sensor



700 series sensor