



Siker Mirror / Polio Blades

English Profile

Laryngoscopes


With Conventional Illumination

SIKER MIRROR

For laryngoscopy of patients with anatomical variations which make intubation with conventional blades difficult or impossible. The basic portions of the blade make an angle of 135° with a stainless steel mirror located at that angle on the flange facing the spatula. The distal portion of the blade is three inches long. Because the mirror inverts the reflected image the operator should exercise caution until familiar with its use. A curved stylette is recommended for use with the Siker blade. Placing the blade in a warm water bath for 10 minutes before use will avoid any mirror fogging problems. Blade made of surgical stainless steel.



REPLACEMENT LAMPS
5-0234-01For Size #4



CAT#	SIZE	DESCRIPTION	A	B
5-3075-04	4	Large Adult	158mm	18mm

POLIO

A modification of the English MacIntosh blade which is offset from the handle at an obtuse angle to allow intubation of patients in respirators, body jackets, and other difficult situations. It is a curved blade type for indirect laryngoscopy. Blade made of surgical stainless steel.



REPLACEMENT LAMPS
5-0234-04For Size #3

CAT#	SIZE	DESCRIPTION	A	B
5-3077-53	3	Medium Adult Angled at 120°	130mm	24mm