

Clinically proven,
noninvasive method
of prophylaxis to
reduce the incidence
of deep vein
thrombosis (DVT)

- For use on patients with high risk conditions such as obesity, advanced age, malignancy or prior thromboembolism.
- Pneumatic controller provides intermittent of pulses of compressed air which sequentially inflate the multiple chambers of the sleeve beginning at the ankle and moving up the leg. Results in a wave-like milking action which empties the veins and results in greatly increased peak blood flow velocity.
- Disposable sleeves available in knee and thigh lengths.
- Thigh-length sleeve available in sterile packaging for use on the draped leg during major hip surgery.
- Equipped with features that assure maximum control, safety and patient comfort.
- Can detect each patient's individual vascular refill rate to deliver the most reliable, efficient, customized compression therapy on the market.

