


Drug Product Contract Manufacturing Services

Robotic Gloveless
GMP Fill & Finish for Sterile Injectables



Formulation: Bottle Mixer

Formulation development and production



Scan or click to learn more



We provide a diverse range of formulation services adapted to small batches, with the capability to work with mixing volumes as low as 50mL, and up to 50L.

Single use approach.

Sterile or non-sterile mixing.

No CapEx for our customers.

Integrated in our IT process reducing human error risk.

High chemical compatibility - Extensive extractable package and very low particle load.


Active substance / Placebo / Diluent / Adjuvant.

Pooling / Dilution.

Powder (API) dissolution.







Specifications & range

Fully closed singe-use system.

Formulation volume range: 50mL to 50L.

Sterility Assurance Level 10-6 (ISO11137)

Animal Derived Component Free (ADCF)

Benefits


Agile mixing: Match your needs without CapEx investment.

Single use: No contamination or mixing error.

Low shear forces: Low degradation and aggregate formation on your sensitive protein.

High chemical compatibility.

Low extractables and particles generation.



#WeRelease