




**Drug Product
Contract Manufacturing Services**

Robotic Gloveless
GMP Fill & Finish for Sterile Injectables




Formulation: Isolator

Formulation development and production



Scan or click to learn more




Key Features

- High containment unit to perform sensitive weighing or formulation activities
- H14 HEPA air filtration.
- H2O2 decontamination system.
- RTP port connection.
- Airlock with dedicated decontamination system.




Specifications & range



- Integrated hydrogen peroxide bio-decontamination device Isosté H+.
- Dynamic aeration and a 20-minute bio-decontamination cycle for the rapid airlock.
- H14 HEPA filters on air intake and exhaust with 99.995% MPPS efficiency.

Benefits

- Superior containment system to protect your product.
- Flexible isolator to match your sterile/contained weighing and formulation activities.



#WeRelease