

Pedigo Procedure Rests

Pedigo brand equipment

The Procedure Rest offers unmatched versatility in a support device. It can be used in a variety of ways depending upon your needs. Features Rotates a full 360 degrees Adjustable height Upholstered for comfort Special density foam provides comfort and support during long procedures Move it to the front as an abdominal support While in front, raise it for use as a wrist rest Lower it for use as an elbow rest Move it to the side for use as an arm rest Move it to the back for use as a back rest Available as an option on the P-6000 Surgeon's Stool and the P-526-GS Exam Stool

FEATURES

- P-6000 Procedure Rest 16 inch diameter seat 21 inch diameter steel base Foot-controlled hydraulic operation Single pedal height adjustment of 20 to 28 inches Powder-coated, 5-caster base with 3 inch casters P-526-GS Procedure Rest 14 inch diameter seat 24 inch diameter aluminum base Hand-operated pneumatic height adjustment 22-1/2 to 24-1/4 inch seat height adjustment 2 inch dual-wheel casters
- Rotates a full 360 degrees Adjustable height Upholstered for comfort Special density foam provides comfort and support during long procedures Move it to the front as an abdominal support While in front, raise it for use as a wrist rest Lower it for use as an elbow rest Move it to the side for use as an arm rest Move it to the back for use as a back rest Available as an option on the P-6000 Surgeon's Stool and the P-526-GS Exam Stool

