

Model	GMM-H01	GMM-H03A	GMM-H05
Weight:			
Measure Range:	5 to 200kg	5 to 200kg	5 to 200kg
Division:	50g	50g	50g
Height:			
Measure Range:	80 to 195 cm	80 to 195 cm	80 to 210 cm
Resolution:	0.5cm	0.5cm	0.5cm
Height Rod:	Able to lower down for easy transport or storage.	Station	Station
BMI Function:	Yes	Yes	Yes
Methods of measurement of Body Fat:	Four points contact electrode, bioelectrical impedance analysis (BIA)	Two points contact electrode at two sides of LCD display, using bioelectrical impedance analysis (BIA)	Four points contact electrode, bioelectrical impedance analysis (BIA)
ON / OFF Voice Function:	Yes	Yes	Yes
Auto zero-Auto off Function:	Yes	Yes	Yes
LCD display size:	102.5 x 67.5mm	70*80mm (blue light)	90*90mm (blue light)
Thermal printer:	No	No	Yes
Scale Platform:	44 x 34 cm	39.5 x 36 cm	42 x 37 cm
Output Function:	No	RS232 port connect to computer with GMM copyright interface software	RS232 port connect to computer with software
Power by:	4*AAA or AC adaptor	AC adaptor	AC adaptor
Material:	ABS Plastic and Aluminium Rod	8mm tempered glass platform, ABS Plastic and Aluminium Rod	ABS Plastic and Aluminium Rod
Total scale weight:	6 KG	15 KG	17 KG