

ANNEX 1

TECHNICAL DATA SHEET

AUTOCLAVABLE BIOHAZARD BAG

INTENDED USE : Autoclavable Biohazard Bags are used for the sterilization of medical wastes before disposal. Bags must be marked with the standard “Biohazard” symbol.

PRODUCT SPECIFICATIONS

PROPERTIES	METHOD/STANDARD
THICKNESS	Minimum 50 µm
MATERIAL	PP (Polypropylene)
TEAR RESISTANCE	ASTM D 1922 or Equivalent standard
DART IMPACT	ASTM D 1709-01 or Equivalent standard
INDICATOR	Indicator patch field turn into dark when autoclaved
TEMPERATURE	Withstand up to 135°C
SIZE Autoclavable Biohazard Bag sizes should be within the ranges mentioned	Small 300x500mm Yellow ± 50mm
	Small 300x500mm Red ± 50mm
	Medium 400x780mm Yellow ± 50mm
	Medium 400x780mm Red ± 50mm
	Large 600x800mm Yellow ± 50mm
	Large 600x800mm Red ± 50mm
	X Large 700x1100mm Yellow ± 50mm
	X Large 700x1100mm Red ± 50mm