







**5. Category 05**

Item no:	Item	Description	Required quantity	
1	<p>Microscope objectives for Philip Harris Microscopes (x 40, high power)</p> <p>MOTIC 40/ 0.65 NA :160/0.17</p>	<p>To be used with Philip harris microscope</p>	<p>10 no's</p>	
2	<p>Microscope objectives for Philip Harris Microscopes (x 100, oil immersion)</p> <p>PLAN 100/1.25, NA: 160/0.17</p>	<p>To be used with Philip harris microscope</p>	<p>10 no's</p>	
3	<p>Test tube (10cm L x 1cm W )</p>	<p>Glass test tubes</p>	<p>100 no's</p>	



4	<p><b>Microscope cleaning and maintenance kit</b></p>	<p><b>Each individual kit should contain:</b></p> <p><b>1.Cleaning kit:</b></p> <ul style="list-style-type: none"> <li>- Instruction manual</li> <li>- Lint-free optical swabs</li> <li>- Lens cleaning solution</li> <li>- Soft Lens cleaning pads</li> <li>- Optical surface cleaning pads</li> <li>- Dust removal compressed air spray</li> </ul> <p><b>2.Maintenance kit:</b></p> <ul style="list-style-type: none"> <li>- Maintenance manual</li> <li>- Lens cleaning solution</li> <li>- Lens paper</li> <li>- Grease</li> <li>- Six-piece jeweler's screwdriver set</li> </ul>	<p><b>2 Kits</b></p>	 
---	---	---	----------------------	---

5	Refrigerator freezer thermometer	<p>Digital Temperature range: -20°C to 60°C Thermometer for refrigerator / freezer with an ISO 17025 A2LA Traceable NIST cert battery included resolution 0.1°C. <b>For an example: Fisher scientific online catalogue no: 11705853</b></p>	7 No's	
6	Bench mounted single hand eye wash- angle head	<p>To be used to wash eyes in accidental spillage of chemicals into the eyes. <b>For an example: Philip harris online catalogue no: B8R01505</b></p>	2 No's	



**Wall-mounted hand-held hose spray**

- Complies with ANSI Standard Z358.1-2009
- For eyes and face or body rinsing
- 12' (3655mm) yellow reinforced thermoplastic hose is self-coiling. Minimum working pressure of 170 psi and burst strength 3 times the maximum working pressure. Universal identification sign and inspection tag included

Plastic perforated spray head for soft-flow spray head for flooding large areas as well as eyes only. Includes protective spray head cover. Integral flow control assures safe, steady flow under varying water supply conditions.

**Valve**

½" stay-open ball valve. Hand-operated by a large highly visible stainless steel push flag handle.

**Wall Bracket**

Stainless steel wall bracket

**For an example:**


**BRAND: BRADLEY®**

BRADLEY ONLINE CATALOGUE



MODEL NO: **S1944022ABC.**

7




**2 No's**






8	<p><b>Apparatus Trolley three tier</b></p>	<p>Fully welded, tubular end frame construction, finished powder coated blue. Shelf heights 210mm, 520mm and 820mm from floor level. The shelves are fixed (i.e. not removable) for with lipped edges Maximum weight limit is 200kg uniform distributed load.</p> <p><b>Size: L1060 x W460 x H900mm.</b></p> <p><b>For an example:</b></p> <p><b>Phillip Harris online catalogue no: B8L02069</b></p>	<p><b>2 trolleys</b></p>	
9	<p><b>Cuvettes ( 1ml/1.5ml)</b></p>	<p>Re-usable 10mm path length cuvettes designed for visible and UV/Vis assays.</p> <p>Clear four-sided cuvettes</p> <p>Molded-in arrow indicates direction of transmission on standard and semi-micro cuvettes</p> <p>Made up of quartz glass.</p>	<p><b>100 No's</b></p>	



10	Slide holder drying holder	<p>To hold slides, holds slides vertically</p> <p>To be used as a receptacle for 3 x 1" or 75 x 25mm slides during drying and cleaning.</p> <ul style="list-style-type: none"> <li>• Light, strong, chemical-resistant high-density polyethylene holder.</li> <li>• 10° angle of slots facilitates clear viewing of desired slide</li> <li>• Beveled edges.</li> </ul> <p><b>Fisher scientific online catalogue no: 11947984</b></p>	4 Racks	
11	Slide folder (Slide storage file)	<p>Slide folder to store 76mm x 26mm slides, 20 slide capacity stackable with transparent cover.</p> <p><b>For an example: Fisher scientific online catalogue no: 12608336</b></p>	10 folders	
12	Grabbit Hand Mitt	<p>Protect hands from hot or cold lab ware. This silicone rubber mitt studded, oversized surface for secure grip. Able to withstand temperatures from -70 to 500°F (-57 to 260°C)</p> <p><b>For an example: Phillip Harris online catalogue no: B8R04579</b></p>	4 No's	



13	<p><b>Foot-switch gas burners</b> (Electrical foot switch, and burner connectable to any type of gas source)</p>	<p>To be used in microbiology laboratory with Stainless steel foot step</p> <p>Stainless steel housing.</p> <p>Connectable to any kind of gas source.</p> <p>Safety Bunsen Burner optimized air-flow system allows a uniform and constant flame even in safety cabinets with horizontal or vertical air flow.</p> <p>Automatically cuts the gas supply after every flame extinguish and checks the burner temperature.</p>	<p>5 No's</p>	
14	<p><b>Multiple inoculator or multipoint inoculator (for inoculation into media)</b></p>	<p>A multipoint inoculator is designed for inoculating agar plates with a maximum of twenty-four strains of bacteria or yeasts simultaneously.</p>	<p>2 No's</p>	