



بِسْمِ اللَّهِ الرَّحْمَنِ الرَّحِيمِ
ދިވެހިރާއްޖޭގެ ސަރުކާރު

Faafu Atoll Hospital

ފ.ނިލަންދޫ، ދިވެހިރާއްޖެ

F.Nilandhoo , Republic of Maldives



Physio items in Faafu Atoll Hospital

#	Equipment Name	QTY	Rate	Total
1	Abdominal Bench	1		
2	Air Splints	1		
3	Basic Ergonomic Hand Exerciser	1		
4	Overhead pulley	1		
5	Diathermy Units, Microwave	1		
6	Diathermy Units, Shortwave	1		



7	Electrotherapy Unit	Electrotherapy Unit (IFT/MS/TENS) LCD/ Programmable	1		
8	Ergometer Bicycles	Ergometer Bicycles (Manual)	1		
9	Exercise Floor mat / Floor Mattress (6 x 2.5)	Exercise Floor mat / Floor Mattress (6 x 2.5)	1		
10	Exercise Chair	Exercise Chair / Quadriceps Chair with Back Rest	1		
11	Shoulder Wheels (Mariner's Wheel)	Shoulder Wheels) 360degrees (Big)	1		
12	Foam Roller	Foam Roller	1		
13	Interferential Therapy	Interferential Therapy Digital With 100 Programs (Metal Cabinet) Operating Voltage: 200/230V 50Hz AC Medium Frequencies: 2 KHz & 4 KHz Base: 0-150 Hz Sweep : 0-100 Hz Sweep Programs: a) Rectangular 1/1 sec. b) Trapezoidal 1/6 sec. c) Triangular 6/6 sec. Therapy Modes : a) 4 Pole Interferential b) 4 Pole Vector c) 2 Pole Modulated Intensity Control: 0-100 MA Continuously Variable Intensity Balance: To equalize 2 Channels Timer : 99 minutes Digital with automatic alarm Output type : 100% isolated	1		



14	Neuro Muscular Stimulator	Neuro Muscular Stimulator (Digital) Main Input :200/230V 50Hz AC Therapy Mode : Faradic, Surge Faradic, Galvanic Int., Galvanic Pulse Width Faradic: 06msec. Int. Galvanic: .01,.03, .1, .3, 1,3,10,30, 100, 300, 1000 msec. Frequency Faradic: 50Hz Galvanic: 0.3,1&3 Pulse/sec.	1		
15	BMI Machine	BMI Machine (Body Composition Analyzer)	1		
16	Light Infrared	Light Infrared / IR Lamp / Tabletop Model	1		
17	Packs Heater (Hot & Cold Packs)	Packs heater (hot & cold packs)	1		
18	Parallel Bar	Parallel Bar (Height and width Adjustable)	1		
19	Shockwave Unit	Shockwave Unit Technology: Electronically Controlled Ballistic Air Pressure Compressor Shockwave Type: Radial Shockwave Shock Intensity: 60MJ -185 MJ Frequency Mode: 1 Hz -16 Hz Pulse Mode: Continuous 4/8/16 Pulses Applicator Heads: Dia 6mm,10,15,20,25mm Power: AC100 V -110V / AC 220-230V, 50Hz / 60Hz	1		
20	Stability Ball	Stability Ball / Physio Ball	1		



21	TENS Therapy Equipment	TENS Therapy Equipment Main Input: 200/230V 50Hz AC Output Channel : Four Timer: 0-99 Min Mode Selection: Continuous, Burst1, Burs t 2, PWM, Sweep, combi. Frequency: 2,10,50,100,125,200 Hz	1		
22	Cervical & Pelvic Traction Unit	Cervical & Pelvic Traction Unit Lumbar cum cervical) separate switch, 1kg. Steps in cervical MAIN SUPPLY 220 V 50 HZ CONSUMPTION 85 Watts LUMBER FORCE 15-45 Kg with Doublers Up to 90Kg CERVICAL FORCE 3,4,5,6,7,8,9,10,11,12,13 kg MODE SELECTION Separate Switch TIMER DIGITAL 0-99 minutes with auto alarm HOLD TIME 20,40,60,80 Secs. REST TIME 1,5,10 secs.	1		



23	Ultrasound Units, Physical Therapy	Ultrasound Units, Physical Therapy 1 MHz (LCD) Power input :90VA 200V 50Hz AC Ultrasonic Frequency :1 MHz Ultrasonic Out Put: Maximum 15 Acoustic Output Mode: Continuous and Pulsed Pulse Ratio: 1:1, 1:3, 1:5, 1:7, 1:9, 1:12,1:16 Treatment time: Fixed according to programs which vary from 3 to 15-minutes Timer: 0-99 min auto alarm	1		
24	Wax Bath	Wax Bath	1		
25	Ankle/ Wrist Weight cuffs	Ankle/ Wrist Weight cuffs 0.5kg	1		
26	Ankle/ Wrist Weight cuffs	Ankle/ Wrist Weight cuffs 1kg	1		
27	Ankle/ Wrist Weight cuffs	Ankle/ Wrist Weight cuffs 2kg	1		
28	Ankle/ Wrist Weight cuffs	Ankle/ Wrist Weight cuffs 3kg.	1		
29	Ankle/ Wrist Weight cuffs	Ankle/ Wrist Weight cuffs 4kg.	1		
30	Ankle/ Wrist Weight cuffs	Ankle/ Wrist Weight cuffs 5kg	1		
31	Staircase with Ramp	Staircase with Ramp (Wooden)	1		
32	Foot stool	Foot stool	1		
33	Racks, Storage for small Exercise	Storage Rack steel	1		
34	Tables , Physical Therapy	Tables, Physical Therapy / Wooden Couch with Mattress / Drawer	1		



35	Tables , Traction W/ Traction Unit	Tables, Traction W/ Traction Unit with 3-Fold Bed Lumbar cum cervical) separate switch, 1kg. Steps in cervical MAIN SUPPLY 220 V 50 HZ CONSUMPTION 85 Watts LUMBER FORCE 15-45 Kg with Doubler's Up to 90Kg CERVICAL FORCE 3,4,5,6,7,8,9,10,11,12,13 kg MODE SELECTION Separate Switch TIMER DIGITAL 0-99 minutes with auto alarm HOLD TIME 20,40,60,80 Secs. REST TIME 1,5,10 secs.	1		
36	Treatment/ Examination Couch-hydraulic	Treatment/ Examination Couch-hydraulic Hi-Low Table (Electrical / with Remote)	1		
37	Hand Dryer	Electric Wall hanging	1		
				Total	
				GST	
				Total	

