



FORMULATION DEVELOPMENT OF OXYCODONE HCL SOFT GELATIN CAPSULES

*S. Gumudavelli

Research Associate, Ascent Pharmaceuticals Inc, 400 S. Technology Dr, Central Islip, NY 11790.

***Corresponding Author: Dr. S. Gumudavelli**

Research Associate, Ascent Pharmaceuticals Inc, 400 S. Technology Dr, Central Islip, NY 11790.

DOI: 10.20959/ejbps20196-7446

Article Received on 26/03/2019

Article Revised on 16/04/2019

Article Accepted on 06/05/2019

ABSTRACT

The objective of this work is to develop simple, elegant, liquid-filled immediate release soft gelatin capsules containing Oxycodone HCl intended for the management of moderate to severe pain in patients who require treatment with an oral opioid analgesic. The present work describes the design of formulation and manufacturing process of Oxycodone HCl 30 mg capsules. Manufacturing of oxycodone HCl Capsules, 30 mg, consists of two concurrent manufacturing processes: gel mass and fill material, which converge during the encapsulation process. Capsule fill material was formulated using citric acid and high molecular polyethylene glycol matrix system. The gel mass was prepared in gelatin as shell forming material, and glycerin was used as plasticizer at a 1: 0.5 ratio. The prepared soft gelatin capsules were evaluated for water content in fill material, loss on drying of shell, weight variation, burst strength, assay, and dissolution. The capsules had a fill weight of 300 mg, disintegration time less than 5 minutes, and drug release rate >85% in 15 minutes. Furthermore, Texture analysis showed the capsules had burst strength of 125.0 N with an elasticity value of 2.23 mm. With this we, conclude that oxycodone HCl immediate release soft gelatin capsules can be manufactured at commercial scale with similar composition and improved quality attributes of other existing solid dosage forms such as tablets.

KEYWORDS: Oxycodone HCl, Soft Gelatin Capsules, Encapsulation.

INTRODUCTION

The one-piece softgel capsule features an easy-to-swallow shell that is tamper-resistant, masks taste and smell, and available in a broad range of unique shapes and colors.

Wide varieties of lipid-based fill excipients are used to improve the bioavailability of poorly water-soluble compounds, reduce food effects, and improve drug stability. The wide range of shapes and colors available for the shell meet objectives for branding, line extension and lifecycle management. Additionally, softgel capsules offer unparalleled dosage-form flexibility: Encapsulation of a wide variety of fill formulations, wide array of sizes, shapes, colors and forms, oral, topical, ophthalmic and vaginal administration, and softgel capsules address a broad array of formulation challenges, often incorporating a variety of pharmaceutical-approved lipid-based excipients. Formulation flexibility is a hallmark of softgel formulation and manufacturing technology. Furthermore, it allows multi-fold improvement in bioavailability and reduction of food effect, assures dose uniformity, provides safe handling of high potency API, and constitutes a suitable dosage form for oxygen and/or light sensitive API.^[1]

MATERIALS AND METHODS

Oxycodone HCl, citric acid monohydrate, polyethylene glycol 600, gelatin, glycerin all are commercial grade and obtained from reputed vendors.

Preparation of Oxycodone HCl Fill Material

Oxycodone HCl fill material was prepared by dissolving in citric acid solution followed by mixing with pre-molten polyethylene glycol 600 and mixed continuously for 15-20 minutes under continuous nitrogen purging. The resultant clear solution was filtered through 10 micron filter and stored in a S.S. vessel until encapsulation. Oxycodone HCl fill material was prepared with varying amount of citric acid and polyethylene glycol 600. Formulation details are given in Table 1.

Table 1: Composition of Oxycodone HCl Fill Material for Oxycodone HCl Softgel Capsules, 30 mg.

Ingredient(s)	Formulations Code			
	OXY-1 (mg/unit)	OXY-2 (mg/unit)	OXY-3 (mg/unit)	OXY-4 (mg/unit)
Oxycodone HCl	30.00	30.00	30.00	30.00
Citric Acid	1.00	1.50	2.00	2.50
P.Water	60.00	60.00	60.00	60.00
PEG 600	209.00	208.50	208.00	207.50
Total	300.00	300.00	300.00	300.00

Preparation of Gel Mass

Gel formulation was prepared with gelatin and glycerin as a shell forming material and plasticizer at a concentration of 42% and 22%, respectively. Gel mass was prepared by cooking gelatin in a mixture of glycerin and purified water for approximately 65 ± 10 minutes at 80°C until a homogenous mass was obtained. Cooked gelatin mass was deaerated to remove air bubbles. Water content in the gel mass was maintained between 20% - 30% (Target 25%). After the desired water content was achieved, the gelatin mass was discharged to preheated gelatin holding tanks maintained at $60 \pm 5^\circ\text{C}$ till encapsulation.^[2]

Encapsulation

Oxycodone HCl soft gelatin capsules, 30 mg were prepared by using rotary die machine (Sky Pharma, Korea). Light mineral oil (Miglyol 812 N) was used for machine and ribbon lubrication as well. Capsules were dried in two stages: primary drying was carried out in tumbles using mineral oil-lecithin solution wetted towels at room temperature for 30-45 minutes in continuous motion and secondary drying was carried out in drying tunnel at controlled temperature ($28^\circ\text{C} \pm 2^\circ\text{C}$) and humidity (NMT 25% RH) for 72 hours to achieve the desired moisture content.^[3]

EVALUATION OF SOFT GELATIN CAPSULES

The soft gelatin capsules filled with oxycodone HCl fill material were subjected to different tests to evaluate for various parameters.^[4]

Capsule Description

Physical appearance of the product was visually determined by spreading ten capsules on a clean white surface and inspected for their physical attributes.

Capsule Dimensions

The dimensions of the capsules were measured using a Vernier Caliper.

Weight Variation of Capsules

Weight variation was determined by cUSP general Test Method <95>.^[5]

Capsule Shell and Fill Weight

The capsule was sectioned into halves, and the fill material was drained. The shell was wiped with lint free paper towels to remove any remaining fill material before weighing. Capsule fill weight was determined by

subtracting the shell weight from the whole capsule weight.

Moisture Content

Moisture content for shell and fill material was determined as below:

Shell Moisture Content

Shell moisture content of softgels was determined according to Test Method FP131G. Softgel capsule shells were sectioned into small pieces, and heated at 105°C for 15 minutes (standard drying) using an O'HAUS Moisture Analyzer. Data was reported as Loss on Drying (LOD), percent weight/weight (%w/w).

Fill Material Moisture Content

The water content in the capsule fill material was determined by Karl Fischer Titration, cUSP <921> Method 1a. Data obtained by this method was noted as "KF".^[5]

Capsule Shell Thickness

The capsule shell thickness was measured using a Vernier Caliper.

Capsule Burst Strength

Capsule burst strength (n=6) was performed with the CS225 Chatillon force measurement (Manufacturer: Ametek) fitted with a flat hardness probe 1.5 inch in diameter. The amount of force (N) and distance (mm) the probe traveled until the capsule burst was recorded.

Drug Content

Drug content was determined by opening ten capsules and the fill solution equivalent to 5.0 mg of oxycodone HCl, to a 50ml volumetric flask, add 10ml methanol and added with solvent mixture (methanol: water=4:1) to volume. The resulting solution was passed through a $0.45\mu\text{m}$ membrane filter. Sample solution: Dilute the 1ml stock solution with solvent mixture to a 10ml volumetric flask. Further, this stock solution 1 mL was diluted to 10 mL with mobile phase. Finally drug content was estimated using HPLC (Agilent Technologies) at 214 nm (Mobile phase: 7mM potassium dihydrogen phosphate aqueous solution (containing 0.1% phosphoric acid, 0.1% n-nonylamine, with potassium hydroxide solution (1 $\rightarrow$ 2), the pH was adjusted to 4.9 ± 0.1 , Column: Welchrom C18 ($5\mu\text{m}$, $250 \times 4.6\text{mm}$), injection volume: 20 μl).

Dissolution Study

Drug release study was carried out using USP dissolution rate test Apparatus I. The study was carried out at $37 \pm 5^\circ\text{C}$ and 50 rpm for 45 minutes in 900 mL of water. After specified time interval, 10 mL of sample was withdrawn from the vessel. The sample was filtered through $0.45 \mu\text{m}$ filter and the drug contents were estimated using HPLC at 214 nm.

RESULTS AND DISCUSSION

Evaluation of Oxycodone HCl Soft Gelatin Capsules, 30 mg revealed that the 13.74×4.65 mm dimensional Oblong, transparent softgel capsule contains the active drug Oxycodone HCl in clear liquid form at a

concentration of 30 mg per average fill weight of 453 mg. Physical evaluation showed that the average weight of 10 capsules was 451 mg; shell weight and thickness were 152 mg and 0.80 mm respectively. The moisture content in the fill and shell was 3.08% and 6.8% respectively. Texture analysis showed that the capsules had a burst strength of 100.8 N with an elasticity value of 2.23 mm. Oxycodone HCl Soft Gelatin Capsules, 30 mg had a potency of 99.6%, 0.34% related compounds. All four Softgel formulations containing varying amount of citric acid showed similar physicochemical properties. The drug release rate at 45 minutes was 98.2%. Results are presented in **Table 2**.

Table 2: Evaluation of Oxycodone HCl Soft Gelatin Capsules.

Attribute	Formulations Code			
	OXY-1	OXY-2	OXY-3	OXY-4
Description	Clear transparent oblong capsules containing clear liquid	Clear transparent oblong capsules containing clear liquid	Clear transparent oblong capsules containing clear liquid	Clear transparent oblong capsules containing clear liquid
Dimensions (mm)	13.74 X 4.65	13.77 X 4.62	13.73 X 4.60	13.71 X 4.64
Fill Weight (mg)	300.01	302.00	300.25	299.98
Shell Weight (mg)	151.02	150.20	149.25	148.52
Total Capsule weight (mg)	451.00	453.01	449.78	449.25
Burst Strength (N)	121	135	128	125
Elasticity (mm)	2.21	2.23	2.20	2.18
Disintegration Time (min)	<7	<7	<7	<7
Drug Content (%)	98.6	99.8	100.25	100.01
Dissolution (in 15 minutes)	>99	>95	>95	>94

CONCLUSION

The proposed formulation showed optimal drug release, implicating that this method can be used for treatment of moderate to severe pain in patients who required oral mediated opioid analgesic. Furthermore, the proposed formulation shows that Oxycodone HCl soft gelatin capsules can be prepared on a commercial scale with great viability and improved patient compliance and bioavailability.

ACKNOWLEDGEMENTS

The author is thankful to M/s. Ascent Pharmaceuticals Inc, USA for providing the opportunity to do research work.

REFERENCES

1. Leon Lachman, Herbert. A. Lieberman. Joseph L. Kanig., The Theory and Practice of Industry Pharmacy, 3rd ed., Varghese Publishing House, Mumbai, 1991; 410.
2. Gabriele R., Formulation and physical properties of soft gelatin capsules, 3rd ed., Whitehouse Station and Co. Inc. NJ., 2003; 201-212.
3. Aulton MV, Pharmaceutics- The Science of Dosage Form Design, 2nd ed., Churchill Livingstone, Elsevier Limited, 2004; 246.
4. HTTP:// WWW.LOCUMUSA.COM. Product development guide pre-formulation-soft gelatin

capsules. International Journal of Generic Drugs- ISSN 0793694 X US/ Canada ISSN 0793 7784 Euro, ISSN 0793 7822. Pacific rim, accessed on 10th September 2010.

5. United States Pharmacopoeia, 30th ed., United States Pharmacopoeial Convention Inc., 2007; 1965; 2325.