

About Bodyshape

TSB 360

Plenty of innovations in one device

An innovative fitness device that fuses cardio with skincare benefits. Closed capsule with a built-in treadmill (optionally with a stepper or a bicycle), in which the training is supported by infrared heating therapy, under pressure vacuum and other functions increasing the effectiveness of exercises, encouraging weight loss and cellulite reduction. Modernist look, intuitive operation, innovative technology - these features make Bodyshape give a completely new quality to slimming training.



Results

Reach your ideal size

Are you constantly trying to reach your ideal size? Is your diet ineffective? You don't have to worry anymore. Body Shape Multivision will help you burn up to **3,000 calories during one training session**. The purpose of all applied technologies is to increase the energy consumption of muscles. The energy generated as a result of burning the fatty tissue. Body Shape Multivision will help you lose excess kilograms without any special diet.



Design



Timeless design

Our rich experience in design of slimming capsules allowed us to create a product with a futuristic design, capable of drawing attention of any customer. Precisely designed elements of the device serve aesthetic purposes and distinguish the Bodyshape Multivision from other models available on the market.

Underpressure

Increase the effectiveness of training

Generated by a specially designed vacu pump, which increases several times the efficiency of each training session.



**Maximum
negative
pressure**

30 mbar

Underpressure

supports **fat burning**

eliminates cellulite **up to 80%**,

reduction in circumference from
8 to 20 cm

slims and shapes the figure

improves blood circulation
(smooth, elastic and firm skin)

improves the functioning of the
circulatory and **lymphatic**
systems



**Vacu pump
power**

1500_w

Infrared IR

Deep thermal healing

Modern radiators provide the full spectrum of infrared radiation IR-A, IR-B and IR-C. They deliver healing thermal radiation directly to the deep subcutaneous layers. In this way, they accelerate the metabolic process, dilate blood vessels, improve blood flow and oxygenate tissues



Infrared radiators power

700_w

Infrared IR

cleansing the body of toxins

treatment of **injuries** and
rheumatism

discharge of stress tensions

reducing the feeling of fatigue

fat burning

skin firming



**Maximum
temperature**

50°C

122°F

Features

Additional functions combined in 1 device



Aromatherapy

Stimulates the body's immune system to function and restores its natural balance. The use of properly selected essential oils supports: regeneration of the epidermis, biological renewal of the body, fighting infections, circulation and oxygenation of the blood.



Ozonotherapy

Uses the natural properties of ozone, i.e. a strong bactericidal and fungicidal stimulant. This element also has significant oxidizing properties. The mixture of oxygen and ozone is the most effective weapon in the fight against microorganisms that attack the human body. Ozone therapy is used in the case of bacterial infections, to stimulate the circulatory system and to improve cell metabolism.

Features



Colortherapy

Based on the use of properly selected colors that affect the functioning of the body. A specific color of light stimulates the brain to specific hormonal reactions, and the hormones produced in the glands are transmitted through the bloodstream to internal organs, activating them and regulating their work. Chromotherapy increases immunity, gives additional energy, restores the body's balance, calms down and increases the overall level of satisfaction.



Collagen lamp

The simulation of skin cellular activity by means of red light with an appropriately selected wavelength. The light penetrates deep into the skin and additionally stimulates the processes of its self-healing and reconstruction. In this way, the healing of wounds and the disappearance of scars is accelerated. The aging process of the skin also slows down - collagen and elastin fibers are rebuilt.



Ionization

Our environment today is full of positive ions produced by electronic equipment. Being in such an environment increases sleepiness, worsens general well-being and increases fatigue. This function enriches the air with health ions, improving the psychophysical functioning of the body, which also receives a portion of additional energy.



Electrostimulation

Stimulates the muscles with low-frequency current. Using special electrodes, it makes it possible to work on selected parts of the body - supplies hard-to-reach places with necessary substances. Electrostimulation reduces tension and improves the process of nourishing the muscles, improves blood supply to tissues, and also has an analgesic effect.



Treadmill

**Maximum
speed**

12 km/h

Manufactured directly by Vacu Activ on special components. Equipped with a powerful **three-phase brushless motor** that guarantees failure-free operation.

**Orthopedic
running belt**

3.2 mm

**Maximum
load**

180 kg

**Maximum
incline**

10°

Multivision

Modern and intuitive interface

During training modify parameters according to your needs. Control vacuum massage level, treadmill incline, infrared heating sensors, speed and other functions.

Information display

Speed and inclination angle of the treadmill, vacuum value, infrared intensity, calories burned counter, time remaining until the end of the training, heart rate level, distance traveled.

Bodyshape



Multivision

Training with virtual track

This innovative feature enriches the training with a virtual view, synchronized with your training. Imagine running along the coast of Los Angeles, through the streets of Venice, among the desert sands or forest wilderness are just a few of the many possibilities **offered by Multivision.**

Bodyshape



Technical specification

Voltage	230VAC / 50Hz / 1PH
Power consumption	2.5kW
Functions	
Infrared	Temperature adjustment: 25 – 50°C, Two infrared radiators: 350W each, Protection against scalding, The range of infrared radiation: 1.4 - 14µm – close to the solar radiation,
Vacuum	5 predefined programs: 0 – 30mBar, A function to adjust the vacuum manually, A vacuum pump, 1500W,
Chromotherapy	Based on LED technology, Wireless remote control allowing to choose the colour, Predefined programs for automatic colour change,
Aromatherapy	An automatic dosing station with essential oils for the THERMAL operation mode, An option to turn on/of the aromatherapy during the training,

Air ionization	Negative ion emission: 9000000 ions/cm ³ ,
Ozone	Ozone emission: 0.05ppm,
Collagen	A single LED lamp, 30W, Wavelength – 620nm, An option to turn on/of during the training,
Electrostimulation	An independent electrostimulation system, Predefined programs, The electrodes have the disinfection function.



Bodyshape


Dimensions



Width
123 cm 48.5 in



Height
159 cm
62.6 in



Length
228 cm 90 in



Device weight
230 kg 508 lbs