





A. -ESI MRM(682.4/466.2)

LpxH minus IPTG





B. -ESI MRM(710.4/466.2)

LpxH minus IPTG





C. -ESI MRM(682.4/466.2)

Standard Lipid X





D. -ESI MRM(710.4/466.2)

Standard Lipid X



E. -ESI Product Ion Scan (682.4) (rt 3.6-3.7 min, CE -50)

LpxH minus IPTG



F. -ESI Product Ion Scan (710.4) (rt 3.7-3.8 min, CE -50)

Standard Lipid X

