

COLD CURED FOAM



- Two Part Polyurethane Components.
- Pre-measured components are poured into the mold cavity and pressurized.
- Consistent density throughout the cushion.
- Cushions are cured for 48 hours before assembly.
- Cells inside the cushion are more defined.
- Moisture resistant outer seal.
- Regains shape almost immediately.

SLAB FOAM (Dual Density)



- Often referred to as “Batch” or “Bun” foam.
- Made in vats 8’ x 8’ x 40’.
- Density can vary throughout the “Bun”.
- The “Bun” is sliced after it’s been made, into the shapes for the cushions.
- Many cells are split and “popped” during the making and cutting of slab foam.
- No moisture resistance.
- Can take time to regain shape.



FREQUENTLY ASKED QUESTIONS

DOES COLD CURED FOAM CONTAIN CHEMICALS?

Yes, all foam is made from a chemical process; however, no harmful chemicals are present in Cold Cured Foam, it is actually a preferred method as no gases are emitted, since the foam is never heated.

WHY IS COLD CURED FOAM NOT USED IN THE USA?

This would be a false statement. Cold Cured Foam is used extensively in automobiles, theatre seating, and in many high end furnishings throughout North America.

IS COLD CURED FOAM JUST A “CHEAP” WAY TO MAKE FOAM?

No, Cold Cured Foam is actually more expensive to produce than slab foam. ComforTek Seating chooses to offer this as an option, as we believe in its merits and longevity.

IS COLD CURED FOAM THE OLD WAY OF MAKING FOAM?

No, Cold Cured Foam is one of the newest methods of producing foam. Slab foam has been around for much longer, as it was a simple process of making large batches of foam, which could be cut to use with many different products.

IF COLD CURED FOAM IS BETTER, WHY DO YOU OFFER SLAB FOAM?

While we encourage the use of Cold Cured Foam, we accept that it is not for everyone. We ensure a high quality in our slab foam as well, providing dual density foam. Which consists of two different density foams placed together to create a comfortable seat.

ADDITIONAL QUESTIONS ON COLD CURED FOAM?

If you still have additional questions regarding Cold Cured Foam or any of its processes, please contact customer service at (888)-678-2060 or by email at info@comfortekseating.com



...better by design.

www.comfortekseating.com
info@comfortekseating.com
1 (888) 678-2060

