

Surgical face mask (BFE \geq 98%)
Children. 6-9 and 10-12 years



Surgical face mask.

Type IIR

(BFE \geq 98%)



Children (6-9 and 10-12 years)

Información técnica

Typology: Disposable face mask for kids
Standard: CE: EN 14683-2019 + AC:2019,
Size:

150 x 65 mm (6-9 year)
170 x 70 mm (10-12 years)


Bacterial filtration efficiency: \geq 98%

Features: Splash resistant

Structure: Mask body, nose clip and mask belt

Validity period: 24 months

Allergy warning: Non-woven textile



Fitting instructions

1. With the blue side facing outwards, place it over the nose and mouth, with the nose clip facing upwards. Then place the loops behind the ears.
2. Gently press on the nose clip in order to adjust it to the nose.
3. Pull the lower part of the mask so that it is beneath the chin and check that it fits well.



Important information


- Hypo-allergenic waterproof polypropylene material, including flexible nasal clip.
- CE Test Report issued by TÜV SÜD.
- Conditions for storage: To be stored in an area with relatively low humidity, free from corrosive gases and well ventilated, avoiding high temperatures ($>40^{\circ}\text{C}$).

**Surgical face mask (BFE ≥ 98%)
Children (6-9 and 10-12 years)**



Surgical face mask. Type IIR (BFE ≥ 98%)

Children (6-9 and 10-12 years)



CHARACTERISTICS	STANDARD	RESULTS
Bacterial Filtration Efficiency (BFE)	EN 14683:2019+AC:2019	BFE ≥ 99%
Differential pressure (ΔP)	EN 14683:2019+AC:2019	ΔP < 38 Pa/cm ²
Microbial cleaning	EN 14683:2019+AC:2019	CFU/g < 20
Synthetic blood penetration	ISO 22609:2004	Not detected

Source: Test Report No.: 200193327 (TÜV SÜD)

10 pcs bag specifications

Quantity: 10 units/bag

SKU: MASK-ST-00025.0010

Dimensions: 215 x 115 mm

Weight: 21 g

Size: 6 - 9 years



Size: 10 - 12 years



Box specifications

Quantity: 100 bags/box

Total quantity: 1000 units

Dimensions: 500 x 370 x 230 mm

Some product features may vary depending on the version available.

