

**Surgical face mask  
Type IIR (BFE  $\geq$  98%)**



## Surgical face mask Type IIR (BFE $\geq$ 98%)



### Technical information

**Typology:** Medical face mask (Type IIR)

**Model:** 175 mm x 94 mm

**Standard:** EN 14683-2019 + AC:2019

**Size:** 175 mm x 94 mm


**Bacterial filtration efficiency:**  $\geq$  98%

**Features:** Type IIR: maximum surgical quality  
Splash resistant

**Structure:** mask body, nose clip and mask belt

**Validity period:** 24 months

**Allergy warning:** Non-woven textile



### Fitting instructions

1. With the blue side facing outwards, place it over the nose and mouth, with the nose clip facing upwards. Then place the loops behind the ears.
2. Gently press on the nose clip in order to adjust it to the nose.
3. Pull the lower part of the mask so that it is beneath the chin and check that it fits well.



### Important information

- Hypoallergenic waterproof polypropylene material, including flexible nasal clip.
- CE Test Report issued by TÜV SÜD.
- Conditions for storage: To be stored in an area with relatively low humidity, free from corrosive gases and well ventilated, avoiding high temperatures ( $>40^{\circ}\text{C}$ ).





Surgical face mask  
Type IIR (BFE  $\geq$  98%)



## Surgical face mask Type IIR (BFE $\geq$ 98%)



CHARACTERISTICS	STANDARD	RESULTS
Bacterial Filtration Efficiency (BFE)	EN 14683:2019+AC:2019	BFE > 99.4%
Differential pressure ( $\Delta P$ )	EN 14683:2019+AC:2019	$\Delta P < 32.4 \text{ Pa/cm}^2$
Microbial cleaning	EN 14683:2019+AC:2019	CFU/g < 20
Synthetic blood penetration	ISO 22609:2004	Not detected

Source: Test Report No.: 60379483001 (TÜV SÜD)

### 10 pcs bag specifications

Quantity: 10 units/bag

Dimensions: 240 x 135 x 23 mm

Weight: 36 g

### Master carton specifications

Quantity: 100 bags/master carton

Total quantity: 1000 units

Dimensions: 540 x 420 x 290 mm

Weight: 5.5 kg

SKU: MASK-ST-00032.0010

