

Surgical Face Mask
Type IIR (BFE $\geq$ 98%)



Surgical face mask. Type IIR (BFE $\geq$ 98%) Kids

From 6-9



Información técnica

Typology: Disposable face mask for kids

Standard: EN 14683-2019 + AC:2019.Type IIR

Size: 150 x 65 mm (6-9 year)


Bacterial filtration efficiency: $\geq$ 98%

Features: Splash resistant

Structure: Mask body, nose clip and mask belt

Validity period: 24 months

Allergy warning: Non-woven textile



Fitting instructions

1. With the blue side facing outwards, place it over the nose and mouth, with the nose clip facing upwards. Then place the loops behind the ears.
2. Gently press on the nose clip in order to adjust it to the nose.
3. Pull the lower part of the mask so that it is beneath the chin and check that it fits well.



Important information

- Hypo-allergenic waterproof polypropylene material, including flexible nasal clip.
- CE Test Report issued by TÜV SÜD.
- Conditions for storage: To be stored in an area with relatively low humidity, free from corrosive gases and well ventilated, avoiding high temperatures ($>40^{\circ}\text{C}$).

Surgical Face Mask
Type IIR (BFE $\geq$ 98%)



Surgical face mask. Type IIR (BFE $\geq$ 98%) Children

From 6-9 years



Surgical Face Mask
Type IIR (BFE ≥ 98%)



Surgical face mask. Type IIR (BFE ≥ 98%) Children

From 6-9 years



CHARACTERISTICS	STANDARD	RESULTS
Bacterial Filtration Efficiency (BFE)	EN 14683:2019+AC:2019	BFE ≥ 99%
Differential pressure (ΔP)	EN 14683:2019+AC:2019	ΔP < 38 Pa/cm ²
Microbial cleaning	EN 14683:2019+AC:2019	CFU/g < 20
Synthetic blood penetration	ISO 22609:2004	Not detected

Source: Test Report No.: 200193327 (TÜV SÜD)

10 pcs bag specifications

Quantity: 10 units/bag

Dimensions: 215 x 115 mm

Weight: 21 g

SKU: MASK-ST-00025.0010

Size: 6 - 9 years



Master carton specifications

Quantity: 100 bags/master carton

Total quantity: 1000 units

Dimensions: 500 x 370 x 230 mm

Some product features may vary depending on the version available.

