

LOOKING AFTER YOUR BGV'S EARS

The inside of your BGV's ears should be clean, with no odour. A simple, regular check will help to keep them clean and infection free.

What you need

- * Smooth-jawed locking surgical scissors
- * Pluck It or Thornit powder to give your fingers some grip
- * Ear cleaner
- * Cotton wool pads



With your fingers, coated with powder for grip, or smooth-jawed scissors with a lock, gradually tease out small strands of hair.



DO NOT probe too deep.



Put a small amount of a simple canine ear cleaning solution direct into the ear canal - or onto some cotton wool and insert this gently on the end of your small finger.

DO NOT use ear cleaning solution on a ruptured ear drum or on any open sore or wound in this area.



Continue to wipe the inside of the ear canal until your cotton wool is clean.



Finally, wipe around the inside of the ear flap with some solution on clean cotton wool until your BGV's entire ear is completely clean.

CHECKING YOUR BGV'S EARS REGULARLY CAN PREVENT PROBLEMS!