



Hill Rom Total Care Bed

6574 N State Road 7 – Suite 228
Coconut Creek, FL 33073
877-446-8700 – alereusa@gmail.com



Hill Rom Total Care Refurbished Bed w/New Mattress
Order #132278 \$3685.00

The TotalCare Bed features both manual and hydraulic controls with which the caregiver can easily transition the patient surface into numerous positions from either side. Voltage switches operate the electric controls and are accessible to both caregiver and patient. The manual controls allow the bed to be operated and adjusted during transport when electric power is not available. Head, knee, foot, and bed up/down functions are easily operated through the hydraulic manual control, and a lock out feature is available to caregivers in order to disable the bed system articulating functions. Due to its versatility and ease of use in emergency situations, this bed is ideal for acute patient care, medical/surgical care.

Hill-Rom TotalCare Hospital Bed Features:

- Side rail controls.
- Removable Headboard.
- Removable Footboard.
- Locking Wheels (4).
- Brake pedal on both sides of bed.
- CPR Release - lowers head down from elevated position to a flat position.
- CGI Interface.
- Manual and electric controls.

Scale System (optional): Can weigh patient and keeps history of patient weights.

Bed Exit Alarm (optional).

Hill-Rom TotalCare Hospital Bed Specifications:

Overall Width (Side Rails Up): 40".

Overall Width (Side Rails Down): 36.5".

Overall Length: 91.5".

Head Elevation (Maximum): 75°.

Knee Elevation (Maximum): 20°.

Trendelenburg (Maximum): 12°.

Reverse Trendelenburg (Maximum): 12°.

Safe Working Load (Patient and Equipment): 500 lbs.

Casters: 5".

Sleep Surface: 35" x 83".

Height Range: 14.5 - 33.5".

Maximum Scale Weight (optional): 400 lbs (181 kg).