

Bed and Stretcher Refurbishing Process

Beds and stretchers are completely disassembled

Bed frames are sandblasted to remove all old paint

Beds are electrostatically painted

Stretchers are steam cleaned

Casters are removed and replaced

Reassembly begins with new plastics

All motors and functions are biomedically tested

All circuit boards are electronically tested

Parts are replaced to look and work like new

Beds and Stretchers are cosmetically restored to “like new” condition

Beds and Stretchers receive final inspection and certified

Beds and stretchers are delivered “patient ready”

Refer to separate download for our industry exclusive 3-yr warranty

WE BACK UP THIS PROCESS WITH A FULL COVERAGE 3-YEAR WARRANTY